



2800T

All-in-one GC Autosampler: Headspace, Liquid and SPME



INSTRUCTION MANUAL

Version: h

Original instructions



NOTICE: For safe and proper use, follow these instructions. Please keep them for future reference.

Copyright © 2025, HTA s.r.l.

All rights reserved in all countries.

Any distribution, alteration, translation or reproduction, partial or total, of this document is strictly prohibited unless with prior authorization in writing from HTA with the exception of the following actions:

- Printing all or part of the document in its original format.
- Transferring the document on websites or other electronic systems.
- Copying contents without any modification and stating HTA as copyright owner.

HTA reserves the right to make modifications or improvements to the relative documentation without prior notice.

This document is for information only. The manufacturer accepts no liability for errors contained herein or for incidental or consequential damages with the furnishing, performance, or use of this material.

Requests for authorization, additional copies of this manual or technical information on the latter, must be addressed to:

HTA s.r.l.
Via del Mella, 21
25131 Brescia (BS)
Italy
enquiry@hta-it.com
www.hta-it.com
+39 030 3582920

Trademarks and registered marks

In this document, registered marks and trademarks are not marked with ® and ™ symbols.

CONTENTS

1. Information on this document	7		
1.1 Document updates	8		
1.2 Provided documentation	10		
1.3 Message concepts	11		
2. Safety	12		
2.1 Warnings	13		
3. Reception	17		
3.1 Unpacking	18		
3.2 Identification	26		
4. Product description	28		
4.1 Introduction to the autosampler	29		
4.2 How to operate the autosampler	46		
4.3 Basic data about syringes	51		
4.4 Structures	52		
4.5 Oven location	57		
5. Installation	58		
5.1 Before proceeding	59		
5.2 Autosampler installation	60		
5.3 Sample rack installation	64		
5.4 Install the headspace syringe	65		
5.5 Install the tray	68		
5.6 Select the analyzer	69		
5.7 Align the injectors	70		
5.8 Set the position of the vials - The Touch and Plunger Zero procedure	72		
5.9 Set the communication with the analyzer	73		
5.10 Test the communication with the analyzer	74		
5.11 Set the autosampler in Headspace mode	75		
5.12 HTA Autosampler Manager installation	76		
5.13 Change to Liquid mode	78		
5.14 Set the autosampler in Liquid mode	88		
5.15 Enable the SPME feature for standard fibers	89		
5.16 Enable the SPME feature for 1.1 / 1.5 fibers	90		
5.17 Change to SPME mode	91		
5.18 Set the autosampler in SPME mode	108		
5.19 Accessories installation	109		
5.20 Re-installation or setup after a long downtime in Headspace mode	110		
5.21 Re-installation or setup after a long downtime in Liquid mode	111		
5.22 Re-installation or setup after a long downtime in SPME mode	112		
5.23 Enable the temporary feature	113		
6. Use in Headspace mode	114		
6.1 Get into Headspace mode	115		
6.2 Prepare the sample vials	122		
6.3 Operating Modes	123		
6.4 Set the method and sequence	126		
6.5 Start the Run	127		
6.6 Run management	128		
6.7 Sequences management	129		
6.8 Racks management	130		
6.9 Syringe management	131		
6.10 Syringe kit management	142		
6.11 Obtain support in case of autosampler problems	143		
6.12 Change the analyzer to be used	144		
6.13 Realign injectors	145		
6.14 Reset the position of the vial	147		
6.15 Get out of Headspace mode	148		
7. Use in Liquid mode	152		
7.1 Get into Liquid mode	153		
7.2 Prepare the solvent vials	160		
7.3 Prepare the sample vials	162		
7.4 Operating Modes	163		
7.5 Set the method and sequence	168		
7.6 Start the Run	169		
7.7 Run management	170		
7.8 Sequences management	171		
7.9 Racks management	172		
7.10 Syringe management	173		
7.11 Obtain support in case of autosampler problems	179		
7.12 Change the analyzer to be used	180		
7.13 Realign injectors	181		
7.14 Reset the position of the vial	183		

CONTENTS

7.15 Get out of Liquid mode	184	10.9 The Method page	298
7.16 Get out of Liquid mode (sealed tip)	184	10.10 The Method menu	300
8. Use in SPME mode	187	11. Screen description in SPME mode	303
8.1 Get into SPME mode (fiber)	188	11.1 Screen elements	304
8.2 Get into SPME mode (sealed tip)	194	11.2 Main menus	305
8.3 Prepare the sample vials	201	11.3 The Run menu (main menu)	308
8.4 Operating Modes	202	11.4 The Settings menu	311
8.5 Set the method and sequence	205	11.5 The Setup menu	315
8.6 Start the Run	206	11.6 The Sequences page	320
8.7 Run management	207	11.7 The Configuration menu	322
8.8 Sequences management	208	11.8 The Run page	327
8.9 Racks management	209	11.9 The Method page	328
8.10 Fiber management	210	11.10 The Method menu	330
8.11 Sealed tip management	220	12. Maintenance	333
8.12 Obtain support in case of autosampler problems	232	12.1 Preventive maintenance	334
8.13 Change the analyzer to be used	233	12.2 Maintain the autosampler	335
8.14 Realign injectors	234	12.3 Clean the autosampler	336
8.15 Reset the position of the vial	236	12.4 Autosampler packing	337
8.16 Get out of SPME mode (fiber)	237	13. Troubleshooting and errors in Headspace mode	341
8.17 Get out of SPME mode (sealed tip)	239	13.1 Troubleshooting	342
9. Screen description in Headspace mode	242	13.2 Start up errors	346
9.1 Screen elements	243	13.3 Run issues / errors	348
9.2 Main menus	244	13.4 System Integrity test errors	351
9.3 The Run menu (main menu)	247	13.5 Mismatch errors	352
9.4 The Settings menu	250	13.6 Preventive maintenance warnings	353
9.5 The Setup menu	254	14. Troubleshooting and errors in Liquid mode	355
9.6 The Sequences page	259	14.1 Troubleshooting	356
9.7 The Configuration menu	261	14.2 Start up errors	362
9.8 The Run page	266	14.3 Run issues / errors	364
9.9 The Method page	267	14.4 Mismatch errors	367
9.10 The Method menu	269	14.5 Preventive maintenance warnings	369
10. Screen description in Liquid mode	272	15. Troubleshooting and errors in SPME mode	371
10.1 Screen elements	273	15.1 Troubleshooting	372
10.2 Main menus	274	15.2 Start up errors	376
10.3 The Run menu (main menu)	277	15.3 Run issues / errors	378
10.4 The Settings menu	280	15.4 Mismatch errors	381
10.5 The Setup menu	284		
10.6 The Sequences page	288		
10.7 The Configuration menu	290		
10.8 The Run page	296		

CONTENTS

15.5 Preventive maintenance warnings	382
16. Manage unsolved errors	383
16.1 Deal with unsolved errors	384
16.2 Send log files to your Customer Representative	385
16.3 Tray emergency release	386
17. Software related errors	387
17.1 HTA Monitor related errors	388
17.2 HTA Autosampler Manager related errors	389
18. Technical features	390
18.1 Technical data	391
18.2 Electrical connections	399
18.3 Ordering information	402
18.4 Disposal	406
19. Conformity	407
19.1 Declarations	408
20. Warranty and service	410
20.1 Warranty	411
20.2 Service	412
21. Glossary of terms	413

1. Information on this document

1.1 Document updates

Version	Publication date	Updates
-	09-2023	First publication
a	11-2023	<ul style="list-style-type: none"> • Correct some typos • Update PC requirements for HTA Monitor and HTA Autosampler Manager
b	12-2023	<ul style="list-style-type: none"> • Update HTA Monitor installation procedure • Revised consumables ordering information
c	04-2024	<ul style="list-style-type: none"> • Add analyzers on which the CDS is applicable through HTA software integration and HTA software connectors • Add purge line specifications and external pressure regulator technical data • Add Carry over issue in Analyzer reproducibility issue chapter • Add green OR ring to identify Vial locator for Liquid technique • Add Needles in Ordering information • Change plunger with PTFE cover indicators in Ordering information (from 'x/-' to 'Yes/No') • Update PC requirements for HTA Autosampler Manager: add account requirement • Update Syringes in Ordering information
d	01-2025	<ul style="list-style-type: none"> • Add information in How to choose desorption depth • Add information in How to choose extraction depth • Add flowchart diagrams for Screen description • Add a Preventive maintenance warning • Add Usage in Troubleshooting chapter • Add two Vial detection problems in HS and SPME Mode • Add instructions on which vial locator use in each mode • Add some internal links to chapters • Correct Change to Liquid/SPME mode procedure • Correct Liquid kit contents • Correct Memory, Touch sensor and Mismatch errors messages • Correct Tray errors in Generic start-up errors • Correct some wrong data • Revise to Get into chapters • Update HTA Monitor installation procedure • Update HTA Autosampler Manager installation procedure and license purchasing code • Update PC requirements for HTA Autosampler Manager: remove Windows XP and update .NET Framework release • Revised the order of chapters in Ordering information
e	01-2025	<ul style="list-style-type: none"> • New HTA address • Update software connectors
f	03-2025	Add Sealed tip management in SPME mode
g	03-2025	<ul style="list-style-type: none"> • Remove references to 1.1 fibers / 1.5 fibers in Use in SPME

Version	Publication date	Updates
		mode
		<ul style="list-style-type: none">• Add SPME kit in Accessories
h	04-2025	<ul style="list-style-type: none">• Correct typos• Rearrange Consumables for Liquid Technique list

1.2 Provided documentation

Manual	Intended users	Code	Date	Distribution format
2800T Instruction Manual	Installer, Final Users	1.97.083	02/04/2025	Electronic

1.3 Message concepts

1.3.1 Safety messages

Warnings related to safety of the user and equipment as envisaged in this document are as follows:



CAUTION Indicates a dangerous situation which, if not avoided, may cause slight injury.

NOTICE: *Indicates obligations that if not observed may cause damage to the equipment.*

1.3.2 Other messages

Note: *Neutral and positive information that emphasize or add information to the main text. They provide information that can only be applied in special cases.*



Indicates the path in the software to access the page described.

2. Safety

2.1 Warnings

2.1.1 Personal protective equipment (PPEs)

During procedures wear the PPEs as prescribed by your laboratory regulations.

2.1.2 Electrical risks

- Live part hazard. Do not touch live parts while the autosampler is switched on (for example, connection power, power supply).
- Do not pour water or any other liquids over the power supply.
- Do not use water, solvents, detergents or spray on electrical parts.
- Ensure the power supply and power cord are not damaged.
- The supplied power cord must be inserted into a power outlet with a protective earth contact (ground).
When using an extension cord, make sure that the cord also has an earth contact.
- Use only the power cord and power supply provided with the autosampler and connect the power supply to the indicated voltage source. See ['Technical data' on page 391](#).
- Do not turn the autosampler on if you suspect that it has incurred any kind of electrical damage or if it shows visible signs of damage. If, for any reason, the safety protections were compromised or if the autosampler shows visible signs of damage, switch off the autosampler, disconnect the power cord and contact your Customer Representative.
- Switch off the autosampler and disconnect the power cord from the rear before cleaning or maintaining it.
- Do not use solvents nor detergents on the touch screen and all fragile parts of the autosampler.

2.1.3 Mechanical risks

- Finger crushing hazard. Do not touch nor place hands near the tray, turret, sledge, needle, oven cover, fan and plunger while they are moving.
- Handle the autosampler and all its consumables and spare parts with care, protecting them from mechanical damage.
- Do not touch nor adjust the injector adapter during injector alignment.

2.1.4 Hazardous substance risks

- The operator must be trained to operate in accordance with the safety standards for the laboratory. The autosampler does not have a specific protection mechanism against corrosive substances, or potentially infectious, toxic, radioactive or other hazardous substances.
- If the autosampler is used with hazardous samples the operator must use individual protection measures such as gloves, masks and goggles where appropriate. The use of such devices must comply with the regulations currently in force regarding the safety of testing laboratories.
- To avoid accidents, observe the rules of good laboratory practice while handling solvents, reagents and consumables. Observe the safety regulations about the used chemicals, as indicated in their specific Material safety data sheets (MSDS). Check carefully the physical and chemical properties of the substances before use.

2.1.5 Thermal risks

- Do not touch oven or what it is inside the oven during operations. The oven surfaces reach temperature of 170°C that cause burns.

- Do not touch, if present, the cleaning device during operations. The cleaning device reaches temperature of 300°C that cause burns.
- Do not touch vials without adequately protected hands for at least 30 minutes after the vial is removed from the oven. The vial reaches temperature of 170°C that cause burns.
- Do not touch syringe:
 - during operations.
 - without adequately protected hands for at least 30 minutes after the end of the operations.The syringe surfaces reach temperature of 150°C that cause burns.
- Do not obstruct the oven fan. The oven surfaces could reach higher temperature that cause burns.

2.1.6 Other risks

- The use of the autosampler for any technique invasive for humans or animals is explicitly prohibited. The autosampler must be used only in the manner described in this documentation. Any other use may affect the autosampler and the operator safety.
- Do not store or use the autosampler outside the indicated temperature and humidity range or outdoors.
- The autosampler is heavy. Lift it carefully and using both hands. Do not lift it from the sliding lid.
- The mounting plate has sharp metal parts. Do not open the locking hook by hand.
- The autosampler tray and syringe locations have sharp metal parts. Do not touch them.
- To avoid damaging the equipment and compromising its safety, only use the interface cable, power cord and power supply provided with the autosampler.
- Do not move the turret by hand.
- Maintain a distance of at least 10 cm (3.9") from the autosampler.
- Only use permitted cleaning products.
- If the barcode reader is enabled, do not look at the LED light within a distance of 1 cm. This could cause eye injury.
- Do not use flammable gas like hydrogen to supply the autosampler.







2.1.7 Disclaimer





- The operator is responsible for damages caused, even only in part, by improper use of the autosampler, by its unauthorized modification, by a different or missing execution of the procedures described in this manual, or by the use of the autosampler violating the applicable laws and regulations.
- The insertion and positioning of sample, sample vials and syringes and the selection of reagents and washing solvents and purge gas are the sole responsibility of the operator. The operator must operate in accordance with the law, in particular the area concerning the safety of testing laboratories and the accident prevention.
- The standard safety measures and precautions regarding use, maintenance and repair, are not intended for applications that require the use of hazardous substances. For these and other applications requiring IVD certification, extra safety measures are necessary. These additional measures are the sole responsibility of the operator and of the supplier of the IVD system in which the autosampler will be integrated.
- The autosampler operator must be trained about the hazards, use and proper disposal of any hazardous substance. This can be found in the Material safety data sheets (MSDS) and in the regulations, laws, directives and safety decrees currently in force for testing laboratories.

- It is responsibility of the operator the handling and use of hazardous substances in accordance with their safety data sheets (MSDS) and with the law about safety in laboratories. Do not use sample vials without sealing caps or closures. Vapors may be hazardous. Acidic vapors may be corrosive for the mechanical parts.
- It is responsibility of the operator to decontaminate the autosampler from any residual hazardous substances. The cleaning and decontamination agents shall not cause hazard due to reaction with the autosampler or other materials contained inside. We recommend to use a solution of 70% isopropanol (isopropyl alcohol) and 30% water as cleaning solution for external instrument surfaces. If there are doubt about the use of cleaning or decontamination agents contact your Customer representative.
- The autosampler can heat the vials in the oven and the sample in the syringe up to 170°C. Verify, according to an analysis of the laboratory's safety standards and safety cards of chemical substances, that substances can be used at this temperature without any risk to the operator and/or the environment.

2.1.8 Safety labels on the instrument

Here are the labels that warn about how to avoid risks:

Symbol	Description	Location
	Electrical risk	
	Generic risk. Risk for moving parts.	
	Finger crushing hazard. Do not touch nor place hands near the tray, turret, sledge, needle and plunger while they are moving.	 

Symbol	Description	Location
	Hot surface. Do not touch oven, syringe and vials during operations.	
	Risk of electrocution. The power supply shall be kept away from liquid sources (liquid containers, liquid spillages and so on). Do not open the power supply.	

3. Reception

3.1 Unpacking


3.1.1 Storage warnings

NOTICE: The package must be stored indoor to avoid exposure to adverse weather conditions. After unpacking, store the packing materials for future use.

3.1.2 Precautions on reception

On reception, check the integrity of the package. If the package is damaged or tampered with, contact your Customer representative. After unpacking, check the presence of all components and their integrity.


3.1.3 Contents



Part	Name	Function
A	2800T All-in-one GC Autosampler (hereinafter "autosampler")	Automate the sampling process.
B	Power supply	Allow to power the autosampler.
C	Power cord	Allow to power the autosampler.
D	Ethernet cable	Connect the autosampler to the network.
E	Syringe pointer	Allow to set the autosampler in operating position, to open the locker and to realign the vial locator to the needle tip.
F	Allen key 2 mm	Allow to change the vial locator.
G	Documentation download link (2x)	<ul style="list-style-type: none"> 1x contain the digital code to download the manual; 1x contain the digital code to download the software.


Part	Name	Function
H	Product paper	List the package contents and provide declarations of conformity and warranty.
I	Activation key for HTA Autosampler Manager	Allow to use the trial license of the software.
J	Safety and conformity information	Provide safety and conformity information.
K	QR code label for installation	Obtain technical assistance on the autosampler installation.
L	Warning label	Warn to read the manual.

Headspace kit:




Part	Name	Function
A	Sample standard rack	Hold 42 sample vials of 20 ml / 10 ml / 6 ml.
B	Vial locator	Allow to transport the vial and guide the needle. For details, see 'Vial locators' on page 55 .
C	Heated syringe assembly	Heat and contain syringe glass and syringe plunger.
D	Allen key 2.5 mm	Tool for syringe replacement.
E	Sample vials 20 ml	Contain the samples.
F	Caps with septa	Close the sample vials and allow perforation.

Liquid kit:



Part	Name	Function
A	Sample standard rack	Hold 121 sample vials.
B	Vial locator	Guide the needle. For details, see 'Vial locators' on page 55 .
C	Syringe holder	Hold and move the syringe.
D	Plunger extension	Hold and move the syringe plunger.
E	Syringe	Perform the necessary procedures on the sample.
F	Tweezers	Remove the spacers from the oven.
G	Sample vials 1,5 ml	Contain the samples.
H	Vials 10 ml	Contain the solvents / Used as waste vials.
I	Spacers	Allow to use 10 ml vials as solvent vials in the oven.
J	Low-maintenance caps with septa	Close the solvent vials and allow perforation.

3.1.4 Optional contents
SPME kit for standard fibers:



Part	Name	Function
A	Vial locator	Allow to transport the vial and guide the needle. For details, see 'Vial locators' on page 55 .


Part	Name	Function
B	Fiber holder	Hold and move the fiber.
C	Sample vials 20 ml	Contain the samples.
D	Caps with septa	Close the sample vials and allow perforation.
E	Product paper	List the package contents and provide declarations of conformity and warranty.
F	Activation code for SPME mode	Provide the feature code to enable the SPME mode.

Note: Two kinds of kit are available:

- SPME kit - standard fibers (w/ fibers);
- SPME kit - standard fibers (w/o fibers).

In case the purchased kit is 'w/ fiber', the kit contains 3 fibers too (PDMS 10 mm, needle size 23 ga).


SPME kit for 1.1 fibers / 1.5 fibers:



Part	Name	Function
A	Vial locator	Allow to transport the vial and guide the needle. For details, see 'Vial locators' on page 55 .
B	Sealed tip holder	Hold the sealed tip.
C	Plunger extension	Move the sealed tip.
D	Sealed tip spacer	Allow to mount the sealed tip on the instrument.
E	Activation code for sealed tip fiber	Provide the feature code to enable the usage of sealed tip fiber.


Note: To use the SPME kit - 1.1 fibers / 1.5 fibers it is necessary to have the SPME fiber - standard fibers (w/ fibers) or the SPME kit - standard fibers (w/o fibers).

HTA Autosampler Manager (Standard license):



Part	Name	Description
A	HTA Autosampler Manager (Standard license)	Allow to use the Standard license of the HTA Autosampler Manager software.


Upgrade to CFR 21 Part 11 for HTA Autosampler Manager:




Part	Name	Description
A	Upgrade to CFR 21 Part 11 for HTA Autosampler Manager	Allow to use the HTA Autosampler Manager software in CFR 21 Part 11 version. It requires 1.90.499 HTA Autosampler Manager (Standard license) to work.

3.1.5 Unpack the autosampler


- 1. Open the box and remove the free package cushioning.
- 2. Hold the autosampler with two hands and extract it from the box.
- 3. Place the autosampler vertically on a flat surface.
- 4. Cut the adhesive tape on the top to open the packaging.
- 5. Remove the upper package cushioning.




6. Remove the accessory boxes and extract the included accessories.




7. Remove the lower package cushioning.



8. Check the product paper to verify if the position of oven is consistent with your analyzer specification. If not, please contact your Customer Representative.
9. Remove the protection belt and the polythene bag from the turret.



10. Remove the turret protection.
NOTICE: Do not turn the turret.




3.1.6 Lift the autosampler

1. Remove all the removable parts, such as sample vials and sample rack.
2. Hold the autosampler base without touching the frontal part and the oven structure.

Note: Parts to be touched are highlighted in green; parts not to be touched are highlighted in red.



CAUTION Heavy product. Lift the autosampler carefully and using both hands.








3.2 Identification

3.2.1 Manufacturer's contacts






HTA s.r.l.
Via del Mella, 21
25131 Brescia (BS)
Italy
enquiry@hta-it.com
www.hta-it.com
+39 030 3582920

3.2.2 Autosampler identification label

HTA s.r.l.
Via del Mella 77/79
25131 Brescia ITALY
www.hta-it.com
Made in Italy



Model: 2800T
P.N.: 2800T
S.N.: 213052

Part	Description
	CE mark, declaration of conformity for Europe
	UK mark, declaration of conformity for the United Kingdom
	China RoHS
	CMIM mark, declaration of conformity for Morocco
	Dispose of the instrument as electronic waste (European Community Directive 2012/19/EU).

Part	Description
Model	Product model
P.N.	Part number
S.N.	Serial number

4. Product description

4.1 Introduction to the autosampler

4.1.1 Description

The autosampler is an automatic sampling device that automates the procedures for injecting liquid samples and headspace sample gas in GC or other chemical analyzers. By getting the optional kits, the autosampler can inject a sample trap into a SPME fiber (standard fibers or 1.1 mm / 1.5 mm fibers).

4.1.2 Intended users

The autosampler is intended to be used by operators with a background in chemistry and knowledge of standard laboratory procedures and instrumentation.

4.1.3 Intended use

This autosampler is designed to be used in chemical analysis laboratories in R&D, contract or control laboratories.

4.1.4 Switch between modalities

The autosampler allows for quick and easy switching from one injection technique to another: a procedure will show which parts must be changed in order to begin to inject in the desired injection technique.

4.1.5 Typical operating principle for Headspace technique

For correctly serving the analyzer, the autosampler must be mounted on the analyzer specific mounting kit. Once installed, prepared and set for running the desired operations on the samples through specific methods and sequences, the autosampler and the analyzer start to communicate autonomously through the interface cable. The autosampler performs the following actions for each sample vial:

Phase	Description
1	Load the vial into the oven.
2	Incubate and shake the vial for a defined time.
3	When the analyzer is ready, aspirate the headspace sample gas.
4	Inject it in the analyzer injector/s.
5	Signal to the analyzer that the injection has been done.
6	If necessary, clean the syringe with the flush gas.

Note: In case of electrostatic discharge, the autosampler can restart automatically without malfunctioning or losing information.

4.1.6 Sample preparation overlap for Heaspace and SPME technique

The autosampler is able to overlap sample preparation in order to improve the high throughput. In case the combination of analysis time, conditioning time and preparation time will allow it, the autosampler will take automatically care to load the subsequent sample vial in a free position of the oven, to start the conditioning of the subsequent vial during the conditioning of the previous sample.

4.1.7 Methods for Headspace technique

A method instructs the autosampler on what to do for each sample vial. The autosampler has ten preset methods. A method is defined by several parameters, whose values can be modified according

to the specific analysis required.

The following table shows the preset parameter values for the ten methods:

List item	Set parameter
General parameters	
Syringe Volume	2.5 ml
Analysis Time	<i>GC Cycle Time: GC run time plus any additional time needed for the GC to execute any post-run program, then return to a Ready state.</i>
Advanced parameters	
Enrichment Cycle(s)	1
Dwell Between Injections	0 min
Pre-Flush	Off
Preparation parameters	
Syringe Temperature	Off
Oven Temperature	Off
Incubation	The incubation time has to be established with method development. For details see 'Typical operating principle' on page 29 .
Shaker On	0.2
Shaker Off	0.1
Sample parameters	
Sample Volume	0.5
Fill Volume	1 ml
Pull Up Strokes	2
Equilibration delay	1 sec
DHS time	<i>To be set for special application only.</i>
Sample speed	6 ml/min
Syringe prefill	No
Injection parameters	
Injection Speed	30 ml/sec
Pre Injection Dwell	1 sec
Post Injection Dwell	3 sec
Flush Time	3 min

4.1.8 How to calculate the Analysis time correctly

The Analysis time is the GC run time plus any additional time needed for the GC to execute any post-run program and then return to a Ready state.

Please consider that the autosampler uses the Analysis time (GC cycle time) to calculate throughput and timing. An accurate Analysis time is crucial to optimizing throughput and for correctly processing samples.

If the Analysis time set into autosampler method is...	this can cause...
longer than the actual one	lowered throughput because vials wait longer than needed before processing.
shorter than the actual one	sequence faults because the vial may sit too long while waiting for the GC to become Ready.

It is better to enter a longer time than needed than to enter a too short time and possibly reduce sample quality.

To determine the Analysis time, program the GC to perform a sequence of two or three blank (no injection) runs. If you are using a Chromatography Data system, you can determine the cycle time from the data system sequence log and compare the start times for the runs. A good Analysis time is the average difference between the start times, plus 0.2 to 0.5 minutes. If you are not using a Data system, watch the GC. Count the time between the start of the first run and when the GC becomes Ready for the second run, then add 1 minute.

You can also estimate the Analysis time without making a run. By adding the GC oven program duration and the duration of any post-run programs, you can get close to the true cycle time.

Also consider time for data processing. While in most cases data processing is not a problem, a very busy Data system may need extra time between samples.

To validate the Analysis time, look at the GC run time. The Analysis time cannot be less than the total run time. Then look at the GC post run program duration. The Analysis time cannot be less than the total run time plus the post run program duration.

Finally, consider that the GC cooling time may vary according to room temperature: the higher the room temperature the more time is needed for the GC cooling time; typically, it implies a shorter time during winter seasons and a longer time during summer seasons. Therefore, the analysis time should be set according to the worst scenario (summer season).

4.1.9 Principle of “Optimized Sequence” to achieve run time and processing optimizations of sequence steps

During a creation of a sequence it is strongly recommended to program the sequence step with accuracy. 2800T autosampler automatically creates optimized sequences in order to optimize the executions and reduce the total run time of the sequence.

This table shows the difference between a non-optimized sequence and an optimized sequence: during the optimized sequence the autosampler takes care to load more than one vial in oven (if possible, depending on the autosampler method parameters and on GC run time) to get the best possible throughput.

No.	Non-optimized	No.	Optimized
1	Load vial 1	1	Load vial 1
2	Incubate vial 1	2	Incubate vial 1
3	Inject vial 1	3	During incubation of vial 1, load vial 2
4	Unload vial 1	4	Incubate vial 1 and 2
5	Wait for analysis time	5	Inject vial 1
6	Load vial 2	6	Unload vial 1
7	Incubate vial 2	7	Load vial 3
8	Inject vial 2	8	Inject vial 2

No.	Non-optimized	No.	Optimized
9	Unload vial 2	9	...
10	Wait for analysis time	10	
11	Load vial 3	11	
12	Incubate vial 3	12	
13	Inject vial 3	13	
14	Unload vial 3	14	
15	Wait for analysis time	15	
	...		

To make the optimizations there are the easy rules to follow:

- The samples to be executed with the same method must be placed in consecutive position;
- Each sequence step must contain all the samples to be executed with the same method and with the same injection mode.

In the following table there are some example of step programming: to simplify the example they refers just to the method, but the same example can be used for the injection mode parameter.

Example	Details	Description
Optimized sequence	Step 1: <ul style="list-style-type: none"> • Method: 1 • First sample: A1 • Last sample: A3 • Injection mode: front ... Step 2: <ul style="list-style-type: none"> • Method: 2 • First sample: A4 • Last sample: A5 • Injection mode: front 	The sample vials to be executed with the same method and injection mode are placed in consecutive position in the rack, and every group of sample vials is contained in one step, to reduce the required execution time.
Not-optimized sequence: one sample per step	Step 1: <ul style="list-style-type: none"> • Method: 1 • First sample: A1 • Last sample: A1 ... Step 2: <ul style="list-style-type: none"> • Method: 1 • First sample: A2 • Last sample: A2 	The autosampler will not optimize the run-time, because the consecutive sample vials to be executed with the same method were not programmed in the same step.
Not-optimized sequence: sample vials not placed in	Step 1: <ul style="list-style-type: none"> • Method: 1 • First sample: A1 	The autosampler will not optimize the run-time at its best, because the sample vials placed in position A5 and A6 could be placed in position A4 and A5, to be programmed in the step 1 with the rest of the sample

Example	Details	Description
order, taking care of the method to be used	<ul style="list-style-type: none"> Last sample: A3 ... Step 2: <ul style="list-style-type: none"> Method: 2 First sample: A4 Last sample: A4 ... Step 3: <ul style="list-style-type: none"> Method: 1 First sample: A5 Last sample: A6 	vials that must be executed with the same method.
Not-optimized sequence: sample vials not in consecutive positions	Step 1: <ul style="list-style-type: none"> Method: 1 First sample: A1 Last sample: A3 ... Step 2: <ul style="list-style-type: none"> Method: 1 First sample: A5 Last sample: A6 	The autosampler will not optimize the run-time at its best, because the sample vials are not consecutive. They could have been placed in consecutive position (from A1 to A5) and be contained in just one step.

4.1.10 Typical operating principle for Liquid technique

For correctly serving the analyzer, the autosampler must be mounted on the analyzer specific mounting kit. Once installed, prepared and set for running the desired operations on the samples through specific methods and sequences, the autosampler and the analyzer start to communicate autonomously through the interface cable. The autosampler waits until the analyzer is ready to receive the samples and then automatically moves the syringe installed in its turret to perform the following actions for each sample vial:

Phase	Description
1	Wash the syringe with the desired solvent/s.
2	Prime syringe by aspirating the sample from the sample vial and emptying the syringe in the waste.
3	Eliminate the bubbles from the sample in the sample vial.
4	Aspirate the sample and inject it in the analyzer injector/s.
5	Signal to the analyzer that the injection has been done.
6	If necessary, wash the syringe with the appropriate solvent/s.
7	Wait until the analyzer is ready to receive another sample and repeat all the phases for another sample vial.

Note: In case of electrostatic discharge, the autosampler can restart automatically without malfunctioning or losing information.

4.1.11 Sample preparation overlap for Liquid technique

The autosampler is able to overlap sample preparation, in order to improve the high throughput in case of usage of "pre-ahead" sample list (available through HTA Autosampler Manager software).

4.1.12 Methods for Liquid technique

A method instructs the autosampler on what to do for each sample vial. The autosampler has ten preset methods. A method is defined by several parameters, whose values can be modified according to the specific analysis required.

The following table shows the preset parameter values for the ten methods:

List item	Set parameter
General parameters	
Syringe Volume	10 µl
Analysis Time	<i>Time of the sample analysis (time from the injection to the moment in which GC will be ready again).</i>
Fill Speed	1 µl/sec
Viscosity Delay	5 sec
Solvent wash parameters	
Pre Wash Cycle (s)	3
Pre Wash Volume	10 µl
Pre Wash Mode	Inject
Post Wash Cycle (s)	3
Post Wash Volume	10 µl
Post Wash Mode	Inject
Internal standard parameters	
Enable	No
Volume	-
Air Gap Mode	-
Air Gap Volume	-
Sample parameters	
Sample Volume	1 µl
Fill Volume	5 µl
Air Volume	1 µl
Sample Washes	2
Pull Up Strokes	5
Injection parameters	
Injection Speed	100 µl/sec
Pre Injection Dwell	1 sec
Post Injection Dwell	3 sec


4.1.13 Pre and post-washing operations

The autosampler performs specific pre- washing operations on the syringe before injecting the sample in the analyzer. After having injected the sample in the analyzer, it performs post- washing operations.

Both washing operations consist of aspirating a solvent with the syringe and then injecting it in the waste vial. When a combination of solvents is indicated (for example, A + B), the second solvent aspiration occurs only after the disposal of the first solvent in the waste vial.

4.1.14 Use of an internal standard

The autosampler allows the use of the internal standard technique to inject in the analyzer a sample with a reference standard at the same time. In such case, the contents of the syringe will be organized as shown below:



4.1.15 Typical operating principle for SPME technique

For correctly serving the analyzer, the autosampler must be mounted on the analyzer specific mounting kit. Once installed, prepared and set for running the desired operations on the samples through specific methods and sequences, the autosampler and the analyzer start to communicate autonomously through the interface cable. The autosampler performs the following actions for each sample vial:

Phase	Description
1	Load the vial into the oven.
2	Incubate and shake the vial for a defined time.
3	Expose the fiber into the liquid sample or the headspace sample gas.
4	If necessary, expose the fiber into a wash vial (pre-wash procedure).
5	Expose the fiber in the analyzer injector and signal to the analyzer that the exposure has begun.
6	If necessary, expose the fiber into a different injector or the cleaning device (post-clean procedure).

Note: In case of electrostatic discharge, the autosampler can restart automatically without malfunctioning or losing information.

4.1.16 How to choose the desorption depth correctly

Desorption depth is the distance from vial locator to the fiber tip. For different depths set, the autosampler exposed:

- only fiber;
- fiber and needle;
- fiber, needle and fiber holder

Note: the desorption depth is set during the alignment of the injectors position (see ['Realign injectors' on page 234](#)). During the run, the instrument will always expose completely the fiber, therefore with 'desorption depth' we mean the sum of the fiber length (10 mm or 20 mm) and the needle fiber penetration depth.

Example 1 - desorption depth: 40 mm (w/ 10 mm fiber) -> 30 mm of needle movement + 10 mm of fiber exposure

Example 2 - desorption depth: 60 mm (w/ 10 mm fiber) -> 49 mm of needle movement (maximum value for the needle) + 11 mm of fiber exposure (10 mm fiber + 1 mm fiber extension)

Note: Be aware that the fiber exposure is 28 mm:

- in case of 10 mm fiber: 10 mm of fiber + 18 mm of fiber extension

- in case of 20 mm fiber: 20 mm of fiber + 8 mm of fiber extension

Therefore, it is highly important to calculate the proper desired depth before to execute the injector alignment, to be sure that the needle penetration depth will be enough to allow the fiber to be exposed below the analyzer septum.

During this alignment procedure, the instrument could move in a different way than when it will reach the desired depth into the run, if the desired depth is lower than 59 mm.






This behaviour is due to the fact that, during the alignment procedure, the instrument does not know which will be the depth chosen, and therefore it will move first the needle fiber for its maximum (49 mm) and then, it will start exposing the fiber.


Example of different behaviour with a 10 mm fiber

Desorption depth	Behaviour during injector alignment procedure	Behaviour during the run
40 mm	The instrument will move the needle fiber for 40 mm.	The instrument will move the needle fiber for 30 mm, then the fiber will be exposed for 10 mm.
60 mm	The instrument will move the needle fiber for 49 mm, then it will expose the fiber for 11 mm.	As for during the injector alignment.




Use this info to set a proper depth, that will assure that the fiber is exposed below the injector septum.

Injection depth setting in Standard SPME mode

Fiber length	Injection depth (mm)	Instrument behaviour	Notes
10 mm	<10	The instrument moves the needle 1 mm under the vial locator, then expose fiber moving only the plunger.	
	10<depth<59	The instrument moves the needle 10 mm over the final position, then moves the plunger to expose 10 mm of fiber.	
	59<depth<77	The instrument moves the needle until 49 mm, moves the plunger to expose 10 mm of fiber and then from 0 mm to 18 mm of fiber holder to reach the right depth.	
20 mm	<20	The instrument moves the needle 1 mm under the vial locator, then expose fiber moving only the plunger.	
	20<depth<69	The instrument moves the needle 20 mm over the final position, then moves the plunger to expose 20 mm of fiber.	

Fiber length	Injection depth (mm)	Instrument behaviour	Notes
	69<depth<77	The instrument moves the needle until 49 mm, moves the plunger to expose 20 mm of fiber and then from 0 mm to 8 mm of fiber holder to reach the right depth.	

Injection depth setting in Sealed Tip SPME mode





Fiber length	Injection depth (mm)	Instrument behaviour	Notes
20 mm	<20	The instrument moves the needle 1 mm under the vial locator, then expose fiber moving only the plunger.	
	20<depth<69	The instrument moves the needle 20 mm over the final position, then moves the plunger to expose 20 mm of fiber.	
	69<depth<77	The instrument moves the needle until 49 mm, moves the plunger to expose 20 mm of fiber and then from 0 mm to 8 mm of fiber holder to reach the right depth.	



4.1.17 How to choose the extraction depth correctly

For Extraction depth and Derivatization depth, the position “Under septum” (depth =0) is between 5 and 7 mm under vial septum. This spread depends on the vials height and the mechanical tolerance of the oven loading.



Therefore, during the extraction and derivatization phase the fiber is always exposed below the vial septum (see '[Sample vial specifications for SPME technique](#)' on page 395). Take care to choose the proper depth in order to have the fiber exposed into the liquid samples or to stay higher than the samples in case of SPME headspace sampling.


Extraction / Derivatization depth setting in Standard SPME mode

Fiber length	Injection depth (mm)	Instrument behaviour	Notes
10 mm	<10	The instrument moves the needle to have the needle tip just under the septum, then expose the fiber moving only the plunger.	 <p>Under septum, depth = 0</p>
	10 < depth < 28	The instrument moves the needle 28mm over the final position, then moves the plunger to expose 10 mm of fiber and if necessary, up to 18mm of fiber holder.	 <p>10mm Fiber</p>
	28 < depth < 47	The instrument moves the needle 28mm over the final position, then moves the plunger to expose 10 mm of fiber and 18mm of fiber holder and, if necessary, it moves the needle down up to 19mm further.	 <p>Under septum, depth = 0</p>
20 mm	<20	The instrument moves the needle to have the needle tip just under the septum, then expose the fiber moving only the plunger.	 <p>Under septum, depth = 0</p>

Fiber length	Injection depth (mm)	Instrument behaviour	Notes
	20<depth<28	The instrument moves the needle 28mm over the final position, then moves the plunger to expose 20 mm of fiber and if necessary, up to 8mm of fiber holder.	 <p>Under septum, depth = 0</p>
	28<depth<47	The instrument moves the needle 28mm over the final position, then moves the plunger to expose 20 mm of fiber and 8mm of fiber holder and, if necessary, it moves the needle down up to 19mm further.	 <p>Under septum, depth = 0</p>

Extraction / Derivatization depth setting in Sealed Tip SPME mode

Fiber length	Injection depth (mm)	Instrument behaviour	Notes
20 mm	<20	The instrument moves the needle to have the needle tip just under the septum, then expose the fiber moving only the plunger.	 <p>Under septum, depth = 0</p>
	20<depth<28	The instrument moves the needle 28mm over the final position, then moves the plunger to expose 20 mm of fiber and if necessary, up to 8mm of fiber holder.	 <p>Under septum, depth = 0</p>

Fiber length	Injection depth (mm)	Instrument behaviour	Notes
	28<depth<47	The instrument moves the needle 28mm over the final position, then moves the plunger to expose 20 mm of fiber and 8mm of fiber holder and, if necessary, it moves the needle down up to 19mm further.	

4.1.18 Methods for SPME technique

A method instructs the autosampler on what to do for each sample vial. The autosampler has ten preset methods. A method is defined by several parameters, whose values can be modified according to the specific analysis required.

The following table shows the preset parameter values for the ten methods:

List item	Set parameter
General parameters	
Vial Type	20 ml
Analysis Time [min]	5
Preparation Time	00:13:30
Wash & Clean parameters	
Pre-Washing	No
Pre-Washing Time [min]	0.0
Post-Cleaning	No
Post-Cleaning Time [min]	0.0
Extraction parameters	
Extraction Time [min]	10.0
Extraction Depth [mm]	25
Derivatization	None
Derivatization Time [min]	0.0
Derivatization Depth [mm]	3
Conditioning parameters	
Oven Temperature [°C]	Off
Incubation Time [min]	5
Shaker On [min]	0.2
Shaker Off [min]	0.1
Desorption parameters	
Desorption [min]	3.0

4.1.19 Sequences

Once the methods have been modified, it is necessary to set a sequence. The sequence represents the procedure that the autosampler must follow to process all the available sample vials. The sequence contains one or more steps. Each step defines, among other things, the position of the sample vials that must be processed with the same method.

Example:

The operator has ten sample vials: five sample vials must be injected in the analyzer with method A and the other five with method B.


Since the operator wants to use two different methods, he/she creates a sequence with two steps. The operator sets the first step so that the autosampler processes the sample vials from number 1 to number 5 with method A, and the second step so that the autosampler processes the sample vials from number 6 to number 10 with method B.

4.1.20 Software used

The autosampler is equipped with HTA Monitor and HTA Autosampler Manager software. Their functions are shown in the following table:

Software name	Is the software required?	Function
HTA Monitor	Yes, Required	<ul style="list-style-type: none"> • Allow to control the autosampler through a virtual screen (requires Ethernet connection). • Save the autosampler log records in a database for receiving technical assistance. • Generate a QR code with a basic log to be sent to the Customer Representative in case help is needed.
HTA Autosampler Manager	No, Optional	<ul style="list-style-type: none"> • Allow to control the autosampler by its dedicated interface • It offers CFR 21 Part 11 compliancy by buying additional license • Allow to backup autosampler data, restore them on critical events and print relevant information.

Depending on the operation mode desired, the autosampler can be managed directly by the analyzer software (also known as CDS – Chromatography Data System or Data Acquisition software). For this to be possible, specific software must be activated. Some software have activation keys that must be provided to the autosampler in order for the communication to be possible.

Read '[How to operate the autosampler](#)' on page  46 or contact your Customer Representative for more information.

Versions and licenses of HTA Autosampler Manager

HTA Autosampler Manager is available in the following versions. All the versions requires a license. All the licenses require an activation key to be activated.

Version	License Type	Description
Standard	Trial	<ul style="list-style-type: none"> • Time limitation: 60 days • Supplied for free with the autosampler

Version	License Type	Description
		<ul style="list-style-type: none"> Linked to the PC used to activate the license and the autosampler serial number
Standard	Full	<ul style="list-style-type: none"> No time limitation Linked to the autosampler serial number
CFR 21 Part 11	Full	<ul style="list-style-type: none"> No time limitation Linked to the PC used to activate the license Requires also a full standard license

4.1.21 System integrity test for Headspace technique

The System integrity kit is sold as an option. The System integrity test is used to verify the syringe integrity.

If this test is enabled, at the beginning of each batch the system performs a check to verify the syringe integrity through a heuristic procedure. This test is carried out using the system integrity tool. The septum of this tool should be replaced periodically.

This check can verify if it is better to replace the syringe (barrel and plunger). Please note that this check does not replace the necessary periodic validations you have to program for the instrument qualification. If the test fails, a warning message appears.

4.1.22 Vial leakage check for Headspace technique

This function can only be enabled by the HTA Autosampler Manager: **Setup > Options**. If this option is enabled, the pressure inside vials is monitored by a heuristic procedure in order to check against anomalous values that are indicative of a vial leakage problem.

The Vial leakage test can discriminate samples correctly sealed vs. bad crimping or missing septa. Each sample is marked with **PASS** or **FAIL** information in HTA Autosampler Manager Windows.

4.1.23 Method development option for Headspace and SPME technique


During method development it can be useful to test methods characterized by a progressive increment of a specific parameter in order to define which are the best settings for this parameter.

The parameters that can be increased in this way are the syringe/oven temperature and the incubation time.


This option is only available with HTA Autosampler Manager, using the **Method** development tab. For the detailed explanation of this option, please refer to the HTA Autosampler Manager online help.

4.1.24 Possible Ethernet configurations


Autosampler directly connected to the PC:



Autosampler connected to a LAN:



Autosampler connected to a wide LAN:



4.2 How to operate the autosampler

4.2.1 By autosampler touch screen

There are two ways to have access to the autosampler touch screen: by using the on-board touch screen display or by using the Virtual Screen, accessible by HTA Monitor.

Both ways allow to perform the same operation.

By autosampler touch screen: all samples are processed in the same way

It will just be necessary to load the samples on the autosampler and then press the **START** button on the CDS and on the autosampler touch screen. No further programming/action is required.

Additional requirement	None
Suggested if	The daily routine analysis is always constant (e.g. routine laboratories, manufacturing control laboratories, etc.) and what changes is just the number of samples to be injected.
Intended use	<p>When the user needs during daily routine are just to choose how many samples to inject, and therefore:</p> <ul style="list-style-type: none"> ○ The autosampler method is always the same; ○ The following autosampler sequence parameters for HS injections are always the same: <ul style="list-style-type: none"> ▪ Tray type ▪ Injection mode ○ The following autosampler sequence parameters for Liquid injections are always the same: <ul style="list-style-type: none"> ▪ Tray type ▪ Cycle(s) ▪ Prewash solvent ▪ Postwash solvent ▪ Injection mode ○ The following autosampler sequence parameters for SPME injections are always the same: <ul style="list-style-type: none"> ▪ Tray type ▪ Injection mode

By autosampler touch screen: you are required to process samples in different ways

Once the programming is done (when the parameters for the specific batch of injections have been set), the user must start the sample list on the CDS and on the autosampler touch screen.

Additional requirement	None
Suggested if	The user must change method and/or sequence more often, due to the variety of the analysis that must be done (e.g. research laboratories, demonstration laboratories, etc.).
Intended use	When the user must change often the autosampler main parameters (for example, to perform different injection with different injected volume, or executing washes from

different solvent vials, or use different incubation time or temperatures), the user can easily have access to the complete set of autosampler parameter to program it according to his needs.

4.2.2 By HTA Autosampler Manager

By HTA Autosampler Manager: Standard license

Once the programming is done (when the parameters for the specific batch of injections have been set), the user must start the sample list on the CDS and on the HTA Autosampler Manager software.

Additional requirement	Purchase 1.90.499 - HTA Autosampler Manager (Standard license)
Suggested if	<ul style="list-style-type: none"> ○ The user wants to back up files about autosampler functioning (e.g. method files, setup files, report, etc.) ○ the user wants to generate report about the autosampler operations ○ the user wants to use multi-language user interface ○ the user does not need to work in compliance to CFR 21 part 11
Intended use	When the user needs one of the above features, the HTA Autosampler Manager software is the best option to control the autosampler.

By HTA Autosampler Manager – CFR 21 part 11 version

Once the programming is done (when the parameters for the specific batch of injections have been set and the modifications have been tracked and confirmed, the user must start the sample list on the CDS and on the HTA Autosampler Manager software.

Additional requirement	Purchase 1.90.499 - HTA Autosampler Manager (Standard license) and 1.90.503 - Upgrade to CFR 21 Part 11 for HTA Autosampler Manager
Suggested if	<ul style="list-style-type: none"> ○ The user wants to back up files about autosampler functioning (e.g. method files, setup files, report, etc.) ○ the user wants to generate report about the autosampler operations ○ the user wants to use multi-language user interface ○ the user needs to work in compliance to CFR 21 part 11
Intended use	When the user wants to use the HTA Autosampler Manager features and to work according to CFR 21 part 11 compliance, the HTA Autosampler Manager with CFR 21 part 11 license is the best option to control the autosampler.

4.2.3 By CDS through HTA software integration

The CDS will be used to start the sample list on the analyzer and to start the same sample list on the autosampler.

Applicable on:	<ul style="list-style-type: none"> • Agilent <ul style="list-style-type: none"> ○ ChemStation GC revision B.04.03-SP1[87] or higher ○ MassHunter GC/MS revisions B.06 or higher ○ MSD ChemStation E.00, E.01, E.02 or higher ○ OpenLAB CDS ChemStation Edit. revision C ○ OpenLAB CDS - EzChrom Edition for GC • Shimadzu <ul style="list-style-type: none"> ○ GCSolution ○ GCMSsolution ○ LabSolutions LC/GC ○ LabSolutions LCMS ○ LCSolution ○ LCMSsolution • SSI <ul style="list-style-type: none"> ○ EzChrom Elite • Thermo <ul style="list-style-type: none"> ○ Chromquest
Additional requirement	Purchase 1.90.499 - HTA Autosampler Manager (Standard license)
Intended use	The method/sequence must be programmed on the autosampler by using HTA Autosampler Manager. This is typically done at autosampler installation or when a change in the method/sequence is required. Then, the sample list is started directly from the CDS

4.2.4 By CDS through HTA software connectors

Applicable on	<ul style="list-style-type: none"> • Agilent <ul style="list-style-type: none"> ○ OpenLAB CDS ChemStation Edit. revision C ○ OpenLAB CDS • Scion <ul style="list-style-type: none"> ○ Compass CDS <p>The list of the CDS which are supported by advanced integration is continuously increasing: contact your Customer Representative to check if your CDS has been recently supported.</p>
Additional requirement	<ul style="list-style-type: none"> • 1.90.465 - Driver for 28x0T autosampler: Agilent OpenLAB CDS ChemStation for GC rev. C • 1.90.456 - Driver for 28x0T autosampler: Agilent OpenLab CDS 2 • 1.90.455 - Driver for 28x0T autosampler: ICF (for Scion Compass CDS release 4.1 or higher). It requires two additional licenses: <ul style="list-style-type: none"> ○ PN 4141003 CompassCDS RC.Net control driver ○ PN 4141004 CompassCDS ICF driver <p>Contact your CDS Representative to get updated information regarding the required license to add HTA driver.</p>

Intended use	The CDS will be used to program the method(s) on the autosampler and start the sample list on the autosampler: the whole daily routine will be processed by using just the CDS.
---------------------	---

4.2.5 By CDS through third party integration

Some of the most popular multi-vendor CDS offer the possibility to control the HTA autosamplers by special control modules. Check with your CDS Representative if the multi-vendor CDS you are using provides control module for HTA autosamplers.

Additional requirement	CDS may require an additional license.
Intended use	Contact your CDS Representative.

4.3 Basic data about syringes

4.3.1 Syringe kit for Headspace technique

Name	Characteristic
Kit for 2,5ml syringe	Standard syringe kit. See ' Syringes for Headspace technique ' on page ⁴⁰² .
Kit for 1ml syringe	Syringe kit for specific syringe code. See ' Syringes for Headspace technique ' on page ⁴⁰² .

4.3.2 Syringe volume for Liquid technique


The following table shows the suggested syringe according to the desired injection volume range.

Injection volume range	Suggested syringe
0.5–5 µl	10 µl
0.2–1 µl	5 µl
0.01–0.3 µl	0.5 µl

For other ranges, consider that the injection volume must be between the 10% and 90% of the syringe nominal volume.

4.4 Structures


4.4.1 Autosampler front structure




Part	Name	Function
A	Turret	Move the syringe.
B	Oven	Heat and shake the vials.
C	Tray	Hold the standard rack with the sample vials.
D	Syringe location	Hold and move the syringe.
E	Injection area	Free space to inject the sample in the injector/s.
F	Touch screen	Interact with the autosampler.
G	(For Headspace) Syringe integrity tool location	Hold the syringe integrity tool.
	(For Liquid) Waste vial location	Hold the waste vial.
	(For SPME) Wash vial location	Hold the wash vial.


4.4.2 Autosampler rear structure

Connection panel:




Part	Name	Function
A	Ethernet port	Connect the autosampler to the network.
B	Service Ethernet reset	Only for Customer Representative.
C	Interface cable connector	Connect the autosampler to the analyzer.
D	Power supply socket	Connect to the power source. <div></div>
E	Gas connector	Allow to purge the syringe.
F	Service Emergency release	Only for Customer Representative.
G	Power switch	Switch the autosampler on and off.

4.4.3 Syringe location structure




Shipping position



Operating position

Part	Name	Function
A	Safety lock	Lock the vial locator during shipment.
B	Trolley	Move the syringe plunger.
C	Sliding lid	Protect the syringe location.
D	Needle height regulator	Allow to regulate the needle height according to the syringe in use.
E	Plunger locker	Lock the syringe plunger in the trolley.
F	Syringe purge connector	Allow to purge the syringe.
G	Syringe electrical connection	Allow to heat the syringe.
H	Vial locator	Allow to transport the vial and find the injector.
I	Finger nuts	Allow to fix the syringe.
J	LED	Illuminate the needle.
K	Screws for vial locator	Allow to fix the removable vial locator.


4.4.4 Oven structure




Part	Function
A	Holes for: <ul style="list-style-type: none">• sample vials (Headspace and SPME mode) and, if required, grey spacers• black spacers and solvent vials (Liquid mode)
B	
C	
D	
E	
F	


4.4.5 Vial locators




For Liquid technique:



O-ring colour	For...
 (green)	Liquid technique

For Headspace and SPME technique:



O-ring colour	For...
 (red)	Headspace technique
 (blue)	Headspace technique and SPME for standard fibers technique
 (yellow)	SPME for 1.1 fibers / 1.5 fibers technique

4.5 Oven location

4.5.1 Position of oven

The oven is positioned to the left or right side of the autosampler, depending on the analyzer. The correct position of the oven for your autosampler is indicated in the Configuration side field of the product paper.

Note: *This manual and the illustrations contained therein have been realized with the oven on the left side of the autosampler, seen from the front.*

4.5.2 Oven positions

The autosampler identifies each oven position with the letter of the position it occupies in the oven location.

4.5.3 Solvent vials

The autosampler identifies each solvent vial with the letter of the position it occupies in the oven location.

5. Installation

5.1 Before proceeding

Before install the autosampler:

- Ensure to have correctly installed the mounting kit following the relevant instruction manual.
- If this is the first installation of the autosampler, follow the next paragraph, otherwise see '[Re-installation in Headspace mode' on page 110](#)', '[Re-installation in Liquid mode' on page 111](#)' or '[Re-installation in SPME mode' on page 112](#)'.

5.2 Autosampler installation


NOTICE: Do not move the turret by hand.

5.2.1 Open the locking hook




CAUTION Sharp metal parts. Do not open the locking hook by hand.


1. Insert the syringe pointer in the locking hook.



2. Push down the handle of the syringe pointer so that the other end lifts.



3. Pull the syringe pointer until the locking hook opens.




5.2.2 Install the autosampler

1. Ensure the locking hook is open. Otherwise, open it.
2. Lift the autosampler by hand. For details, see ['Lift the autosampler' on page 24](#).



CAUTION Heavy product. Lift the autosampler carefully and using both hands.

3. Place the autosampler on the mounting plate so that the four feet and the central part of the autosampler perfectly fit into the four support holes and into the central hole of the mounting plate respectively.



4. To lock the autosampler, close the locking hook.

5.2.3 Connect the autosampler



CAUTION To avoid damaging the equipment and compromising its safety, only use the interface cable, power cord and power supply provided with the autosampler.

1. Ensure the power switch is to off.
2. If present, insert the interface cable provided with the mounting kit in the autosampler and in the analyzer. For details, see '[Autosampler rear structure](#)' on page 53.
3. Connect the power cord to the power supply.
4. Insert the power supply plug in the power supply socket of the autosampler.
5. Connect the power cord to the wall socket.
6. Connect the Ethernet cable to the autosampler and to the PC or wall jack/switch.
7. Connect the purge gas line at the purge connector, max pressure: 2 bar (29psi or 200kPa).

5.2.4 Start the autosampler

1. Press the autosampler power switch to switch the autosampler on.


5.2.5 Install HTA Monitor


1. Ensure the Ethernet cable is connected.
2. Download the **HTA_Monitor.exe** following the instructions included in the packaging.
3. Double click **HTA_Monitor.exe** and authorize the device modification.

Note: In case you want to use the HTA Monitor software from any user of the PC, it is necessary to choose "Install for all users".


4. Select the desired language and click **OK**: the **Setup - HTAMonitor** window appears.
5. Click **Next**, accept the agreement and click **Next**.
6. To choose a different installation folder, click **Browse....** Otherwise, click **Next > Next**.


- 7. To confirm the automatic search of the software updates, select **Next**.
- 8. Click **Install**: the software installation begins.



- 9. Click **Finish**:  appears in the taskbar and a popup opens.
- 10. If after the installation appears a Windows Firewall message asking on which network you want to work:

If...	Then...
Private is selected	click on Allow access .
Public is selected	select also Private and then click on Allow access .



- 11. Right click  and click **Install**: the **HTA Monitor** window opens.
- 12. To confirm the sharing of useful technical data with HTA, select **Yes**.
- 13. Click **Add new...**: the installation wizard appears.
Note: You can connect the autosampler only to one PC at the time.
Note: Close all the application that can communicate with the autosampler (i.e. HTA Autosampler Manager, third party software,...) during this procedure.
- 14. Follow the wizard instructions to install the autosampler on HTA Monitor.
- 15. If the autosampler IP is not unique in the network, follow the wizard instructions to assign a new IP to the autosampler. Contact your network administrator if necessary.


5.2.6 Open the HTA Monitor virtual screen

In the taskbar, right click and click **Show Virtual Screen**: the virtual screen opens.


5.2.7 Prepare the syringe location

NOTICE: Do not move the turret by hand.

- 1. If this is the first installation, the **Move the left block from the shipping position to the operating position as described in user manual** page appears. Unlock the safety lock proceeding with the next step.
- 2. Raise the sliding lid above the syringe location.
- 3. Put your finger below the vial locator: it will be necessary to sustain the vial locator in the following points of this procedure.
- 4. Take the syringe pointer and use its edge to loosen the safety lock.



5. With your finger, guide the right block down until it rests on the sledge. Remove your finger from below the vial locator.




6. To lock the syringe location in position, screw the left block 2-3 mm above the needle height regulator level using the syringe pointer.
7. Lower the sliding lid.
8. If this is the first installation, select **Continue**.


5.3 Sample rack installation

5.3.1 Insert the sample standard rack in the autosampler tray

1. Insert the front side of sample rack.





2. Insert the other side of the sample rack and push it so that the label with the letters is parallel to the base of the tray.



3. If this is the first installation, select **Continue**: the autosampler tray closes and the screen displays the autosampler model and the firmware version.

5.3.2 Rack correct positioning




5.4 Install the headspace syringe


1. When the screen displays the **Remove and install new syringe** page, do what follows:

If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the syringe installation without moving the turret	proceed to point 2.


2. Raise the sliding lid and unscrew the two finger nuts.




3. Take the syringe pointer, insert its edge in the plunger locker and pull it out.




4. Take the syringe.
 5. Push the vial locator up and insert the syringe needle in the central hole of the vial locator.



6. Insert the syringe body in the syringe holder and the plunger in the syringe plunger holder.




7. Lower the syringe plunger until it rests on the syringe plunger holder.



8. Screw the two finger nuts.

9. Insert the plunger locker with the knurled side upward in the syringe plunger holder.



10. Lower the sliding lid.
11. Select **Continue**: the **Select syringe Volume** page appears.
12. Choose the syringe volume and select **Save**: the automatic alignment procedure starts.

5.5 Install the tray

1. At the end of the alignment procedure, select **Continue**.
2. Set the depth below the septum where the needle must be positioned to aspirate the sample and select **Continue** > **Save** > **Continue**.

Note: *If you set the sample depth to 0 mm (the minimum value), the needle tip will be just below the septum. If you set the sample depth to 25 mm (the maximum value), the needle tip will be 25 mm below the septum.*

5.6 Select the analyzer

1. Select the analyzer to be used and then **OK** to confirm.
2. Select the number of injectors to be used and then **OK** to confirm and **Load**.
3. Select **Continue**.

5.7 Align the injectors

1. To align the autosampler to the front injector, select **Front Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the front injector adapter.



CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

2. Do what follows:

If you want to...	Then...
confirm the injector alignment	select Confirm .
align the autosampler to the injector manually	<ol style="list-style-type: none"> 1. Select Do Manual. 2. Select the arrows to move the vial locator to the correct position (Touch Sensor must be ON). 3. Select Store to confirm.

3. Set the desired injection depth: the autosampler tries to insert the needle in the injector.
4. Select **Save** to confirm.
5. Select the desired injection speed (from **Very Low** to **Very High**) and select **Save**: the turret returns to its original position and the **Alignment** menu appears.
6. If the rear injector is present, select **Rear Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the rear injector adapter.







CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

7. Repeat from step 2 to 7.

5.7.1 Vial locator aligned to the injector adapter

To be aligned, the vial locator must lean against the top of the injector adapter completely.





5.8 Set the position of the vials - The Touch and Plunger Zero procedure


1. Select **Touch & Plunger Zero**: the turret moves and the tray opens.
2. Insert a crimped sample vial in position A4 of the rack and in C position of the oven (if using 6ml or 10 ml vial, ensure to insert first the proper spacer in C position). In case syringe integrity tool is available, insert it in its position.
3. Select **Continue**: the vial locator touches the syringe integrity tool, the sample vial in the rack, the sample vial in the oven, the cover and the injector adapter to memorize their positions.
4. Remove the vial from oven and press **Continue**.
5. At the end of the procedure, when the **Setup menu** appears, select **Exit > Exit**.

5.9 Set the communication with the analyzer

✱ *Run menu > Settings > Setup > Configuration > Input Output*




1. Set **Ready On** and **Start On** as specified in the Mounting kit manual.
2. Select **Exit**: the **Configuration** menu appears.

5.10 Test the communication with the analyzer

1. Ensure the autosampler and the analyzer are switched on and the interface cable is properly connected.
2. In the **Run** menu select **Settings > Setup > Manual Operation**.
3. Set the analyzer to be not ready to receive the injection and check that the **Ready** value on the screen is **Off**.
4. Set the analyzer to be ready to receive the injection and check that the **Ready** value on the screen is **On**.
5. Select  to select **Start**.
6. Select **OK** and check that the **Ready** value on the screen goes to **Off** and that analyzer starts the analysis.

5.11 Set the autosampler in Headspace mode

1. Choose how do you want to operate the autosampler. See '[How to operate the autosampler](#)' on page [46](#) for details.
2. Do what follows:


If you want to operate the autosampler...	Then...
by autosampler touch screen: all samples are processed in the same way	Set the method and the sequence selecting A.1 as First Sample and G.6 as Last Sample . See ' Set the method and sequence ' on page 126 .
by autosampler touch screen: you are required to process samples in different ways	<ol style="list-style-type: none"> 1.  <i>Run menu > Settings > Setup > Configuration > Run</i> 2. Select User Interface and set the parameter to Normal.
by HTA Autosampler Manager	Install HTA Autosampler Manager and connect the autosampler. Full license is required. See ' HTA Autosampler Manager installation ' on page 76 .
by HTA Autosampler Manager – CFR 21 part 11 version	Install HTA Autosampler Manager and connect the autosampler. Full license and CFR 21 part 11 license are required. See ' HTA Autosampler Manager installation ' on page 76 .
by CDS through HTA software integration	<ol style="list-style-type: none"> 1. See 'HTA Autosampler Manager installation' on page 76. Full license is required. 2. Perform the setup and preliminary operations on your CDS to connect with the autosampler. For further details, see your CDS documentation and HTA Autosampler Manager help.
by CDS through HTA software connectors	<ol style="list-style-type: none"> 1. See 'Install HTA Autosampler Manager' on page 76. 2.  <i>Run menu > Settings > Setup > Configuration > Activation</i> 3. Select the feature indicated in the Password for Optional Features purchased and select OK. 4. Select Full and insert the password indicated in the Password for Optional Features purchased. 5.  <i>Run menu > Settings > Setup > Configuration > Run</i> 6. Select Protocol version and set the parameter to -. 7. Perform the setup and preliminary operations on your CDS to connect with the autosampler. For further details, see your CDS documentation and HTA Autosampler Manager help. <p>Dedicated feature is required.</p>
by CDS through third party integration	Perform the setup and preliminary operations on your third party integration software to connect with the autosampler. See your third party software documentation for further details.

5.12 HTA Autosampler Manager installation

5.12.1 Install HTA Autosampler Manager

1. Ensure the Ethernet cable is connected.
2. Download the **HTA_Autosampler_Manager.exe** following the instruction included in the packaging.
3. Double click **HTA_Autosampler_Manager.exe** and authorize the device modification.
4. Select the desired language of the installer and click **OK**: the **Setup - HTA Autosampler Manager** window appears.
5. Accept the agreement and click **Next**.
6. To choose a different installation folder, click **Browse....** Otherwise, click **Next**.
7. Do what follows:

If you have...	Then...	And...
the standard version of HTA Autosampler Manager (trial or full)	select Standard version	click Next > Next > Next > Install .
the CFR 21 Part 11 version of HTA Autosampler Manager	select CFR 21 Part 11 version	

8. Click **Next**, automatically USB driver installation procedure start, in case of Windows 10 or Windows 8 this procedure could fail, ignore the status information reported at the end of the procedure then click **Finish > Finish: START - HTA Autosampler Manager** opens and  appears in the taskbar.
9. To proceed, see ['Connect HTA Autosampler Manager to the autosampler' on page 76](#).

5.12.2 Connect HTA Autosampler Manager to the autosampler

1. Do what follows:

If you have...	Then...
the standard version of HTA Autosampler Manager (trial or full)	Proceed with step 2.
the CFR 21 Part 11 version of HTA Autosampler Manager	<ol style="list-style-type: none"> 1. Click Install > Licenses > Add...: the CFR 21 Part 11 Licenses wizard appears. 2. Follow the wizard instructions. When required, insert the CFR 21 Part 11 activation key of your license. 3. Click Install > Users > Add...: the CFR 21 Part 11 Users wizard appears. 4. Follow the wizard instructions adding at least one Administrator user. 5. Proceed with step 2.

2. In **START - HTA Autosampler Manager**, to add the autosampler in use to the **Instruments** list click **Install > Add new**: the Instrument wizard appears.
3. Follow the wizard instructions. When required, insert the Instrument activation key of your license (trial or full).
4. If the online activation fails, for example because the Internet connection is missing, select **Activate offline > I need a license key > Next**, then choose the way to get the license key you need (**By your**

smartphone is the fastest and recommended option). When you receive the license key, repeat the entire operation but take care to select **I have a license key** and to insert the license key just received.

5. If the autosampler is connected to a LAN, ensure the default connection parameters of the autosampler are not in conflict with other devices.
6. If the autosampler is connected to a wide LAN, provide the default connection parameters of the autosampler to your network administrator to ensure they are not in conflict with other devices.
7. If the **Connection troubleshooting** windows appears, do what follows:

If your autosampler is...	Then...						
directly connected to the PC	<ol style="list-style-type: none"> 1. Open the Control Panel of the PC and click Network and Internet > Network and Sharing Center > Change adapter settings. 2. Right click Local Area Connection > Properties > Internet Protocol Version 4 (TCP/IPv4) > Properties. 3. Select Use the following IP address and insert what follows: <table border="1"> <tr> <td>IP address</td><td>192.168.0.208</td></tr> <tr> <td>Subnet mask</td><td>255.255.0.0</td></tr> <tr> <td>Default gateway</td><td>192.168.0.1</td></tr> </table> 4. Click OK > Close. 5. In the Connection troubleshooting windows select Proceed to retry the autosampler connection. 6. Follow the wizard instructions to complete the installation. 	IP address	192.168.0.208	Subnet mask	255.255.0.0	Default gateway	192.168.0.1
IP address	192.168.0.208						
Subnet mask	255.255.0.0						
Default gateway	192.168.0.1						
connected to a LAN or a WLAN	<ol style="list-style-type: none"> 1. Ask your network administrator for new IP address, Subnet mask and Gateway parameters. 2. Select Modify network configuration of autosampler and click Proceed. 3. In Ethernet insert the new parameters provided by your network administrator and click Next. 4. Restart the autosampler and wait until it is in standby, then click Next. 5. Follow the wizard instructions to complete the installation. 						



7. Click **Run**, select the desired autosampler and click **Open**: HTA Autosampler Manager connects to the autosampler.

5.13 Change to Liquid mode



5.13.1 Get out of Headspace mode

✱ *Run menu > Change Mode*

1. Select **Liquid Mode** and select **Continue**.
2. Lift the sliding lid.
3. Remove the vial locator.




4. Select **Continue**.
5. Take the syringe pointer, insert its edge in the plunger locker and pull it out.



6. Unscrew the two finger nuts.


7. Remove the syringe warmer assembly.




8. Select **Continue**.

5.13.2 Get into Liquid mode


1. Select the Syringe volume and select **Save**.
2. If necessary, raise the sliding lid and remove the two finger nuts.
3. Install the syringe holder.



4. Fix completely the two finger nuts.
5. If necessary, take the syringe pointer, insert its edge in the plunger locker and pull it out.
6. Insert the plunger extension in the plunger holder, taking care to keep frontwards the side with the hole for the syringe plunger.



7. Insert the plunger lock.




Install the syringe



1. When the screen displays the **Remove and install new syringe** page, do what follows:

If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the syringe installation without moving the turret	proceed to point 2.


2. Open the syringe locker.





3. Take the syringe.
4. Insert the syringe body in the syringe holder and the plunger in the plunger extension.




5. Close the syringe locker.




6. Select **Continue**: the **Select syringe Volume [μl]** page appears.
7. Choose the syringe volume and select **Save**: the automatic alignment procedure starts.
8. Wait for the automatic positing.
9. Install the vial locator for Liquid technique (green o-ring: see ['Vial locators' on page 55](#)), taking care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars.




10. Fasten completely the two fixing screws, taking care to have the small bars completely inside the vial locator bars.



11. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the syringe.



12. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.




13. Lower the sliding lid.


14. Press **Continue**.

Insert the sample standard rack in the autosampler tray

1. Insert the front side of sample rack.





2. Insert the other side of the sample rack and push it so that the label with the letters is parallel to the base of the tray.



3. If this is the first installation, select **Continue**: the autosampler tray closes and the screen displays the autosampler model and the firmware version.

Rack correct positioning



Install the tray


1. At the end of the alignment procedure, select **Continue**.
2. Select **OK**.
3. If necessary, modify the vial type, then select **Continue**.
4. Set the depth below the septum where the needle must be positioned to aspirate the sample and select **Continue** > **Save** > **Continue**.

Note: If you set the sample depth to 0 mm (the minimum value), the needle tip will be just below the septum. If you set the sample depth to 32 mm (the maximum value), the needle tip will be close to the vial bottom.

Note: The autosampler shows automatically the vial depth used the last time that the autosampler has been used in Liquid mode.

Change the analyzer to be used

✱ *Run menu > Settings > Setup > Alignment > Analyzer Default*

1. If the name of the analyzer to be used is already listed, select it. Otherwise, select **Generic**.
2. Select **OK** to confirm.
3. Select the number of injectors to be used and then **OK** to confirm and **Load**.
4. Select **Continue**.
5. Proceed with the ['Align injectors' procedure on the following section](#). 

Realign injectors

✱ *Run menu > Settings > Setup > Alignment*

1. To align the autosampler to the front injector, select **Front Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the front injector adapter.



CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

2. Do what follows:

If you want to...	Then...
confirm the injector alignment	select Confirm .
align the autosampler to the injector manually	<ol style="list-style-type: none"> 1. Select Do Manual. 2. Select the arrows to move the vial locator to the correct position (Touch Sensor must be ON). 3. Select Store to confirm.

3. Set the desired injection depth: the autosampler tries to insert the needle in the injector.
4. Select **Save** to confirm.
5. Select the desired injection speed (from **Very Low** to **Very High**) and select **Save**: the turret returns to its original position and the **Alignment** menu appears.
6. If the rear injector is present, select **Rear Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the rear injector adapter.





CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

7. Repeat from step 2 to 7.

Vial locator aligned to the injector adapter

To be aligned, the vial locator must lean against the top of the injector adapter completely.








Set the position of the vials - The Touch and Plunger Zero procedure

1. Select **Touch & Plunger Zero**: the turret moves and the tray opens.
2. Insert a sample vial with cap in position A6 of the standard rack, a black spacer and a solvent vial in position C of the oven and ensure the waste vial is in its proper position.
3. Select **Continue**: the vial locator touches the waste vial, the sample vial, the solvent vial and the injector adapter to memorize their positions.
4. At the end of the procedure, when the **Setup menu** appears, select **Exit > Exit**.

5.14 Set the autosampler in Liquid mode

1. Choose how do you want to operate the autosampler. See ['How to operate the autosampler' on page 46](#) for details.
2. Do what follows:

If you want to operate the autosampler...	Then...
by autosampler touch screen: all samples are processed in the same way	Set the method and the sequence selecting A.1 as First Sample and K.11 as Last Sample . See 'Set the method and sequence' on page 168 .
by autosampler touch screen: you are required to process samples in different ways	<ol style="list-style-type: none"> 1.  <i>Run menu > Settings > Setup > Configuration > Run</i> 2. Select User Interface and set the parameter to Normal.
by HTA Autosampler Manager	Install HTA Autosampler Manager and connect the autosampler. Full license is required. See 'HTA Autosampler Manager installation' on page 76 .
by HTA Autosampler Manager – CFR 21 part 11 version	Install HTA Autosampler Manager and connect the autosampler. Full license and CFR 21 part 11 license are required. See 'HTA Autosampler Manager installation' on page 76 .
by CDS through HTA software integration	<ol style="list-style-type: none"> 1. See 'HTA Autosampler Manager installation' on page 76. Full license is required. 2. Perform the setup and preliminary operations on your CDS to connect with the autosampler. For further details, see your CDS documentation and HTA Autosampler Manager help.
by CDS through HTA software connectors	<ol style="list-style-type: none"> 1. See 'Install HTA Autosampler Manager' on page 76. 2.  <i>Run menu > Settings > Setup > Configuration > Activation</i> 3. Select the feature indicated in the Password for Optional Features purchased and select OK. 4. Select Full and insert the password indicated in the Password for Optional Features purchased. 5.  <i>Run menu > Settings > Setup > Configuration > Run</i> 6. Select Protocol version and set the parameter to -. 7. Perform the setup and preliminary operations on your CDS to connect with the autosampler. For further details, see your CDS documentation and HTA Autosampler Manager help. <p>Dedicated feature is required.</p>
by CDS through third party integration	Perform the setup and preliminary operations on your third party integration software to connect with the autosampler. See your third party software documentation for further details.

5.15 Enable the SPME feature for standard fibers

✱ *Run menu > Settings > Setup > Configuration > Activation*

1. Select the feature number 1 and select **OK**.
2. Select **FULL**.
3. Insert the code you find in the envelope of the SPME kit - standard fibers and press **OK**.

5.16 Enable the SPME feature for 1.1 / 1.5 fibers

✱ *Run menu > Settings > Setup > Configuration > Activation*



1. Select the feature number 9 and select **OK**.
2. Select **FULL**.
3. Insert the code you find in the envelope of the SPME kit - 1.1 fibers / 1.5 fibers and press **OK**.

5.17 Change to SPME mode



5.17.1 Get out of Liquid mode


✱ *Run menu > Change Mode*

1. Select **SPME mode** and select **Continue**.
2. Wait for the automatic positioning.
3. Remove the two fixing screws to remove the vial locator.




4. Select **Continue**.
5. Wait for the automatic positioning, open the locker and remove the syringe.







6. Take the syringe pointer, insert its edge in the plunger locker and pull it out.



7. Remove the plunger extension.



8. Unscrew the two fixing nuts and remove the syringe holder.



9. Select **Continue**.


5.17.2 Get into SPME mode (fiber)

- 1. Select the **Type** of fiber you are going to use.
- 2. When the screen displays the **Install new fiber holder or Sealed tip holder** page, do what follows:


If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the fiber holder installation without moving the turret	proceed to point 3.

- 3. Wait for the automatic positioning.


4. Unscrew the fiber guide and remove it.




5. Screw the fiber fixing point to the fiber holder plunger.



6. Lift up the fiber holder plunger until the fiber is completely inside the fiber holder.




7. Insert the fiber guide along the fiber and fix it at the fiber holder.




Note: Take care to use the standard fiber guide in case of standard fiber or the sealed tip guide in case of sealed tip.


8. Install the fiber holder.



9. Fix completely the two finger nuts.




10. Insert the plunger locker.




11. Install the vial locator for Headspace and SPME standard fibers technique (blue o-ring): see '[Vial locators](#)' on page 55. Take care to let the needle tip to enter the central hole on the vial locator and


to insert the two bars along the two long bars and then select **Continue**.



12. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the fiber.




13. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.




Insert the sample standard rack in the autosampler tray

1. Remove the sample vials and the 121 positions sample rack.
2. Insert the front side of the 42 positions sample rack.




3. Insert the other side of the sample rack and push it so that the label with the letters is parallel to the base of the tray.



4. If this is the first installation, select **Continue**: the autosampler tray closes and the screen displays the autosampler model and the firmware version.

Rack correct positioning



Change the analyzer to be used

✱ *Run menu > Settings > Setup > Alignment > Analyzer Default*

1. If the name of the analyzer to be used is already listed, select it. Otherwise, select **Generic**.
2. Select **OK** to confirm.
3. Select the number of injectors to be used and then **OK** to confirm and **Load**.
4. Select **Continue**.
5. Proceed with the ['Align injectors' procedure on the following section](#). 70

Realign injectors

✱ *Run menu > Settings > Setup > Alignment*

1. To align the autosampler to the front injector, select **Front Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the front injector adapter.



CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

2. Do what follows:

If you want to...	Then...
confirm the injector alignment	select Confirm .
align the autosampler to the injector manually	<ol style="list-style-type: none"> 1. Select Do Manual. 2. Select the arrows to move the vial locator to the correct position (Touch Sensor must be ON). 3. Select Store to confirm.

3. Set the desired injection depth: the autosampler tries to insert the needle in the injector.
4. Select **Save** to confirm.
5. Select the desired injection speed (from **Very Low** to **Very High**) and select **Save**: the turret returns to its original position and the **Alignment** menu appears.
6. If the rear injector is present, select **Rear Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the rear injector adapter.




CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

7. Repeat from step 2 to 7.

Vial locator aligned to the injector adapter

To be aligned, the vial locator must lean against the top of the injector adapter completely.



Set the position of the vials - The Touch and Plunger Zero procedure


1. Select **Touch & Plunger Zero**: the turret moves and the tray opens.
2. Insert a crimped sample vial in position A4 of the rack and in C position of the oven (if using 6ml or 10 ml vial, ensure to insert first the proper spacer in C position). In case wash vial is used, insert it in its position.
3. Select **Continue**: the vial locator touches the wash vial, the sample vial in the rack, the sample vial in the oven, the cover and the injector adapter to memorize their positions.
4. Remove the vial from oven and press **Continue**.
5. At the end of the procedure, when the **Setup menu** appears, select **Exit > Exit**.

5.17.3 Get into SPME mode (sealed tip)


1. Select the **Type** of sealed tip you are going to use.
2. When the screen displays the **Install new fiber holder or Sealed tip holder** page, do what follows:

If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the fiber holder installation without moving the turret	proceed to point 3.


- 3. Wait for the automatic positioning.
- 4. Insert the sealed tip arrow into the sealed tip holder.




- 5. Insert the sealed tip spacer in the position shown in the following pictures.




6. Install the fiber holder.




7. Fix completely the two finger nuts.




8. Insert the plunger extension.




9. Lower the plunger extension until it rests on the plunger holder.




10. Insert the plunger locker.




11. Install the vial locator for SPME for 1.1 fibers / 1.5 fibers technique (yellow o-ring): see ['Vial locators' on page 55](#). Take care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars and then select **Continue**.



12. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the fiber.




13. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.




Insert the sample standard rack in the autosampler tray

1. Remove the sample vials and the 121 positions sample rack.
2. Insert the front side of the 42 positions sample rack.





3. Insert the other side of the sample rack and push it so that the label with the letters is parallel to the base of the tray.




4. If this is the first installation, select **Continue**: the autosampler tray closes and the screen displays the autosampler model and the firmware version.

Rack correct positioning



Change the analyzer to be used

✱ *Run menu > Settings > Setup > Alignment > Analyzer Default*

1. If the name of the analyzer to be used is already listed, select it. Otherwise, select **Generic**.
2. Select **OK** to confirm.
3. Select the number of injectors to be used and then **OK** to confirm and **Load**.
4. Select **Continue**.
5. Proceed with the ['Align injectors' procedure on the following section](#). 

Realign injectors

✱ *Run menu > Settings > Setup > Alignment*

1. To align the autosampler to the front injector, select **Front Injector**: the turret moves and the vial

locator touches a known part of the autosampler and then the front injector adapter.



CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

2. Do what follows:

If you want to...	Then...
confirm the injector alignment	select Confirm .
align the autosampler to the injector manually	<ol style="list-style-type: none"> 1. Select Do Manual. 2. Select the arrows to move the vial locator to the correct position (Touch Sensor must be ON). 3. Select Store to confirm.

3. Set the desired injection depth: the autosampler tries to insert the needle in the injector.
4. Select **Save** to confirm.
5. Select the desired injection speed (from **Very Low** to **Very High**) and select **Save**: the turret returns to its original position and the **Alignment** menu appears.
6. If the rear injector is present, select **Rear Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the rear injector adapter.






CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

7. Repeat from step 2 to 7.

Vial locator aligned to the injector adapter

To be aligned, the vial locator must lean against the top of the injector adapter completely.






Set the position of the vials - The Touch and Plunger Zero procedure

1. Select **Touch & Plunger Zero**: the turret moves and the tray opens.
2. Insert a crimped sample vial in position A4 of the rack and in C position of the oven (if using 6ml or 10 ml vial, ensure to insert first the proper spacer in C position). In case wash vial is used, insert it in its position.
3. Select **Continue**: the vial locator touches the wash vial, the sample vial in the rack, the sample vial in the oven, the cover and the injector adapter to memorize their positions.
4. Remove the vial from oven and press **Continue**.
5. At the end of the procedure, when the **Setup menu** appears, select **Exit > Exit**.

5.18 Set the autosampler in SPME mode

1. Choose how do you want to operate the autosampler. See '[How to operate the autosampler](#)' on page ⁴⁶ for details.
2. Do what follows:

If you want to operate the autosampler...	Then...
by autosampler touch screen: all samples are processed in the same way	Set the method and the sequence selecting A.1 as First Sample and G.6 as Last Sample . See ' Set the method and sequence ' on page ²⁰⁵ .
by autosampler touch screen: you are required to process samples in different ways	<ol style="list-style-type: none"> 1.  <i>Run menu > Settings > Setup > Configuration > Run</i> 2. Select User Interface and set the parameter to Normal.
by HTA Autosampler Manager	Install HTA Autosampler Manager and connect the autosampler. Full license is required. See ' HTA Autosampler Manager installation ' on page ⁷⁶ .
by HTA Autosampler Manager – CFR 21 part 11 version	Install HTA Autosampler Manager and connect the autosampler. Full license and CFR 21 part 11 license are required. See ' HTA Autosampler Manager installation ' on page ⁷⁶ .
by CDS through HTA software integration	<ol style="list-style-type: none"> 1. See 'HTA Autosampler Manager installation' on page ⁷⁶. Full license is required. 2. Perform the setup and preliminary operations on your CDS to connect with the autosampler. For further details, see your CDS documentation and HTA Autosampler Manager help.
by CDS through third party integration	Perform the setup and preliminary operations on your third party integration software to connect with the autosampler. See your third party software documentation for further details.

5.19 Accessories installation

External pressure regulator (optional)

This accessory is used to regulate the pressure of the gas used for the syringe flushing (e.g. Nitrogen). The regulator input pressure must be lower than 10 bar (145 psi or 1000 kPa).

1. Set the pressure regulator output at 2 bar maximum (29 psi or 200 kPa).
2. Insert the pressure regulator between the gas source and the gas inlet of the autosampler in order to regulate the inlet pressure.
3. Connect the 1/8" tube coming out from the regulator to the gas inlet 1/8" located on the rear panel of the autosampler.
4. Connect the other end of the regulator to the gas source.

Swagelok Adapter

1. Unscrew the Rapid fitting terminal using an Hexagon key (2mm).
2. Put the O.R. on the new fitting.
3. Screw the fitting and the washer on the female adapter using an open-end wrench (8mm).






System integrity tool for Headspace technique

1. Remove the upper part of the Syringe integrity tool.
2. Insert the septum in the lower part of the Syringe integrity tool.
3. Screw by hands the upper part of the Syringe integrity tool to the lower part.
4. Place the Syringe integrity tool in its proper location.
5. Set the position of the vial. For more details, see ['Set the position of the vials - The Touch and Plunger Zero procedure' on page 72](#).

5.20 Re-installation or setup after a long downtime in Headspace mode

This procedure must be followed if:






- the autosampler has not been used for a long time;
- the autosampler must be installed on the same analyzer where it was previously installed;
- the autosampler must be installed on a different analyzer.

Step	Installation task	Procedure	
1	Start the autosampler	See 'Start the autosampler' on page 61 .	
2	Prepare the syringe location	If...	Then...
		the autosampler is in shipping position	See 'Prepare the syringe location' on page 62 .
		the autosampler is in operating position	Proceed with the next step.
3	Sample rack installation	 <i>Run menu > Load Sample Tray</i> Insert the sample standard rack in the autosampler tray. See 'Sample rack installation' on page 64 .	
4	Install the syringe	See 'Syringe management' on page 131 .	
5	Install the tray	 <i>Run menu > Settings > Setup > Alignment > Tray Install</i> See 'Racks management' on page 130 .	
6	Select the analyzer	 <i>Run menu > Settings > Setup > Alignment > Analyzer Default</i> See 'Select the analyzer' on page 69 .	
7	Align the injectors	 <i>Run menu > Settings > Setup > Alignment</i> See 'Align the injectors' on page 70 .	
8	Set the position of the vials - The Touch and Plunger Zero procedure	 <i>Run menu > Settings > Setup > Alignment</i> See 'Set the position of the vials - The Touch and Plunger Zero procedure' on page 121 .	
9	HTA Monitor installation	See 'Install HTA Monitor' on page 61 .	
10	Set the communication with the analyzer	See 'Set the communication with the analyzer' on page 73 .	
11	Test the communication with the analyzer	See 'Test the communication with the analyzer' on page 74 .	
12	Set the autosampler	See 'Set the autosampler' on page 75 .	

5.21 Re-installation or setup after a long downtime in Liquid mode

This procedure must be followed if:






- the autosampler has not been used for a long time;
- the autosampler must be installed on the same analyzer where it was previously installed;
- the autosampler must be installed on a different analyzer.

Step	Installation task	Procedure	
1	Start the autosampler	See 'Start the autosampler' on page 61 .	
2	Prepare the syringe location	If...	Then...
		the autosampler is in shipping position	See 'Prepare the syringe location' on page 62 .
		the autosampler is in operating position	Proceed with the next step.
3	Sample rack installation	 <i>Run menu > Load Sample Tray</i> Insert the sample standard rack in the autosampler tray. See 'Sample rack installation' on page 157 .	
4	Install the syringe	See 'Syringe management' on page 173 .	
5	Install the tray	 <i>Run menu > Settings > Setup > Alignment > Tray Install</i> See 'Racks management' on page 172 .	
6	Select the analyzer	 <i>Run menu > Settings > Setup > Alignment > Analyzer Default</i> See 'Select the analyzer' on page 69 .	
7	Align the injectors	 <i>Run menu > Settings > Setup > Alignment</i> See 'Align the injectors' on page 70 .	
8	Set the position of the vials - The Touch and Plunger Zero procedure	 <i>Run menu > Settings > Setup > Alignment</i> See 'Set the position of the vials - The Touch and Plunger Zero procedure' on page 159 .	
9	HTA Monitor installation	See 'Install HTA Monitor' on page 61 .	
10	Set the communication with the analyzer	See 'Set the communication with the analyzer' on page 73 .	
12	Test the communication with the analyzer	See 'Test the communication with the analyzer' on page 74 .	
12	Set the autosampler	See 'Set the autosampler' on page 88 .	

5.22 Re-installation or setup after a long downtime in SPME mode

This procedure must be followed if:

- the autosampler has not been used for a long time;
- the autosampler must be installed on the same analyzer where it was previously installed;
- the autosampler must be installed on a different analyzer.

Step	Installation task	Procedure	
1	Start the autosampler	See 'Start the autosampler' on page 61 .	
2	Prepare the syringe location	If...	Then...
		the autosampler is in shipping position	See 'Prepare the syringe location' on page 62 .
		the autosampler is in operating position	Proceed with the next step.
3	Sample rack installation	 <i>Run menu > Load Sample Tray</i> Insert the sample standard rack in the autosampler tray. See 'Sample rack installation' on page 64 .	
4	Install the fiber	See 'Fiber management' on page 210 .	
5	Install the tray	 <i>Run menu > Settings > Setup > Alignment > Tray Install</i> See 'Racks management' on page 209 .	
6	Select the analyzer	 <i>Run menu > Settings > Setup > Alignment > Analyzer Default</i> See 'Select the analyzer' on page 69 .	
7	Align the injectors	 <i>Run menu > Settings > Setup > Alignment</i> See 'Align the injectors' on page 70 .	
8	Set the position of the vials - The Touch and Plunger Zero procedure	 <i>Run menu > Settings > Setup > Alignment</i> See 'Set the position of the vials - The Touch and Plunger Zero procedure' on page 193 .	
9	HTA Monitor installation	See 'Install HTA Monitor' on page 61 .	
10	Set the communication with the analyzer	See 'Set the communication with the analyzer' on page 73 .	
11	Test the communication with the analyzer	See 'Test the communication with the analyzer' on page 74 .	
12	Set the autosampler	See 'Set the autosampler' on page 108 .	

5.23 Enable the temporary feature

✦ *Run menu > Settings > Setup > Configuration > Activation*

1. Select the desired feature and select **OK**.
2. Select **Temporary**: a random number appears.
3. Do what follows:



If...	Then...
you do not have the temporary activation key	<ol style="list-style-type: none"> 1. Send the random number to your Customer representative to obtain the key for the temporary activation. 2. While waiting for the temporary activation key, select Continue > Cancel.
you have the temporary activation key	<ol style="list-style-type: none"> 1. Select ACTIVATE. 2. Insert the received temporary activation key. 3. Read the instructions on the HTA Autosampler Manager help to correctly use the feature.


6. Use in Headspace mode

6.1 Get into Headspace mode


To get into Headspace mode, you need first to get out of the mode you are currently using: therefore, please follow the ['Get out of Liquid mode' on page 184](#) or ['Get out of SPME mode' on page 237](#).

- 1. Select the syringe volume and select **Save**.
- 2. Wait for the automatic positioning and install the complete HS syringe warmer.


If the vial locator is...	Then...
already mounted	<div>1. Push the vial locator up and insert the syringe needle in the central hole of the vial locator.</div> <div></div> <div>2. Insert the syringe body in the syringe holder and the plunger in the syringe plunger holder.</div> <div></div>
not mounted	Insert the syringe body in the syringe holder and the plunger in the syringe plunger holder.

If the vial locator is...	Then...
	

3. Lower the syringe plunger until it rests on the syringe plunger holder.








4. Fix the finger nuts.
5. Insert the plunger locker with the knurled side upward in the syringe plunger holder.






6. Select **Continue**.




7.


If...	Then...
the autosampler asks to mount the vial locator and you are provided with the vial locator for Headspace and SPME standard fibers technique (blue o-ring: see 'Vial locators' on page 55)	<div>1. Install the vial locator for Headspace and SPME standard fibers technique (blue o-ring: see 'Vial locators' on page 55), taking care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars. Note: Both the blue o-ring vial locator and the red o-ring vial locator can be used in Headspace technique. See 'Vial locators' on page 55.</div> <div></div> <div>2. Fasten completely the two fixing screws, taking care to have the small bars completely inside the vial locator bars.</div> <div></div> <div>3. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the syringe.</div>

If...	Then...
	 <p data-bbox="403 786 1276 887">4. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.</p>  <p data-bbox="403 1417 632 1447">5. Press Continue.</p>
<p data-bbox="148 1458 389 1843">the autosampler asks to mount the vial locator and you are not provided with the vial locator for Headspace and SPME standard fibers technique (blue o-ring: see 'Vial locators' on page 55)</p>	<p data-bbox="403 1458 1331 1559">1. Install the vial locator for Headspace technique (red o-ring: see 'Vial locators' on page 55), taking care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars.</p> 

If...	Then...
	<p data-bbox="403 259 1355 327">2. Fasten completely the two fixing screws, taking care to have the small bars completely inside the vial locator bars.</p>  <p data-bbox="403 887 1355 954">3. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the syringe.</p>  <p data-bbox="403 1476 1355 1576">4. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.</p> 

If...	Then...
	5. Press Continue .
the autosampler asks to remove the sample vials	proceed to point 8.

8.	
If the autosampler asks to...	Then...
remove the sample vials and the 121 positions rack is installed	<p>proceed to sample rack installation:</p> <ol style="list-style-type: none">1. Remove the sample vials and the 121 positions sample rack.2. Insert the front side of sample rack.  <ol style="list-style-type: none">3. Insert the other side of the sample rack and push it so that the label with the letters is parallel to the base of the tray.  <ol style="list-style-type: none">4. Check that rack is correctly positioned. 

If the autosampler asks to...	Then...
	
remove the sample vials and the 42 positions rack is installed	remove the sample vials.

6.1.1 Install the tray

1. Select **Continue**.
2. Set the depth below the septum where the needle must be positioned to aspirate the sample and select **Continue** > **Save** > **Continue**.

Note: If you set the sample depth to 0 mm (the minimum value), the needle tip will be just below the septum. If you set the sample depth to 25 mm (the maximum value), the needle tip will be 25 mm below the septum.

Note: The autosampler shows automatically the vial depth used the last time that the autosampler has been used in Headspace mode.

6.1.2 Set the position of the vials - The Touch and Plunger Zero procedure

1. Insert a crimped sample vial in position A4 of the rack and in C position of the oven (if using 6ml or 10 ml vial, ensure to insert first the proper spacer in C position). In case syringe integrity tool is available, insert it in its position.
2. Select **Continue**: the vial locator touches the syringe integrity tool, the sample vial in the rack, the sample vial in the oven, the cover and the injector adapter to memorize their positions.
3. Remove the vial from oven and press **Continue**.

6.2 Prepare the sample vials

6.2.1 Basic information on sample vials


- Vials must be approved for Headspace techniques;
- Rounded end;
- Clear or amber glass (amber glass is suitable for light sensitive samples). For more details see ['Sample vial specifications' on page 395](#).

6.2.2 Close a sample vial with a crimp cap

You must have a crimper available. Fit the septum on the cap and then:

1. Clean the inside surfaces of the crimper jaws.
2. Place the crimp cap over the top of the vial.
3. Hold the vial and place the crimper on top of it. Squeeze the handle until it touches the adjuster screw.

Note: Following picture shows a not correct crimping.



4. Check each vial for proper crimping:
 - Check the aluminum cap is not folded or wrinkled.
 - The cap should be tight and not easy to move or rotate.
 - If either of the above problems occur crimp the cap again. Regulate the crimper to adjusting the tightness if needed.
 - Check that each cap has a flat septum and that it is centered over the top of the vial:

If...	Then...
the septum is not flat	remove the cap, turn the crimper adjusting screw and try again.
the cap is not centered	remove the cap and make sure the new cap is flat on the top of the vial before the crimp operation.

6.3 Operating Modes

6.3.1 Operate by autosampler touch screen

To...	Then...						
Set the method and sequence	Do what follows: <table> <tr> <th>If...</th><th>Then...</th></tr> <tr> <td>all samples are processed in the same way</td><td>you don't need to change method and sequence at each run</td></tr> <tr> <td>you are required to process samples in different way</td><td>see 'Set the method and sequence' on page 126.</td></tr> </table>	If...	Then...	all samples are processed in the same way	you don't need to change method and sequence at each run	you are required to process samples in different way	see 'Set the method and sequence' on page 126 .
If...	Then...						
all samples are processed in the same way	you don't need to change method and sequence at each run						
you are required to process samples in different way	see 'Set the method and sequence' on page 126 .						
Run sample(s) - Before proceeding	See 'Run sample(s) - Before proceeding' on page 127 .						
Run sample(s) - Start a sequence	See 'Run sample(s) - Start a sequence' on page 127 .						
Run sample(s) - Stop a sequence	See 'Run sample(s) - Stop a sequence' on page 128 .						
Run sample(s) - Process a single sample	<table> <tr> <th>If...</th><th>Then...</th></tr> <tr> <td>all samples are processed in the same way</td><td>see 'Run sample(s) - Start a sequence' on page 127.</td></tr> <tr> <td>you are required to process samples in different way</td><td>see 'Run sample(s) - Process a single sample' on page 127.</td></tr> </table>	If...	Then...	all samples are processed in the same way	see 'Run sample(s) - Start a sequence' on page 127 .	you are required to process samples in different way	see 'Run sample(s) - Process a single sample' on page 127 .
If...	Then...						
all samples are processed in the same way	see 'Run sample(s) - Start a sequence' on page 127 .						
you are required to process samples in different way	see 'Run sample(s) - Process a single sample' on page 127 .						
Perform additional syringe purge	See 'Perform additional syringe purges' on page 128 .						
Modify the autosampler general settings	See 'Modify the autosampler general settings' on page 128 .						
Manage racks	See 'Racks management' on page 130 .						

6.3.2 Operate by HTA Autosampler Manager

To...	Then...
Set the method and sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Method , Sequence and Run topics.
Run sample(s) - Before proceeding	See 'Run sample(s) - Before proceeding' on page 127 .
Run sample(s) - Start a sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Run topic.
Run sample(s) - Stop a sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Run topic.
Run sample(s) - Process a single sample	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Auxiliary > Single Injection topic.
Perform additional syringe purge	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Auxiliary > Purge topic.
Modify the autosampler general settings	Open the autosampler with HTA Autosampler Manager and see

To...	Then...
	the help by pressing F1 . See in particular Setup topic.
Manage racks	See 'Racks management' on page 130 .

6.3.3 Operate by CDS through HTA software integration

To...	Then...	
Set the method and sequence	If you want to work from...	Then...
	touch screen	see 'Set the method and sequence' on page 126 .
	HTA Autosampler Manager	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Method , Sequence and Run topics.
Run sample(s) - Before proceeding	'Run sample(s) - Before proceeding' on page 127	
Run sample(s) - Start a sequence	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.	
Run sample(s) - Stop a sequence	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.	
Run sample(s) - Set a priority	If supported by your CDS, use your CDS. See your CDS documentation for further details.	
Run sample(s) - Process a single sample	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.	
Perform additional syringe purge	If supported by your CDS, use your CDS. See your CDS documentation for further details.	
Modify the autosampler general settings	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Setup topics.	
Manage racks	See 'Racks management' on page 130 .	

6.3.4 Operate by CDS through HTA software connectors

To...	Then...
Set the method and sequence	Use your CDS. Check the installation and usage guide of the driver.
Run sample(s) - Before proceeding	Check the installation and usage guide of the driver.
Run sample(s) - Start a sequence	Use your CDS. Check the installation and usage guide of the driver.

To...	Then...	
Run sample(s) - Stop a sequence	If supported by your CDS, use your CDS. See your CDS documentation for further details. Check the installation and usage guide of the driver.	
Run sample(s) - Process a single sample	If supported by your CDS, use your CDS. Check the installation and usage guide of the driver.	
Perform additional syringe purge	If supported by your CDS, use your CDS. See your CDS documentation for further details. Check the installation and usage guide of the driver.	
Modify the autosampler general settings	If...	Then...
	supported by your CDS	Use your CDS.
	not supported by your CDS	see ' Modify the autosampler general settings ' on page 128.
	Check the installation and usage guide of the driver.	
Manage racks	See ' Racks management ' on page 130.	


6.3.5 Operate by CDS through third party integration

To...	Then...	
Set the method and sequence	Use your third party software. See your third party software documentation for further details.	
Run sample(s) - Before proceeding	See ' Run sample(s) - Before proceeding ' on page 127.	
Run sample(s) - Start a sequence	Use your third party software. See your third party software documentation for further details.	
Run sample(s) - Stop a sequence	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Run sample(s) - Process a single sample	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Perform additional syringe purge	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Modify the autosampler general settings	If...	Then...
	supported by your third party software	Use your third party software. See your third party software documentation for further details.
	not supported by your third party software	see ' Modify the autosampler general settings ' on page 128.
Manage racks	See ' Racks management ' on page 130.	

6.4 Set the method and sequence

6.4.1 Modify a method

 *Run menu > Settings > Methods*

1. Select the desired method and then **OK** to confirm.
2. Select **General** and, if necessary, modify the parameters. Press **Exit**.
3. If desired, modify the remaining parameters.
4. Select **Exit > Save > Continue**: the **Settings** menu appears.
5. To proceed, see '[Set a sequence](#)' on page  126.

6.4.2 Set a sequence

 *Run menu > Settings > Sequences*

1. Select the first step and **OK** to open it.
2. Set the step parameters as desired.
3. Select **Exit** to exit from the first step menu.
4. To add a new sequence step, select **Settings > Sequences > Tools > OK > Add new > Continue > Exit** and set the parameters of the new step as desired.
5. Select **Exit > Save > Continue**: the **Settings** menu appears.

6.4.3 Restore the default parameters of a method

 *Run menu > Settings > Methods > Tools*

1. Select **Restore Default** and select the method to be modified.
2. Select **Continue**.
3. Select **Load** to confirm.
4. Select **Continue**.

6.4.4 Copy the parameters of a method to another method

 *Run menu > Settings > Methods > Tools*

1. Select **Copy Methods**.
2. Select the source method and the destination method.
3. Select **Copy** to confirm.
4. Select **Continue**.

6.5 Start the Run

6.5.1 Run sample(s) - Before proceeding


Before starting the run:

- Ensure to have correctly prepared the sample vials and the sample list in the analyzer software.
- Ensure methods and, if desired, sequences, have already been set.
- Verify that the parameter values of the desired method are correct

To modify a method or set a sequence, see ['Set the method and sequence' on page 126](#).

6.5.2 Run sample(s) - Start a sequence

This procedure is useful when there are many sample vials to be processed.


 *Run menu > Load Sample Tray*

1. Insert the sample vials in the tray.
2. Select **Close**.
3. Start the sample list on the analyzer software.
4. In the **Run** menu select **Start**: the autosampler starts the sequence and the summary page appears.
5. At the end of the run, select **Continue**: the **Run** menu appears.

6.5.3 Run sample(s) - Process a single sample

This procedure is useful when there is only one sample vial to be processed or to perform some test.


Note: *This procedure is applicable only if the parameter **User Interface** has been set to **Standard**.*

 *Run menu > Load Sample Tray*

1. Insert the sample vial in the tray.
2. Start the sample list on the analyzer software.
3. In the Standard version of the Run menu select **Single Injection**.
4. Set the parameters as desired.
5. Select **Start**: the autosampler starts the single injection and the summary page appears.

6.6 Run management

6.6.1 Perform additional syringe purge

 *Run menu > Syringe Purge*


1. Set **Flush time** and **Purge temperature**.
2. Set the desired parameters and select **Start**: the autosampler performs the defined additional purge.

6.6.2 Run sample(s) - Stop a sequence

1. Select **Stop**: the autosampler pauses.
2. Do what follows:

If you want...	Then...
to put the autosampler into standby	select Abort : the autosampler washes the syringe and goes into standby.
to stop the autosampler	select Terminate Current : the autosampler stops when the operations on the current sample vial end.
to continue the sequence	select Resume .

6.6.3 Modify the autosampler general settings

 *Run menu > Settings > Setup > Configuration*

Select the desired parameters and modify them. For details, see '[Configuration](#)' on page ²⁵⁷.

6.7 Sequences management

6.7.1 Add a sequence step copying from another one

 *Run menu > Settings > Sequences > Tools*


1. Select **Copy as new** and select the step to copy from.
2. Select **Copy**.
3. Select **Continue**.

6.7.2 Reorder sequence step

 *Run menu > Settings > Sequences > Tools*

1. Select **Reorder**.
2. Select the step to move and select **OK**.
3. Select the target position where you want to move the step and select **Move**.
4. Select **Continue**.

6.7.3 Remove sequence step

 *Run menu > Settings > Sequences > Tools*

1. Select **Remove Step**.
2. Select the step to remove and select **OK**.
3. Select **Continue**.

6.8 Racks management

6.8.1 Set the characteristics of the rack to be used

This procedure adds, to the list of the available racks to be used, a rack with a different vial type than the previous one.

 *Run menu > Settings > Setup > Alignment > Tray Install*

1. Select **Add New** and select **OK**.
2. Select the vial type, the vial depth and press **Continue**.
3. The new list item appears.
4. Select **Exit** > **Exit**.

6.8.2 Choose the type of rack to be used

 *Run menu > Settings > Tray Type*

1. Select the desired type of rack and select **OK** > **Mount** to confirm: the tray opens.
2. Follow the screen instructions to place the required sample vial, then select **Continue**: the vial locator touches the sample vial.
3. Select **Continue**.

6.8.3 Remove a rack from the list

 *Run menu > Settings > Setup > Alignment > Tray Install*

1. Select **Remove** and select **OK**.
2. Select the tray to be removed and select **Remove**.
3. Select **Continue** > **Exit** > **Exit**.

6.9 Syringe management


6.9.1 Workflow for replacing the syringe glass and plunger

For replacing the syringe, follow these procedures in the stated order:


1. ['Remove the syringe warmer' on page 131](#).
2. ['Remove the syringe plunger' on page 132](#).
3. ['Install the syringe plunger' on page 133](#).
4. ['Remove the syringe glass' on page 133](#).
5. ['Inspect the syringe before installation' on page 134](#).
6. ['Install the syringe glass' on page 135](#).
7. ['Install the syringe warmer' on page 136](#).

Each procedure is described in detail below.


6.9.2 Remove the syringe warmer

 *Run menu > Settings > Setup > Syringe Installation*


1. Raise the sliding lid and unscrew the two finger nuts.
2. Take the syringe pointer, insert its edge in the plunger locker and pull it out.
3. Pull the syringe to disconnect the gas connection and the electrical connection from the sledge.



4. Carefully remove the syringe needle from the vial locator.





5. Remove the two screws which fix the syringe warmer top cover and remove the top cover.



6.9.3 Remove the syringe plunger

1. Unscrew the syringe plunger top part and remove the syringe plunger from the top cover





6.9.4 Install the syringe plunger

1. Fix the new syringe plunger in its position, taking care to insert the PTFE washer in the proper position, and fasten the syringe plunger top part firmly.
2. Fix the syringe warmer top cover at the syringe warmer, taking care not to overtight the two screws. In case it will be necessary to change also the syringe glass, skip this point and follow the next chapter ['Remove the syringe glass' on page 133](#).


6.9.5 Remove the syringe glass

Syringes should be replaced according to sample throughput and sample dirtiness/matrix. To replace the syringe, it is necessary to disassemble and reassemble the syringe warmer assembly that contains the syringe.

1. Pull the gas connector to remove the syringe glass, the spacer and the aluminum spacer from the syringe warmer.



2. Holding the gas connector, pull the syringe glass to remove it.
Note that on the gas connector there are two o-rings to avoid leaking, so it is necessary to apply some force to make them to move along the syringe glass.




3. Remove the aluminium spacer from the syringe glass.



6.9.6 Inspect the syringe before installation

1. Roll the syringe over a clean flat surface and observe the needle tip:

If the needle tip...	Then...
does not move in a circle	Proceed with step 2.
moves in a circle	<div>1. Straighten the needle by bending it carefully near the connection point to the barrel.  2. Roll the syringe again to ensure the needle is straight. Otherwise, straighten it again and repeat.</div>

2. Observe the needle under a 10X magnifier:

If ridges are...	Then...
absent	Proceed with step 3.

If ridges are...	Then...
present	<ol style="list-style-type: none"> 1. Rub carefully the needle with a fine emery paper to remove all ridges. NOTICE: Do not damage the needle tip. 2. Proceed with step 3.


3. Check the syringe plunger by sliding it up and down a few times.
4. If the syringe plunger sticks, remove and clean it using a suitable solvent. Then replace it and repeat step 3.
5. Check that the needle tip is suitable for the septa to pierce.

6.9.7 Install the syringe glass


1. Insert the new syringe glass in the aluminium spacer, taking care to keep the spacer as shown in the following picture.




2. Insert the gas connector on the syringe glass until it touches the spacer, and insert the syringe glass into the syringe warmer taking care to let the needle tip to pass through the hole on syringe warmer.



NOTICE: In case the two insulating gaskets have been removed, take care to insert them along the syringe.



NOTICE: If you use Merlin Microseal™ septum, then take care to correctly orient the side hole of the needle. See Merlin Microseal™ documentation for further details.




3. Fix the syringe warmer top cover at the syringe warmer, taking care not to overtight the two screw.

6.9.8 Install the syringe warmer


1. When the screen displays the **Remove and install new syringe** page, do what follows:

If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the syringe installation without moving the turret	proceed to point 2.


2. Raise the sliding lid and unscrew the two finger nuts.




3. Take the syringe pointer, insert its edge in the plunger locker and pull it out.




4. Take the syringe.
5. Push the vial locator up and insert the syringe needle in the central hole of the vial locator.



6. Insert the syringe body in the syringe holder and the plunger in the syringe plunger holder.




7. Lower the syringe plunger until it rests on the syringe plunger holder.



8. Screw the two finger nuts.

9. Insert the plunger locker with the knurled side upward in the syringe plunger holder.



10. Lower the sliding lid.

11. Select **Continue**: the **Select syringe Volume** page appears.

12. Choose the syringe volume and select **Save**: the automatic alignment procedure starts.


Note: The automatic conditioning procedure will start, with a count-down which shows an estimation of the time required. This procedure must be done if:


- the syringe glass has been changed.
- the syringe plunger has been changed.
- advised from service support.

In case of installation of a brand new instrument, it is not necessary to perform it.


6.9.9 Align the vial locator to the needle tip

1. Remove the syringe warmer. For details, see ['Remove the syringe warmer' on page 131](#).
2. Unscrew the safety lock and the needle higher regulator and let the vial locator to rest on the frontal panel.







3. Insert the syringe pointer as indicated in the following pictures:



4. Lift up the vial locator until the syringe pointer is completely against the sledge.



5. Using a flat screwdriver, screw the safety lock and the needle height regulator to block the vial locator in position.
6. Remove the syringe pointer.




7. Install the syringe warmer.

6.10 Syringe kit management


6.10.1 Remove the syringe kit

1. Remove the syringe warmer. See ['Remove the syringe warmer' on page 131](#).
2. Remove the syringe plunger. See ['Remove the syringe plunger' on page 132](#).
3. Remove the syringe glass. See ['Remove the syringe glass' on page 133](#).
4. Remove the two spacers inside the syringe warmer.



6.10.2 Install the syringe kit

1. Install the two spacers inside the syringe warmer.



2. Install the syringe glass. See ['Install the syringe glass' on page 135](#).
3. Install the syringe plunger. See ['Install the syringe plunger' on page 133](#).
4. Install the syringe warmer. See ['Install the syringe warmer' on page 136](#).

6.11 Obtain support in case of autosampler problems

This procedure allows to send an auto-generated email with useful technical information about autosampler.


Note: *This procedure is applicable only if the parameter **Log** has been set to **Enable** or **Enable w/o warnings**.*

 *Run menu*

1. Select **Generate QR code**. The Generate QR code wizard appears on the PC screen.
2. Follow the wizard instructions.
3. Select **Continue**.

6.12 Change the analyzer to be used

✱ *Run menu > Settings > Setup > Alignment > Analyzer Default*

1. If the name of the analyzer to be used is already listed, select it. Otherwise, select **Generic**.
2. Select **OK** to confirm.
3. Select the number of injectors to be used and then **OK** to confirm and **Load**.
4. Select **Continue**.
5. Proceed with the ['Align injectors' procedure on the following section.](#) 

6.13 Realign injectors

✱ *Run menu > Settings > Setup > Alignment*

1. To align the autosampler to the front injector, select **Front Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the front injector adapter.



CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

2. Do what follows:

If you want to...	Then...
confirm the injector alignment	select Confirm .
align the autosampler to the injector manually	<ol style="list-style-type: none"> 1. Select Do Manual. 2. Select the arrows to move the vial locator to the correct position (Touch Sensor must be ON). 3. Select Store to confirm.

3. Set the desired injection depth: the autosampler tries to insert the needle in the injector.
4. Select **Save** to confirm.
5. Select the desired injection speed (from **Very Low** to **Very High**) and select **Save**: the turret returns to its original position and the **Alignment** menu appears.
6. If the rear injector is present, select **Rear Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the rear injector adapter.







CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

7. Repeat from step 2 to 7.

Vial locator aligned to the injector adapter

To be aligned, the vial locator must lean against the top of the injector adapter completely.





6.14 Reset the position of the vial

 *Run menu > Settings > Setup > Alignment*

1. Select **Touch & Plunger Zero**: the turret moves and the tray opens.
2. Insert a crimped sample vial in position A4 of the rack and in C position of the oven (if using 6ml or 10 ml vial, ensure to insert first the proper spacer in C position). In case syringe integrity tool is available, insert it in its position.
3. Select **Continue**: the vial locator touches the syringe integrity tool, the sample vial in the rack, the sample vial in the oven, the cover and the injector adapter to memorize their positions.
4. Remove the vial from oven and press **Continue**.
5. At the end of the procedure, when the **Setup menu** appears, select **Exit > Exit**.

6.15 Get out of Headspace mode


✱ *Run menu > Change Mode*

1. Select the mode you are going to use and select **Continue**.


6.15.1 Going to Liquid mode

✱ *Run menu > Change Mode > Liquid*


1. Lift the sliding lid.
2. Remove the vial locator.



3. Select **Continue**.
4. Take the syringe pointer, insert its edge in the plunger locker and pull it out.



5. Unscrew the two finger nuts.
6. Remove the syringe warmer assembly.



- 7. Select **Continue**.
- 8. Follow the instructions of the ['Get into Liquid mode' chapter on page 153](#).


6.15.2 Going to SPME mode


✦ *Run menu > Change Mode > SPME*

1.



If the autosampler...	Then...
asks Which fiber type are you going to use	choose if Sealed tip or Standard fiber.
asks if the installed vial locator o-ring is blue	select Yes or No and proceed to point 2.

2.


If the autosampler...	Then...
asks to Remove the vial locator	<div>1. Remove the two finger screws to remove vial locator.</div> <div></div>

If the autosampler...	Then...
	 <p data-bbox="727 813 963 842">2. Select Continue.</p>
asks to Remove the syringe warmer assembly	proceed to point 3.


3. Raise the sliding lid and unscrew the two finger nuts.
4. Take the syringe pointer, insert its edge in the plunger locker and pull it out.



5. Pull the syringe to disconnect the gas connection and the electrical connection from the sledge.



6. Carefully remove the syringe needle from the vial locator.




7. Select **Continue**.
8. Follow the instructions of the ['Get into SPME mode' chapter on page 93](#).

7. Use in Liquid mode


7.1 Get into Liquid mode

To get into Liquid mode, you need first to get out of the mode you are currently using: therefore, please follow the '[Get out of Headspace mode](#)' on page [148](#) or '[Get out of SPME mode](#)' on page [237](#).


1. Select the Syringe volume and select **Save**.
2. If necessary, raise the sliding lid and remove the two finger nuts.
3. Install the syringe holder.



4. Fix completely the two finger nuts.
5. If necessary, take the syringe pointer, insert its edge in the plunger locker and pull it out.
6. Insert the plunger extension in the plunger holder, taking care to keep frontwards the side with the hole for the syringe plunger.



7. Insert the plunger locker.




7.1.1 Install the syringe



1. When the screen displays the **Remove and install new syringe** page, do what follows:

If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the syringe installation without moving the turret	proceed to point 2.


2. Open the syringe locker.



- 3. Take the syringe.
- 4. Insert the syringe body in the syringe holder and the plunger in the plunger extension.




5. Close the syringe locker.




6. Select **Continue**: the **Select syringe Volume [μl]** page appears.

7. Choose the syringe volume and select **Save**: the automatic alignment procedure starts.


8. Wait for the automatic positing.
9. Install the vial locator for Liquid technique (green o-ring: see ['Vial locators' on page 55](#)), taking care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars.




10. Fasten completely the two fixing screws, taking care to have the small bars completely inside the vial locator bars.



11. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the syringe.



12. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.




13. Lower the sliding lid.

14. Press **Continue**.


7.1.2 Sample rack installation

Insert the sample standard rack in the autosampler tray

1. Remove the sample vials and the 42 positions sample rack.
2. Insert the front side of sample rack.





3. Insert the other side of the sample rack and push it so that the label with the letters is parallel to the base of the tray.



4. If this is the first installation, select **Continue**: the autosampler tray closes and the screen displays the autosampler model and the firmware version.

Rack correct positioning



7.1.3 Install the tray

1. At the end of the alignment procedure, select **Continue**.
2. Select **OK**.
3. If necessary, modify the vial type, then select **Continue**.
4. Set the depth below the septum where the needle must be positioned to aspirate the sample and select **Continue** > **Save** > **Continue**.

Note: If you set the sample depth to 0 mm (the minimum value), the needle tip will be just below the

septum. If you set the sample depth to 32 mm (the maximum value), the needle tip will be close to the vial bottom.

Note: *The autosampler shows automatically the vial depth used the last time that the autosampler has been used in Liquid mode.*

7.1.4 Set the position of the vials - The Touch and Plunger Zero procedure

1. Insert a sample vial with cap in position A6 of the standard rack, a black spacer and a solvent vial in position C of the oven and ensure the waste vial is in its proper position.
2. Select **Continue**: the vial locator touches the waste vial, the sample vial, the solvent vial and the injector adapter to memorize their positions.

7.2 Prepare the solvent vials


7.2.1 Prepare the solvent vials

✱ *Run menu > Load Sample Tray*


1. Take a solvent vial and fill it with the desired solvent.
2. Close the solvent vial with the cap.
3. If necessary, insert the black spacer in the desired oven position.
4. Insert the solvent vial on the black spacer.
5. Repeat for all the solvent vials required.

7.2.2 Prepare the solvent and waste vial caps


1. Remove the cap from the solvent/waste vial.
2. Remove the septum.
3. Place the insert on the white washer.



4. Insert the low maintenance septa into the cap.



5. Insert the vial in the desired location.



Note: Take care not to fill completely the solvent vial (solvent in the vial must not be in contact with any part of the septa).

7.2.3 Reuse a solvent vial

1. Extract the solvent vial from its position in the autosampler with care.
2. Open the solvent vial carefully and empty it.
3. Wash the vial and the cap.
4. If necessary, change the septum.
5. Fill the vial with the desired solvent.
6. Close the solvent vial with the cap.
7. Insert the solvent vial in the desired position in the solvent and waste vial location.

7.2.4 Reuse the waste vial

1. Open the waste vial with care and empty it according to your laboratory instructions.
2. Wash the waste vial and its cap.
3. If necessary, change the septum.
4. Close the waste vial with the cap.
5. Insert the waste vial in the waste position.

7.3 Prepare the sample vials

7.3.1 Basic information on sample vials

For preparing the sample vials correctly, consider what follows:


- The sample volume must be maximum half of the vial nominal volume (for example, 1 ml for a 2 ml vial and 50 μ l for a 100 μ l vial). The air contained in the vial avoids vacuum formation during sample aspiration.
- Large samples should be divided into different vials to obtain more reliable data.
- Once the sample vial has been closed, the septum must be flat and not twisted.
- Contaminants from the previous sample or solvent washes will have a bigger impact on a low sample volume.
- Reassess the method for preparing the sample vials if you change the vial supplier. Different vial manufacturing practices can cause variances in analytical data.
- The autosampler can also be used with clear or amber glass sample vials with snap, crimp or screw cap. Amber glass vials can be used for light-sensitive samples.

7.3.2 Close a sample vial with a crimp cap

You must have a crimper available. Fit the septum on the cap and then:

1. Clean the inside surfaces of the crimper jaws.
2. Place the crimp cap over the top of the vial.
3. Hold the vial and place the crimper on top of it. Squeeze the handle until it touches the adjuster screw.

Note: Following picture shows a not correct crimping.



4. Check each vial for proper crimping:
 - Check the aluminum cap is not folded or wrinkled.
 - The cap should be tight and not easy to move or rotate.
 - If either of the above problems occur crimp the cap again. Regulate the crimper to adjusting the tightness if needed.
 - Check that each cap has a flat septum and that it is centered over the top of the vial:

If...	Then...
the septum is not flat	remove the cap, turn the crimper adjusting screw and try again.
the cap is not centered	remove the cap and make sure the new cap is flat on the top of the vial before the crimp operation.

7.4 Operating Modes

7.4.1 Operate by autosampler touch screen

To...	Then...	
Set the method and sequence	Do what follows:	
	If...	Then...
	all samples are processed in the same way	you don't need to change method and sequence at each run
	you are required to process samples in different way	see ' Set the method and sequence ' on page 168.
Run sample(s) - Before proceeding	See ' Run sample(s) - Before proceeding ' on page 169.	
Run sample(s) - Start a sequence	See ' Run sample(s) - Start a sequence ' on page 169.	
Run sample(s) - Stop a sequence	See ' Run sample(s) - Stop a sequence ' on page 170.	
Run sample(s) - Set a priority	See ' Run sample(s) - Set a priority ' on page 170.	
Run sample(s) - Process a single sample	If...	Then...
	all samples are processed in the same way	see ' Run sample(s) - Start a sequence ' on page 169.
	you are required to process samples in different way	see ' Run sample(s) - Process a single sample ' on page 169.
Perform additional syringe washes	See ' Perform additional syringe washes ' on page 170.	
Modify the autosampler general settings	See ' Modify the autosampler general settings ' on page 170.	
Manage racks	See ' Racks management ' on page 172.	

7.4.2 Operate by HTA Autosampler Manager

To...	Then...
Set the method and sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Method , Sequence and Run topics.
Run sample(s) - Before proceeding	See 'Run sample(s) - Before proceeding' on page 169 .
Run sample(s) - Start a sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Run topic.
Run sample(s) - Stop a sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Run topic.
Run sample(s) - Set a priority	Not applicable.
Run sample(s) - Process a single sample	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Auxiliary > Single Injection topic.
Perform additional syringe washes	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Auxiliary > Wash topic.
Modify the autosampler general settings	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Setup topic.
Manage racks	See 'Racks management' on page 172 .

7.4.3 Operate by CDS through HTA software integration

To...	Then...	
Set the method and sequence	If you want to work from...	Then...
	touch screen	see 'Set the method and sequence' on page 168 .
	HTA Autosampler Manager	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Method , Sequence and Run topics.
Run sample(s) - Before proceeding	See 'Run sample(s) - Before proceeding' on page 169 .	
Run sample(s) - Start a sequence	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.	
Run sample(s) - Stop a sequence	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.	
Run sample(s) - Set a priority	If supported by your CDS, use your CDS. See your CDS documentation for further details.	

To...	Then...
Run sample(s) - Process a single sample	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.
Perform additional syringe washes	If supported by your CDS, use your CDS. See your CDS documentation for further details.
Modify the autosampler general settings	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Setup topics.
Manage racks	See 'Racks management' on page 172 .

7.4.4 Operate by CDS through HTA software connectors

To...	Then...	
Set the method and sequence	Use your CDS. Check the installation and usage guide of the driver.	
Run sample(s) - Before proceeding	Check the installation and usage guide of the driver.	
Run sample(s) - Start a sequence	Use your CDS. Check the installation and usage guide of the driver.	
Run sample(s) - Stop a sequence	If supported by your CDS, use your CDS. See your CDS documentation for further details. Check the installation and usage guide of the driver.	
Run sample(s) - Set a priority	If supported by your CDS, use your CDS. See your CDS documentation for further details. Check the installation and usage guide of the driver.	
Run sample(s) - Process a single sample	If supported by your CDS, use your CDS. Check the installation and usage guide of the driver.	
Perform additional syringe washes	If supported by your CDS, use your CDS. See your CDS documentation for further details. Check the installation and usage guide of the driver.	
Modify the autosampler general settings	If...	Then...
	supported by your CDS	Use your CDS.
	not supported by your CDS	see ' Modify the autosampler general settings ' on page 170.
	Check the installation and usage guide of the driver.	
Manage racks	See ' Racks management ' on page 172.	


7.4.5 Operate by CDS through third party integration

To...	Then...	
Set the method and sequence	Use your third party software. See your third party software documentation for further details.	
Run sample(s) - Before proceeding	See 'Run sample(s) - Before proceeding' on page 169 .	
Run sample(s) - Start a sequence	Use your third party software. See your third party software documentation for further details.	
Run sample(s) - Stop a sequence	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Run sample(s) - Set a priority	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Run sample(s) - Process a single sample	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Perform additional syringe washes	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Modify the autosampler general settings	If...	Then...
	supported by your third party software	Use your third party software. See your third party software documentation for further details.
	not supported by your third party software	see 'Modify the autosampler general settings' on page 170 .
Manage racks	See 'Racks management' on page 172 .	

7.5 Set the method and sequence

7.5.1 Modify a method

 *Run menu > Settings > Methods*

1. Select the desired method and then **OK** to confirm.
2. Select **General** and, if necessary, modify the parameters. Press **Exit**.
3. If desired, modify the remaining parameters.
4. Select **Exit > Save > Continue**: the **Settings** menu appears.
5. To proceed, see '[Set a sequence](#)' on page  168.

7.5.2 Set a sequence

 *Run menu > Settings > Sequences*

1. Select the first step and **OK** to open it.
2. Set the step parameters as desired.
3. Select **Exit** to exit from the first step menu.
4. To add a new sequence step, select **Settings > Sequences > Tools > OK > Add new > Continue > Exit** and set the parameters of the new step as desired.
5. Select **Exit > Save > Continue**: the **Settings** menu appears.

7.5.3 Restore the default parameters of a method

 *Run menu > Settings > Methods > Tools*

1. Select **Restore Default** and select the method to be modified.
2. Select **Continue**.
3. Select **Load** to confirm.
4. Select **Continue**.

7.5.4 Copy the parameters of a method to another method

 *Run menu > Settings > Methods > Tools*

1. Select **Copy Methods**.
2. Select the source method and the destination method.
3. Select **Copy** to confirm.
4. Select **Continue**.

7.6 Start the Run

7.6.1 Run sample(s) - Before proceeding


Before starting the run:

- Ensure to have correctly prepared the sample vials and the sample list in the analyzer software.
- Ensure methods and, if desired, sequences, have already been set.
- Verify that the parameter values of the desired method are correct

To modify a method or set a sequence, see ['Set the method and sequence' on page 168](#).

7.6.2 Run sample(s) - Start a sequence

This procedure is useful when there are many sample vials to be processed.


 *Run menu > Load Sample Tray*

1. Insert the sample vials in the tray.
2. Select **Close**.
3. Start the sample list on the analyzer software.
4. In the **Run** menu select **Start**: the autosampler starts the sequence and the summary page appears.
5. At the end of the run, select **Continue**: the **Run** menu appears.

7.6.3 Run sample(s) - Process a single sample

This procedure is useful when there is only one sample vial to be processed or to perform some test.


Note: *This procedure is applicable only if the parameter **User Interface** has been set to **Standard**.*

 *Run menu > Load Sample Tray*

1. Insert the sample vial in the tray.
2. Start the sample list on the analyzer software.
3. In the Standard version of the Run menu select **Single Injection**.
4. Set the parameters as desired.
5. Select **Start**: the autosampler starts the single injection and the summary page appears.

7.7 Run management

7.7.1 Perform additional syringe wash

 *Run menu > Syringe Wash*

1. Set **Solvent Vial** and **Wash Volume [µl]**.
2. Set the number of desired washes and select **Start Wash**: the autosampler performs the defined additional washes.

7.7.2 Run sample(s) - Stop a sequence

1. Select **Stop**: the autosampler pauses.
2. Do what follows:


If you want...	Then...
to put the autosampler into standby	select Abort : the autosampler washes the syringe and goes into standby.
to stop the autosampler	select Terminate Current : the autosampler stops when the operations on the current sample vial end.
to continue the sequence	select Resume .

7.7.3 Run sample(s) - Set a priority

This procedure allows to process a sample vial immediately, even if an automatic run is already in process and other vials should be processed first. Using this procedure you give priority to a specific sample vial that will be processed first.

1. In the summary page select **Menu**: the **Automatic Run** menu appears.
2. Select **Priority**: the sample tray opens when the operations on the current sample vial end.
3. Insert the priority vial in the desired position.
4. Select **Continue**.
5. Set the position of the sample vial to be processed with priority and set the other parameters as desired. Select **Start Priority**: the autosampler starts the operations on the priority vial and then automatically resume the original sample list.

7.7.4 Modify the autosampler general settings

 *Run menu > Settings > Setup > Configuration*

Select the desired parameters and modify them. For details, see ['Configuration' on page 286](#).

7.8 Sequences management

7.8.1 Add a sequence step copying from another one

 *Run menu > Settings > Sequences > Tools*


1. Select **Copy as new** and select the step to copy from.
2. Select **Copy**.
3. Select **Continue**.

7.8.2 Reorder sequence step

 *Run menu > Settings > Sequences > Tools*

1. Select **Reorder**.
2. Select the step to move and select **OK**.
3. Select the target position where you want to move the step and select **Move**.
4. Select **Continue**.

7.8.3 Remove sequence step

 *Run menu > Settings > Sequences > Tools*

1. Select **Remove Step**.
2. Select the step to remove and select **OK**.
3. Select **Continue**.

7.9 Racks management

7.9.1 Set the characteristics of the rack to be used

This procedure adds, to the list of the available racks to be used, a rack with a different draw depth than the previous one.

 *Run menu > Settings > Setup > Alignment > Tray Install*

1. Select **Add New** and select **OK**.
2. Select **Continue** > **Continue**.
3. Set the draw depth and select **Save** > **Continue**: the new list item appears.
4. Select **Exit** > **Exit**.

7.9.2 Choose the type of rack to be used

 *Run menu > Settings > Tray Type*

1. Select the desired type of rack and select **OK** > **Mount** to confirm: the tray opens.
2. Follow the screen instructions to place the required sample vial, then select **Continue**: the vial locator touches the sample vial.
3. Select **Continue**.

7.9.3 Remove a rack from the list

 *Run menu > Settings > Setup > Alignment > Tray Install*

1. Select **Remove** and select **OK**.
2. Select the tray to be removed and select **Remove**.
3. Select **Continue** > **Exit** > **Exit**.

7.10 Syringe management

7.10.1 Workflow for replacing the syringe

For replacing the syringe, follow these procedures in the stated order:


1. ['Remove the syringe' on page 173](#).
2. Check if the syringe holder is suitable for the new syringe. For details see ['Syringes for Liquid technique' on page 403](#).
3. If necessary, replace the syringe holder.
4. ['Inspect the syringe before installation' on page 174](#).
5. ['Install the syringe' on page 175](#).

Each procedure is described in detail below.


7.10.2 Remove the syringe

✱ *Run menu > Settings > Setup > Syringe Installation*

1. Raise the sliding lid above the syringe location.
2. Open the syringe locker.





3. Turn the syringe until it comes out from the syringe holder.



- Carefully remove the syringe from the syringe holder and the syringe needle from the vial locator.


7.10.3 Prepare the syringe

This procedure is only for a syringe with o-ring or reinforcement.

If the syringe has...	Then...
an o-ring	put the syringe through the o-ring from the needle until the o-ring reaches the top of the syringe. 
a reinforcement	completely put the syringe needle through the reinforcement. 

7.10.4 Inspect the syringe before installation

- Roll the syringe over a clean flat surface and observe the needle tip:

If the needle tip...	Then...
does not move in a circle	Proceed with step 2.
moves in a circle	<ol style="list-style-type: none"> Straighten the needle by bending it carefully near the connection point to the barrel.  <ol style="list-style-type: none"> Roll the syringe again to ensure the needle is straight. Otherwise, straighten it again and repeat.

- Observe the needle under a 10X magnifier:

If ridges are...	Then...
absent	Proceed with step 3.
present	<ol style="list-style-type: none"> Rub carefully the needle with a fine emery paper to remove all ridges. NOTICE: Do not damage the needle tip. Proceed with step 3.


- Check the syringe plunger by sliding it up and down a few times.
- If the syringe plunger sticks, remove and clean it using a suitable solvent. Then replace it and repeat step 3.
- Check that the needle tip is suitable for the septa to pierce.

7.10.5 Install the syringe


1. When the screen displays the **Remove and install new syringe** page, do what follows:

If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the syringe installation without moving the turret	proceed to point 2.



2. Open the syringe locker.



3. Take the syringe.
 4. Push the vial locator up and insert the syringe needle in the central hole of the vial locator.




5. Insert the syringe body in the syringe holder and the plunger in the in the plunger extension.




6. Close the syringe locker.
7. Select **Continue**: the **Select syringe Volume [μl]** page appears.
8. Choose the syringe volume and select **Save**: the automatic alignment procedure starts.
9. Wait for the automatic positing.
10. Lower the sliding lid.


7.10.6 Align the vial locator to the needle tip

1. Remove the syringe and the syringe holder. For details, see ['Remove the syringe' on page 173](#).
2. Unscrew the safety lock and the needle higher regulator and let the vial locator to rest on the frontal panel.




3. Insert the syringe pointer as indicated in the following pictures:






4. Lift up the vial locator until the syringe pointer is completely against the sledge.



5. Using a flat screwdriver, screw the safety lock and the needle height regulator to block the vial locator in position.
6. Remove the syringe pointer.



7. Install the syringe holder and the syringe.

7.11 Obtain support in case of autosampler problems

This procedure allows to send an auto-generated email with useful technical information about autosampler.


Note: *This procedure is applicable only if the parameter **Log** has been set to **Enable** or **Enable w/o warnings**.*

 *Run menu*

1. Select **Generate QR code**. The Generate QR code wizard appears on the PC screen.
2. Follow the wizard instructions.
3. Select **Continue**.

7.12 Change the analyzer to be used

✱ *Run menu > Settings > Setup > Alignment > Analyzer Default*

1. If the name of the analyzer to be used is already listed, select it. Otherwise, select **Generic**.
2. Select **OK** to confirm.
3. Select the number of injectors to be used and then **OK** to confirm and **Load**.
4. Select **Continue**.
5. Proceed with the ['Align injectors' procedure on the following section.](#) 

7.13 Realign injectors

✱ *Run menu > Settings > Setup > Alignment*

1. To align the autosampler to the front injector, select **Front Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the front injector adapter.



CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

2. Do what follows:

If you want to...	Then...
confirm the injector alignment	select Confirm .
align the autosampler to the injector manually	<ol style="list-style-type: none"> 1. Select Do Manual. 2. Select the arrows to move the vial locator to the correct position (Touch Sensor must be ON). 3. Select Store to confirm.

3. Set the desired injection depth: the autosampler tries to insert the needle in the injector.
4. Select **Save** to confirm.
5. Select the desired injection speed (from **Very Low** to **Very High**) and select **Save**: the turret returns to its original position and the **Alignment** menu appears.
6. If the rear injector is present, select **Rear Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the rear injector adapter.





CAUTION Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

7. Repeat from step 2 to 7.

Vial locator aligned to the injector adapter

To be aligned, the vial locator must lean against the top of the injector adapter completely.





7.14 Reset the position of the vial


 *Run menu > Settings > Setup > Alignment*

1. Select **Touch & Plunger Zero**: the turret moves and the tray opens.
2. Insert a sample vial with cap in position A6 of the standard rack, a black spacer and a solvent vial in position C of the oven and ensure the waste vial is in its proper position.
3. Select **Continue**: the vial locator touches the waste vial, the sample vial, the solvent vial and the injector adapter to memorize their positions.
4. At the end of the procedure, when the **Setup menu** appears, select **Exit > Exit**.


7.15 Get out of Liquid mode


✱ *Run menu > Change Mode*

1. Select the mode you are going to use and select **Continue**.
2. Wait for the automatic positioning.
3. Remove the two fixing screws to remove the vial locator.




4. Select **Continue**.
5. Wait for the automatic positioning, open the locker and remove the syringe.







6. Take the syringe pointer, insert its edge in the plunger locker and pull it out.



7. Remove the plunger extension.



8. Unscrew the two fixing nuts and remove the syringe holder.



9. Select **Continue**.

10. Follow the instructions of the 'Get in' chapter of the mode you are going to use.

8. Use in SPME mode


8.1 Get into SPME mode (fiber)



To get into SPME mode, you need first to get out of the mode you are currently using: therefore, please follow the ['Get out of Headspace mode' on page 148](#) or ['Get out of Liquid mode' on page 184](#).


- 1. Select the **Type** of fiber / sealed tip you are going to use.
- 2. When the screen displays the **Install new fiber holder or Sealed tip holder** page, do what follows:

If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the fiber holder installation without moving the turret	proceed to point 3.


- 3. If necessary, raise the sliding lid and unscrew the two finger nuts.
- 4. If necessary, take the syringe pointer, insert its edge in the plunger locker and pull it out.
- 5. Wait for the automatic positioning.

If the vial locator...	Then...
is already mounted	<div>1. Push the vial locator up and insert the fiber in the central hole of the vial locator.</div> <div></div>


If the vial locator...	Then...
	<div data-bbox="496 286 890 835"></div> <div data-bbox="416 880 735 913">2. Install the fiber holder.</div> <div data-bbox="496 952 890 1512"></div>
is not mounted yet	install the fiber holder.

If the vial locator...	Then...
	

6. Fix completely the two finger nuts.








7. Insert the plunger locker.





8. Press **Continue**.

9.

If...	Then...
<p>the autosampler asks to mount the vial locator and the standard fiber is installed</p>	<ol style="list-style-type: none"> <li data-bbox="430 353 1351 492">1. Install the vial locator for Headspace and SPME standard fibers technique (blue o-ring: see 'Vial locators' on page 55), taking care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars.  <ol style="list-style-type: none"> <li data-bbox="430 1070 962 1104">2. Fasten completely the two fixing screws. <li data-bbox="430 1111 1351 1182">3. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the fiber.  <ol style="list-style-type: none"> <li data-bbox="430 1720 1303 1825">4. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.

If...	Then...
	
	5. Select Continue .
the autosampler executes the Automatic Fiber / Sealed tip Zero procedure	proceed to point 10.
10. Wait the end of the of the Automatic Fiber / Sealed tip Zero procedure.	
11.	
If the autosampler asks to...	Then...
remove the sample vials and the 121 positions sample rack is installed	<p>proceed to sample rack installation:</p> <ol style="list-style-type: none">1. Remove the sample vials and the 121 positions sample rack.2. Insert the front side of sample rack.  <ol style="list-style-type: none">3. Insert the other side of the sample rack and push it so that the label with the letters is parallel to the base of the tray.  <ol style="list-style-type: none">4. Check that rack is correctly positioned.

If the autosampler asks to...	Then...
	 
remove the sample vials and the 42 positions rack is installed	remove the sample vials.

8.1.1 Set the position of the vials - The Touch and Plunger Zero procedure

1. Insert a crimped sample vial in position A4 of the rack and in C position of the oven (if using 6ml or 10 ml vial, ensure to insert first the proper spacer in C position). In case wash vial is used, insert it in its position.
2. Select **Continue**: the vial locator touches the wash vial, the sample vial in the rack, the sample vial in the oven, the cover and the injector adapter to memorize their positions.
3. Remove the vial from oven and press **Continue**.


8.2 Get into SPME mode (sealed tip)

To get into SPME mode, you need first to get out of the mode you are currently using: therefore, please follow the ['Get out of Headspace mode' on page 148](#) or ['Get out of Liquid mode' on page 184](#).

- 1. Select the **Type** of fiber / sealed tip you are going to use.
- 2. When the screen displays the **Install new fiber holder or Sealed tip holder** page, do what follows:




If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the sealed tip holder installation without moving the turret	proceed to point 3.


- 3. If necessary, raise the sliding lid and unscrew the two finger nuts.
- 4. If necessary, take the syringe pointer, insert its edge in the plunger locker and pull it out.
- 5. Insert the sealed tip spacer in the position shown in the following pictures.




6.

If the vial locator...	Then...
is already mounted	1. Push the vial locator up and insert the sealed tip in the central hole of the vial locator.


If the vial locator...	Then...
	 <p data-bbox="421 797 798 831">2. Install the sealed tip holder.</p> 
is not mounted yet	install the sealed tip holder. 

If the vial locator...	Then...
	


7. Fix completely the two finger nuts.




8. Insert the plunger extension.



9. Lower the plunger extension until it rests on the plunger holder.



10. Insert the plunger locker.








11. Press **Continue**.

12.

If...	Then...
the autosampler asks to mount the vial locator	1. Install the vial locator for SPME for 1.1 fibers / 1.5 fibers technique (yellow o-ring: see 'Vial locators' on page 55), taking care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars.

If...	Then...
	<div data-bbox="502 286 1257 779"> </div> <div data-bbox="427 831 1337 943"> <p>2. Fasten completely the two fixing screws.</p> <p>3. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the sealed tip.</p> </div> <div data-bbox="512 981 1302 1424"> </div> <div data-bbox="427 1480 1302 1585"> <p>4. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.</p> </div>

If...	Then...
	
	5. Select Continue .
the autosampler executes the Automatic Fiber / Sealed tip Zero procedure	proceed to point 10.
10. Wait the end of the of the Automatic Fiber / Sealed tip Zero procedure.	
11.	
If the autosampler asks to...	Then...
remove the sample vials and the 121 positions sample rack is installed	<p>proceed to sample rack installation:</p> <ol style="list-style-type: none">1. Remove the sample vials and the 121 positions sample rack.2. Insert the front side of sample rack.  <ol style="list-style-type: none">3. Insert the other side of the sample rack and push it so that the label with the letters is parallel to the base of the tray.  <ol style="list-style-type: none">4. Check that rack is correctly positioned.

If the autosampler asks to...	Then...
	 
remove the sample vials and the 42 positions rack is installed	remove the sample vials.

8.2.1 Set the position of the vials - The Touch and Plunger Zero procedure

1. Insert a crimped sample vial in position A4 of the rack and in C position of the oven (if using 6ml or 10 ml vial, ensure to insert first the proper spacer in C position). In case wash vial is used, insert it in its position.
2. Select **Continue**: the vial locator touches the wash vial, the sample vial in the rack, the sample vial in the oven, the cover and the injector adapter to memorize their positions.
3. Remove the vial from oven and press **Continue**.

8.3 Prepare the sample vials

8.3.1 Basic information on sample vials


- Vials must be approved for SPME techniques;
- Rounded end;
- Clear or amber glass (amber glass is suitable for light sensitive samples). For more details see ['Sample vial specifications' on page 395](#).

8.3.2 Close a sample vial with a crimp cap

You must have a crimper available. Fit the septum on the cap and then:

1. Clean the inside surfaces of the crimper jaws.
2. Place the crimp cap over the top of the vial.
3. Hold the vial and place the crimper on top of it. Squeeze the handle until it touches the adjuster screw.

Note: Following picture shows a not correct crimping.



4. Check each vial for proper crimping:
 - Check the aluminum cap is not folded or wrinkled.
 - The cap should be tight and not easy to move or rotate.
 - If either of the above problems occur crimp the cap again. Regulate the crimper to adjusting the tightness if needed.
 - Check that each cap has a flat septum and that it is centered over the top of the vial:

If...	Then...
the septum is not flat	remove the cap, turn the crimper adjusting screw and try again.
the cap is not centered	remove the cap and make sure the new cap is flat on the top of the vial before the crimp operation.

8.4 Operating Modes

8.4.1 Operate by autosampler touch screen

To...	Then...	
Set the method and sequence	Do what follows:	
	If...	Then...
	all samples are processed in the same way	you don't need to change method and sequence at each run
	you are required to process samples in different way	see ' Set the method and sequence ' on page 205.
Run sample(s) - Before proceeding	See ' Run sample(s) - Before proceeding ' on page 206.	
Run sample(s) - Start a sequence	See ' Run sample(s) - Start a sequence ' on page 206.	
Run sample(s) - Stop a sequence	See ' Run sample(s) - Stop a sequence ' on page 207.	
Run sample(s) - Process a single sample	If...	Then...
	all samples are processed in the same way	see ' Run sample(s) - Start a sequence ' on page 206.
	you are required to process samples in different way	see ' Run sample(s) - Process a single sample ' on page 206.
Perform additional fiber conditioning	See ' Perform additional fiber cleaning conditioning ' on page 207.	
Modify the autosampler general settings	See ' Modify the autosampler general settings ' on page 207.	
Manage racks	See ' Racks management ' on page 209.	

8.4.2 Operate by HTA Autosampler Manager

To...	Then...
Set the method and sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Method , Sequence and Run topics.
Run sample(s) - Before proceeding	See ' Run sample(s) - Before proceeding ' on page 206.
Run sample(s) - Start a sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Run topic.
Run sample(s) - Stop a sequence	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Run topic.
Run sample(s) - Process a single sample	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Auxiliary > Single Injection topic.
Perform additional syringe purge	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Auxiliary > Purge topic.
Modify the autosampler general settings	Open the autosampler with HTA Autosampler Manager and see

To...	Then...
	the help by pressing F1 . See in particular Setup topic.
Manage racks	See ' Racks management ' on page 209.

8.4.3 Operate by CDS through HTA software integration

To...	Then...	
Set the method and sequence	If you want to work from...	Then...
	touch screen	see ' Set the method and sequence ' on page 205.
	HTA Autosampler Manager	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Method , Sequence and Run topics.
Run sample(s) - Before proceeding	See ' Run sample(s) - Before proceeding ' on page 206.	
Run sample(s) - Start a sequence	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.	
Run sample(s) - Stop a sequence	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.	
Run sample(s) - Set a priority	If supported by your CDS, use your CDS. See your CDS documentation for further details.	
Run sample(s) - Process a single sample	Use your CDS. For further details, see your CDS documentation and HTA Autosampler Manager help.	
Perform additional syringe purge	If supported by your CDS, use your CDS. See your CDS documentation for further details.	
Modify the autosampler general settings	Open the autosampler with HTA Autosampler Manager and see the help by pressing F1 . See in particular Setup topics.	
Manage racks	See ' Racks management ' on page 209.	

8.4.4 Operate by CDS through software connectors

To...	Then...
Set the method and sequence	Use your CDS. Check the installation and usage guide of the driver.
Run sample(s) - Before proceeding	Check the installation and usage guide of the driver.
Run sample(s) - Start a sequence	Use your CDS. Check the installation and usage guide of the driver.
Run sample(s) - Stop a sequence	If supported by your CDS, use your CDS. See your CDS documentation for further details. Check the

To...	Then...	
	installation and usage guide of the driver.	
Run sample(s) - Process a single sample	If supported by your CDS, use your CDS. Check the installation and usage guide of the driver.	
Perform additional syringe purge	If supported by your CDS, use your CDS. See your CDS documentation for further details. Check the installation and usage guide of the driver.	
Modify the autosampler general settings	If...	Then...
	supported by your CDS	Use your CDS.
	not supported by your CDS	see 'Modify the autosampler general settings' on page 207 .
	Check the installation and usage guide of the driver.	
Manage racks	See 'Racks management' on page 209 .	

8.4.5 Operate by CDS through third party integration

To...	Then...	
Set the method and sequence	Use your third party software. See your third party software documentation for further details.	
Run sample(s) - Before proceeding	See 'Run sample(s) - Before proceeding' on page 206 .	
Run sample(s) - Start a sequence	Use your third party software. See your third party software documentation for further details.	
Run sample(s) - Stop a sequence	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Run sample(s) - Set a priority	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Run sample(s) - Process a single sample	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Perform additional syringe purge	If supported by your third party software, use your third party software. See your third party software documentation for further details.	
Modify the autosampler general settings	If...	Then...
	supported by your third party software	Use your third party software. See your third party software documentation for further details.
	not supported by your third party software	see 'Modify the autosampler general settings' on page 207 .
Manage racks	See 'Racks management' on page 209 .	

8.5 Set the method and sequence

8.5.1 Modify a method

 *Run menu > Settings > Methods*

1. Select the desired method and then **OK** to confirm.
2. Select **General** and, if necessary, modify the parameters. Press **Exit**.
3. If desired, modify the remaining parameters.
4. Select **Exit > Save > Continue**: the **Settings** menu appears.
5. To proceed, see '[Set a sequence](#)' on page 205.

8.5.2 Set a sequence

 *Run menu > Settings > Sequences*

1. Select the first step and **OK** to open it.
2. Set the step parameters as desired.
3. Select **Exit** to exit from the first step menu.
4. To add a new sequence step, select **Settings > Sequences > Tools > OK > Add new > Continue > Exit** and set the parameters of the new step as desired.
5. Select **Exit > Save > Continue**: the **Settings** menu appears.

8.5.3 Restore the default parameters of a method

 *Run menu > Settings > Methods > Tools*

1. Select **Restore Default** and select the method to be modified.
2. Select **Continue**.
3. Select **Load** to confirm.
4. Select **Continue**.

8.5.4 Copy the parameters of a method to another method

 *Run menu > Settings > Methods > Tools*

1. Select **Copy Methods**.
2. Select the source method and the destination method.
3. Select **Copy** to confirm.
4. Select **Continue**.

8.6 Start the Run

8.6.1 Run sample(s) - Before proceeding


Before starting the run:

- Ensure to have correctly prepared the sample vials and the sample list in the analyzer software.
- Ensure methods and, if desired, sequences, have already been set.
- Verify that the parameter values of the desired method are correct

To modify a method or set a sequence, see ['Set the method and sequence' on page 205](#).

8.6.2 Run sample(s) - Start a sequence

This procedure is useful when there are many sample vials to be processed.


 *Run menu > Load Sample Tray*

1. Insert the sample vials in the tray.
2. Select **Close**.
3. Start the sample list on the analyzer software.
4. In the **Run** menu select **Start**: the autosampler starts the sequence and the summary page appears.
5. At the end of the run, select **Continue**: the **Run** menu appears.

8.6.3 Run sample(s) - Process a single sample

This procedure is useful when there is only one sample vial to be processed or to perform some test.


Note: *This procedure is applicable only if the parameter **User Interface** has been set to **Standard**.*

 *Run menu > Load Sample Tray*

1. Insert the sample vial in the tray.
2. Start the sample list on the analyzer software.
3. In the Standard version of the Run menu select **Single Injection**.
4. Set the parameters as desired.
5. Select **Start**: the autosampler starts the single injection and the summary page appears.

8.7 Run management

8.7.1 Perform additional fiber cleaning conditioning

 *Run menu > Fiber cleaning conditioning*


1. Set **Position** and **Remaining time**.
2. Set the desired parameters and select **Start**: the autosampler performs the defined additional fiber cleaning conditioning.

8.7.2 Run sample(s) - Stop a sequence

1. Select **Stop**: the autosampler pauses.
2. Do what follows:

If you want...	Then...
to put the autosampler into standby	select Abort : the autosampler washes the syringe and goes into standby.
to stop the autosampler	select Terminate Current : the autosampler stops when the operations on the current sample vial end.
to continue the sequence	select Resume .

8.7.3 Modify the autosampler general settings

 *Run menu > Settings > Setup > Configuration*

Select the desired parameters and modify them. For details, see '[Configuration](#)' on page [318](#).

8.8 Sequences management

8.8.1 Add a sequence step copying from another one

 *Run menu > Settings > Sequences > Tools*


1. Select **Copy as new** and select the step to copy from.
2. Select **Copy**.
3. Select **Continue**.

8.8.2 Reorder sequence step

 *Run menu > Settings > Sequences > Tools*

1. Select **Reorder**.
2. Select the step to move and select **OK**.
3. Select the target position where you want to move the step and select **Move**.
4. Select **Continue**.

8.8.3 Remove sequence step

 *Run menu > Settings > Sequences > Tools*

1. Select **Remove Step**.
2. Select the step to remove and select **OK**.
3. Select **Continue**.

8.9 Racks management

8.9.1 Set the characteristics of the rack to be used

This procedure adds, to the list of the available racks to be used, a rack with a different vial type than the previous one.

✱ *Run menu > Settings > Setup > Alignment > Tray Install*

1. Select **Add New** and select **OK**.
2. Select the vial type, the vial depth and press **Continue**.
3. The new list item appears.
4. Select **Exit > Exit**.

8.9.2 Choose the type of rack to be used

✱ *Run menu > Settings > Tray Type*

1. Select the desired type of rack and select **OK > Mount** to confirm: the tray opens.
2. Follow the screen instructions to place the required sample vial, then select **Continue**: the vial locator touches the sample vial.
3. Select **Continue**.

8.9.3 Remove a rack from the list

✱ *Run menu > Settings > Setup > Alignment > Tray Install*

1. Select **Remove** and select **OK**.
2. Select the tray to be removed and select **Remove**.
3. Select **Continue > Exit > Exit**.

8.10 Fiber management


8.10.1 Workflow for replacing the fiber

For replacing the fiber, follow these procedures in the stated order:


- 1. ['Remove the fiber holder' on page 210](#).
- 2. ['Remove the fiber' on page 212](#).
- 3. ['Install the fiber' on page 212](#).
- 4. ['Install the fiber holder' on page 213](#).


Each procedure is described in detail below.

8.10.2 Remove the fiber holder



 Run menu > Settings > Setup > Fiber Installation

- 1. Raise the sliding lid.
- 2. Select the **Fiber type** you want to install and press **Save**.
- 3. Do what follows:

If the autosampler asks to...	Then...
remove the fiber holder or the sealed tip holder	proceed to point 4.
remove the vial locator	<div>1. Remove the two fixing screws to remove the vial locator.</div> <div></div>

If the autosampler asks to...	Then...
	 <p data-bbox="735 837 970 871">2. Select Continue.</p>


4. Unscrew the two finger nuts.
5. Take the syringe pointer, insert its edge in the plunger locker and pull it out.
6. Carefully remove the fiber holder.




7. Press **Continue**.

8.10.3 Remove the fiber


1. Unscrew the fiber guide and remove it.



2. Lower the fiber holder plunger until the fiber fixing point is visible.




3. Unscrew the fiber fixing point to remove the fiber.




8.10.4 Install the fiber


1. Screw the fiber fixing point to the fiber holder plunger.



2. Lift up the fiber holder plunger until the fiber is completely inside the fiber holder.



3. Insert the fiber guide along the fiber and fix it at the fiber holder.



Note: Take care to use the standard fiber guide in case of standard fiber or the sealed tip guide in case of sealed tip.




8.10.5 Install the fiber holder


1. When the screen displays the **Install new fiber holder or Sealed tip holder** page, do what follows:

If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the fiber holder installation without moving the turret	proceed to point 2.


- 2. If necessary, raise the sliding lid and unscrew the two finger nuts.
- 3. If necessary, take the syringe pointer, insert its edge in the plunger locker and pull it out.
- 4. Wait for the automatic positioning.

If the vial locator...	Then...
is already mounted	1. Push the vial locator up and insert the fiber in the central hole of the vial locator.


If the vial locator...	Then...
	<div data-bbox="475 293 842 792">  </div> <div data-bbox="922 293 1289 792">  </div> <p data-bbox="395 840 710 873">2. Install the fiber holder.</p> <div data-bbox="443 920 778 1442">  </div>
is not mounted yet	install the fiber holder.

If the vial locator...	Then...
	

5. Fix completely the finger nuts.






6. Insert the plunger locker.



7. Press **Continue**.

8.


If...	Then...
<p>the autosampler asks to mount the vial locator</p>	<ol style="list-style-type: none"> <li data-bbox="397 309 1355 443">1. Install the vial locator for Headspace and SPME standard fibers technique (blue o-ring: see 'Vial locators' on page 55), taking care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars.  <ol style="list-style-type: none"> <li data-bbox="397 987 932 1021">2. Fasten completely the two fixing screws. <li data-bbox="397 1028 1307 1095">3. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the fiber.  <ol style="list-style-type: none"> <li data-bbox="397 1653 1355 1753">4. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.


If...	Then...
	
the autosampler executes the Automatic Fiber / Sealed tip Zero procedure	5. Select Continue . proceed to point 9.

9. Wait the end of the of the Automatic Fiber / Sealed tip Zero procedure.


8.10.6 Align the vial locator to the needle tip

- 1. Remove the fiber holder. For details, see ['Remove the fiber holder' on page 210](#).
- 2. Unscrew the safety lock and the needle higher regulator and let the vial locator to rest on the frontal panel.







3. Insert the syringe pointer as indicated in the following pictures:



4. Lift up the vial locator until the syringe pointer is completely against the sledge.



5. Using a flat screwdriver, screw the safety lock and the needle height regulator to block the vial locator in position.
6. Remove the syringe pointer.



7. Install the fiber holder.

8.11 Sealed tip management


8.11.1 Workflow for replacing the sealed tip

For replacing the sealed tip, follow these procedures in the stated order:


1. ['Remove the sealed tip holder' on page 220](#).
2. ['Remove the sealed tip' on page 222](#).
3. ['Install the sealed tip' on page 223](#).
4. ['Install the sealed tip holder' on page 223](#).

Each procedure is described in detail below.


8.11.2 Remove the sealed tip holder

 *Run menu > Settings > Setup > Fiber Installation*

1. Raise the sliding lid.
2. Select the **Fiber type** you want to install and press **Save**.
3. Do what follows:



If the autosampler asks to...	Then...
remove the fiber holder or the sealed tip holder	proceed to point 4.
remove the vial locator	<div>1. Remove the two fixing screws to remove the vial locator.</div> <div></div>

If the autosampler asks to...	Then...
-------------------------------	---------




2. Select **Continue**.

4. Take the syringe pointer, insert its edge in the plunger locker and pull it out.
5. Remove the plunger extension.



6. Unscrew the two finger nuts.

7. Carefully remove the sealed tip holder.




Note: In case you are going to mount the fiber holder for standard fiber, remove the sealed tip spacer too.

8. Press **Continue**.


8.11.3 Remove the sealed tip

1. Remove the sealed tip fiber from the sealed tip holder.



8.11.4 Install the sealed tip

1. Insert the sealed tip fiber into the sealed tip holder.




8.11.5 Install the sealed tip holder



1. When the screen displays the **Install new fiber holder or Sealed tip holder** page, do what follows:


If you want to...	Then...
move the turret to a more comfortable position	1. Select Manual . 2. Select the arrows for moving the syringe location as you desire.
proceed with the sealed tip holder installation without moving the turret	proceed to point 2.

2. If necessary, raise the sliding lid and unscrew the two finger nuts.
3. If necessary, take the syringe pointer, insert its edge in the plunger locker and pull it out.
4. If necessary, insert the sealed tip spacer in the position shown in the following picture.




5. Wait for the automatic positioning.


If the vial locator...	Then...
is already mounted	<div data-bbox="400 253 1353 322"><p>1. Push the vial locator up and insert the sealed tip in the central hole of the vial locator.</p></div> <div data-bbox="475 360 786 819"></div> <div data-bbox="400 869 772 902"><p>2. Install the sealed tip holder.</p></div> <div data-bbox="440 943 956 1458"></div>
is not mounted yet	install the sealed tip holder.

If the vial locator...	Then...
	


5. Fix completely the finger nuts.




6. Insert plunger extension.



7. Lower the plunger extension until it rests on the plunger holder.






8. Insert the plunger locker.



9. Press **Continue**.

10.


If...	Then...
the autosampler asks to mount the vial locator	<div>1. Install the vial locator for SPME for 1.1 fibers / 1.5 fibers technique (yellow o-ring: see 'Vial locators' on page 55), taking care to let the needle tip to enter the central hole on the vial locator and to insert the two bars along the two long bars.</div> <div></div> <div>2. Fasten completely the two fixing screws.</div>

If...	Then...
	<p>3. By pushing from the indicated point (green circles), manually lift the vial locator until it reaches the sealed tip.</p>  <p>4. Keeping your finger under the vial locator, guide the vial locator down until it returns to its original position, checking that it is free to lower, without the need to pull it.</p>  <p>5. Select Continue.</p>
the autosampler executes the Automatic Fiber / Sealed tip Zero procedure	proceed to point 9.



9. Wait the end of the of the Automatic Fiber / Sealed tip Zero procedure.

8.11.6 Align the vial locator to the needle tip


1. Remove the sealed tip holder. For details, see ['Remove the sealed tip holder' on page 220](#).
2. Unscrew the safety lock and the needle higher regulator and let the vial locator to rest on the frontal panel.



3. Insert the syringe pointer as indicated in the following pictures:




4. Lift up the vial locator until the syringe pointer is completely against the sledge.



5. Using a flat screwdriver, screw the safety lock and the needle height regulator to block the vial locator in position.

6. Remove the syringe pointer.



7. Install the sealed tip holder.

8.12 Obtain support in case of autosampler problems

This procedure allows to send an auto-generated email with useful technical information about autosampler.


Note: *This procedure is applicable only if the parameter **Log** has been set to **Enable** or **Enable w/o warnings**.*

 *Run menu*

1. Select **Generate QR code**. The Generate QR code wizard appears on the PC screen.
2. Follow the wizard instructions.
3. Select **Continue**.

8.13 Change the analyzer to be used


✱ *Run menu > Settings > Setup > Alignment > Analyzer Default*

1. If the name of the analyzer to be used is already listed, select it. Otherwise, select **Generic**.
2. Select **OK** to confirm.
3. Select the number of injectors to be used and then **OK** to confirm and **Load**.
4. Select **Continue**.
5. Proceed with the ['Align injectors' procedure on the following section.](#) 

8.14 Realign injectors

☛ *Run menu > Settings > Setup > Alignment*


- 1. To align the autosampler to the front injector, select **Front Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the front injector adapter.

 **CAUTION** Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

- 2. Do what follows:

If you want to...	Then...
confirm the injector alignment	select Confirm .
align the autosampler to the injector manually	<ul style="list-style-type: none">1. Select Do Manual.2. Select the arrows to move the vial locator to the correct position (Touch Sensor must be ON).3. Select Store to confirm.


- 3. Set the desired injection depth: the autosampler tries to insert the needle in the injector.
- 4. Select **Save** to confirm.
- 5. Select the desired injection speed (from **Very Low** to **Very High**) and select **Save**: the turret returns to its original position and the **Alignment** menu appears.
- 6. If the rear injector is present, select **Rear Injector**: the turret moves and the vial locator touches a known part of the autosampler and then the rear injector adapter.



 **CAUTION** Vial locator with needle in motion. Do not touch nor adjust the injector adapter.

- 7. Repeat from step 2 to 7.

Vial locator aligned to the injector adapter

To be aligned, the vial locator must lean against the top of the injector adapter completely.





8.15 Reset the position of the vial


✱ *Run menu > Settings > Setup > Alignment*

1. Select **Touch & Plunger Zero**: the turret moves and the tray opens.
2. Insert a crimped sample vial in position A4 of the rack and in C position of the oven (if using 6ml or 10 ml vial, ensure to insert first the proper spacer in C position). In case wash vial is used, insert it in its position.
3. Select **Continue**: the vial locator touches the wash vial, the sample vial in the rack, the sample vial in the oven, the cover and the injector adapter to memorize their positions.
4. Remove the vial from oven and press **Continue**.
5. At the end of the procedure, when the **Setup menu** appears, select **Exit > Exit**.


8.16 Get out of SPME mode (fiber)

☀ *Run menu > Change Mode*

- 1. Select the mode you are going to use and select **Continue**.
- 2. Lift the sliding lid.
- 3.

If the autosampler...	Then...
asks to Remove the vial locator	<div>1. Remove the two fixing screws to remove the vial locator.</div> <div></div> <div>2. Select Continue.</div>
asks to Remove the fiber holder	proceed to point 4.

- 4. Unscrew the finger nuts.
- 5. Take the syringe pointer, insert its edge in the plunger locker and pull it out.
- 6. Carefully remove the fiber holder.




7. Press **Continue**.

8. Follow the instructions of the 'Get in' chapter of the mode you are going to use.


8.17 Get out of SPME mode (sealed tip)

- 1. Select the mode you are going to use and select **Continue**.
- 2. Lift the sliding lid.
- 3. Do what follows:

If the autosampler asks to...	Then...
remove the fiber holder or the sealed tip holder	proceed to point 4.
remove the vial locator	<div><div>1. Remove the two fixing screws to remove the vial locator.</div><div></div><div>2. Select Continue.</div></div>


- 4. Take the syringe pointer, insert its edge in the plunger locker and pull it out.

5. Remove the plunger extension.



6. Unscrew the two finger nuts.

7. Carefully remove the sealed tip holder.



9. Screen description in Headspace mode

9.1 Screen elements

Note: The interface may be updated in future. In this case, please follow the screen instructions.




9.1.1 Interacting with the touch screen

The touch screen can be tapped with a finger or a capacitive touch screen stylus.

9.1.2 Types of screens


The screens of the touch screen are mainly of two types: menus, with icons, or pages, with a list.

9.1.3 Common buttons

Part	Function
	Move up and down in a list and select an item.
	<ul style="list-style-type: none"> • Move the needle/plunger up and down. • Open/Close the tray. • Increase and decrease the visualized parameter.
	Move the turret left and right.
OK	<ul style="list-style-type: none"> • Confirm the selection. • Confirm a parameter value. • Commute a signal status.
Save	Save the modifications.
Next	Go to the next page.
Exit	Exit the menu.
Cancel	Exit the page without saving the modifications.
Continue	Confirm and go to the next page.
Abort	Stop the autosampler immediately.
Stop	Pause the running operation.
Menu	Access additional options during the automatic run.
Resume	Continue a paused sequence.
Terminate current	Stop the sequence. The sequence stops only after having completed the injection in progress in the analyzer.


9.2 Main menus

9.2.1 Screen flowchart



9.2.2 Run menu

Quick Start version




Standard version



List item	Function
Start	Only in the quick start version. Process all the sample vials listed in all the sequence steps stored. After two consecutive empty positions, the autosampler stops.
Single Injection	Only in the standard version. Process a single sample vial. For details, see 'Single Injection' on page 247 .
Automatic Run	Only in the standard version. Perform a sequence of injections from different sample vials. For details, see 'Automatic Run' on page 247 .
Syringe Purge	Allow an extra purge of the syringe using the flush gas. For details, see 'Syringe Purge' on page 248 .
Load Sample/Solvent	Open the tray to load/unload the sample vials and racks.
Settings	Access the Settings menu. For details, see 'Settings' on page 248 .
Injector Septa Test	Check if the syringe can pierce the injector septum easily and correctly.
Generate QR code	Open the Generate QR code wizard on the PC in the HTA Monitor.
Change Mode	Change the operation mode (Liquid or SPME).

9.2.3 Automatic Run menu

 (during the automatic run) Menu





List item	Function
Methods	Show the methods menu.
Sequences	Show the sequences menu.
Load Sample Tray	Open the tray to add sample vials.
Immediate Injection	Skip the remaining conditioning time and inject immediately.
Skip Next Vial	Skip the remaining conditioning time and unload the vial, without to perform the injection.
Sample Info	View the remaining conditioning time of the samples in the oven.
Run Info	View important details of the run.
Exit	Exit the menu and return to the summary page.

9.3 The Run menu (main menu)

9.3.1 Single Injection


✱ *Run menu > Single Injection*



List item	Description
Start Injection	Start the single injection set.
Method	From 0 to 9.
Sample	Position of the sample vial to be injected in the analyzer.
Exit	Exit the page and return to the Run menu.

9.3.2 Automatic Run

✱ *Run menu > Automatic Run*




List item	Function
Start Autorun	Start the set sequence.

List item	Function
First Step	Set the first step of the sequence.
Last Step	Set the last step of the sequence.
Run-All-Samples	Set the behaviour of the autosampler when it detects two consecutive empty positions: <ul style="list-style-type: none"> • Yes: the autosampler stops. • No: the autosampler searches for other sample vials.
Exit	Exit the page and return to the Run menu.

9.3.3 Syringe Purge


✱ *Run menu > Syringe Purge*



List item	Description
Start Purge	Run the syringe purge.
Flush Time [min]	Duration of the purge procedure.
Purge Temperature [°C]	Temperature of the purge procedure.
Exit	Exit the page and return to the Run menu.

9.3.4 Settings

✱ *Run menu > Settings*





List item	Function
Methods	Set the method parameters.

List item	Function
	For details, see 'Methods' on page 250 .
Sequences	Set the sequence parameters. For details, see 'Sequences' on page 251 .
Tray Type	Set the rack type to be used. For details, see 'Choose the type of rack to be used' on page 130 .
View Parameters	Set the autosampler parameters. For details, see 'View parameters' on page 251 .
Setup	Access the menu for autosampler setup. For details, see 'Setup' on page 253 .
Eco Savings	Set the oven stand by temperature. <ul style="list-style-type: none"> • St-by Oven: allow to set the oven temperature in stand-by conditions. • St-by Syringe: allow to set the syringe temperature in stand-by condition. • St-by Flush: allow to set if syringe purge is continuously active in stand-by condition.
Exit	Exit the menu and return to the Run menu.

9.4 The Settings menu

9.4.1 Methods


✱ Run menu > Settings > Methods



List item	Function
Method	Select the method to be modified. For details, see 'The method menu' on page 269 .
Tools	Manage the method. For details, see 'Tools' on page 250 .
Exit	Exit the page and return to the Settings menu.

9.4.2 Tools

✱ Run menu > Settings > Methods > Tools




List item	Function
Restore Default	Load the default settings for a method.

List item	Function
Copy Methods	Copy the parameters from one method to another.
Exit	Exit the page and return to the Methods menu.

9.4.3 Sequences


 *Run menu > Settings > Sequences*



List item	Function
Step	Select the desired step. For details, see 'Step' on page 259 .
Tools	Manage the sequence steps. For details, see 'Tools' on page 260 .
Exit	Exit the page and return to the Settings menu.

9.4.4 View Parameters

 *Run menu > Settings > View Parameters*



Run: View Parameters 2/3

Front Injector Insertion Speed: **Low**

Rear Injector Depth [mm]: **0**

Rear Injector Insertion Speed: **Low**

IP Address: **192.168.0.235**

Subnet Mask: **255.255.0.0**

Exit

Next

Run: View Parameters 3/3


Gateway: **192.168.0.1**


Exit

Next

List item	Description
Syringe Volume [ml]	Volume of the installed syringe.
Tray Type	Tray type installed, with number of vial positions, dimension of sample vials and needle draw depth inside the vial.
Analyzer Brand	Brand of the analyzer to which the autosampler is connected.
Analyzer Model	Model of the analyzer to which the autosampler is connected.
Front Injector Depth [mm]	Needle insertion depth into the front injector.
Front Injector Insertion Speed	Needle insertion speed into the front injector.
Rear Injector Depth [mm]	Only with rear injector. Needle insertion depth into the rear injector.
Rear Injector Insertion Speed	Only with rear injector. Needle insertion speed into the rear injector.
IP Address	LAN connection data.
Subnet Mask	
Gateway	
Exit	Exit the page and return to the Settings menu.

9.4.5 Setup

 *Run menu > Settings > Setup*




List item	Function
Alignment	Access the autosampler alignment procedures to adjust the working positions. For details, see 'Alignment' on page 254 .
Syringe Installation	Install the syringe.
Manual Operation	Move the plunger, needle, turret and tray for diagnostic or emergency needs. Status: <ul style="list-style-type: none">• RUN: moving• IDLE: not moving• BUSY: an error impedes the movement• NEED ENC ALIGN: the motor has not been aligned yet (only for service purpose) For details, see 'Manual operation' on page 254 .
Configuration	Edit the autosampler general settings. For details, see 'Configuration' on page 257 .
Service	Password protected. Allow service operations on the autosampler.
Exit	Exit the menu and return to the Run menu.

9.5 The Setup menu

9.5.1 Alignment

✱ *Run menu > Settings > Setup > Alignment*








List item	Function
Analyzer Default	Select the analyzer on which the autosampler is mounted.
Front Injector	Enabled only if the analyzer is selected. Align the vial locator to the front injector.
Rear Injector	Enabled only if the analyzer is selected. Align the vial locator to the rear injector.
Touch Plunger Zero	Perform the autosampler alignment procedure.
Tray Install	List the racks needed for the sequence with the related draw depth (mm) and manage the list.
Exit	Exit the menu and return to the Setup menu.

9.5.2 Manual Operation

✱ *Run menu > Settings > Setup > Manual Operation*

This page is useful for troubleshooting.








List item	Description
Plunger	<ul style="list-style-type: none"> • Show the status sent by the plunger motor. • If selected, it allows to move the plunger motor.
Needle	<ul style="list-style-type: none"> • Show the status sent by the needle motor. • If selected, it allows to move the needle motor.
Turret	<ul style="list-style-type: none"> • Show the status sent by the turret motor.

List item	Description
	<ul style="list-style-type: none"> If selected, it allows to move the turret motor.
Tray	<ul style="list-style-type: none"> Show the status sent by the tray motor. If selected, it allows to move the tray motor.
Touch sensor	Status of the touch sensor: <ul style="list-style-type: none"> On: the vial locator has detected the object. Off: the vial locator has not detected the object.
Cover	<ul style="list-style-type: none"> Show the status sent by cover motor. If selected, it allows to move the cover motor.
Shaker	Allow to select a defined speed and, by pressing OK button, to activate it at the defined speed. It is used for service purpose.
Shaker sensor	Status of the shaker sensor: <ul style="list-style-type: none"> On: the shaker sensor is on. Off: the shaker sensor is off. It is used for service purpose.
Fan	Allow to select a defined speed and, by pressing OK button, to activate it at the defined speed. It is used for service purpose.
Ready	Signal sent by the analyzer to communicate that it has finished the analysis. <ul style="list-style-type: none"> On Off
Start	Signal sent by the autosampler to communicate that it has started the run. <ul style="list-style-type: none"> On Off Select OK to switch the signal.
Flush	Status of the Flush valve. Select OK to switch the Flush valve status.
Sync-out	Status of the sync- out signal. Select OK to switch the signal.
Sync-in	Status of the sync- in signal. <ul style="list-style-type: none"> On Off
Prep	Signal sent by the autosampler to communicate to quit the PREP status and to go in READY status. <ul style="list-style-type: none"> On Off Select OK to switch the signal.
Exit	Exit the page and return to the Setup menu.

9.5.3 Configuration

 *Run menu > Settings > Setup > Configuration*

This menu contains general settings of the autosampler.



List item	Function
Sound	Change the sound effects of the autosampler: <ul style="list-style-type: none"> • Level 2: emit a sound when the screen is selected • Level 1: emit a sound in the following cases: <ul style="list-style-type: none"> ○ The single injection is finished. ○ The sample has been injected in the analyzer. ○ Save or Store have been selected. ○ The autosampler memory has been erased. • Level 0: no sound
Light	Change the syringe illumination settings: <ul style="list-style-type: none"> • On: syringe LED on • Off: syringe LED off • Injection: syringe LED on during injection only
Run	Change the settings for the injection sequence. For details, see 'Run' on page 261 .
Input Output	Change the input/output parameters to/from the analyzer. For details, see 'Input Output' on page 262 .
Maintenance	Show information and allow to change the settings of the preventive maintenance counters. For details, see 'Maintenance' on page 263 .
Advanced	Password protected. Edit the advanced parameters of the autosampler.

List item	Function
Activation	Insert the password for the optional features.
BCR	Show information about the sample vial barcode and change the barcode reader settings. For details, see 'BCR' on page 265 .
Exit	Exit the menu and return to the Setup menu.

9.6 The Sequences page

9.6.1 Step

✱ *Run menu > Settings > Sequences > Step*

Run: Sequences

Step: 1	↑
Tools	
Exit	OK
	↓

Step: 1

Tray Type: 42 positions	↑
Method: 0	
First Sample: A.1	OK
Last Sample: A.1	
Injection Mode: Default	↓
Name:	


Step: 1


Method: 0	↑
First Sample: A.1	
Last Sample: A.1	OK
Injection Mode: Default	
Name:	↓
Exit	

List item	Description
Tray Type	Type of rack to be used.
Method	Method to be used.
First Sample	Position of the first sample vial to be processed with the selected method.
Last Sample	Position of the last sample vial to be processed with the selected method. All the sample vials from the position set in First Sample to this one will be processed.
Injection Mode	<ul style="list-style-type: none">Default: the autosampler injects in the injector defined in 'Run' on page 261.Front: the autosampler injects only in the front injector.Rear: the autosampler injects only in the rear injector.Confirmation: the autosampler injects the same sample in both injectors.

List item	Description
Name	Name of the sequence step. Editable from HTA Autosampler Manager.
Exit	Exit the page and return to the Sequences page.

9.6.2 Tools

 *Run menu > Settings > Sequences > Tools*





List item	Function
Add New	Add a new sequence step.
Copy as New	Create a new sequence step by copying the parameters of an existing ones.
Reorder	Change the order of the sequence steps.
Remove	Remove a sequence step.
Exit	Exit the page and return to the Sequences page.

9.7 The Configuration menu

9.7.1 Run

✱ *Run menu > Settings > Setup > Configuration > Run*

This page allows to change the autosampler behaviour during the run.




List item	Function
Tray Behaviour	Set if the tray must stay open or closed during the run.
Power-on Restart	Set the autosampler restart after a power failure during the run: <ul style="list-style-type: none"> • On: the run restarts. • Off: the run does not restart and the Run menu appears.
Injection Mode	Only if two injectors are present. Set the default injector in which to inject the sample. <ul style="list-style-type: none"> • Front: the autosampler injects only in the front injector. • Rear: the autosampler injects only in the rear injector. • Confirmation: the autosampler injects the same sample in both injectors.
Injection Synchro	Set the autosampler behaviour in relation to the analyzer. For details, see 'Injection Synchro' on page 266 .
Shaker Speed	Define the shaker speed used during the run. <ul style="list-style-type: none"> • Very Low • Low • Normal • High • Very High
Maintenance Warning	Set when the maintenance warning messages can appear: <ul style="list-style-type: none"> • Normal: at Start up and at the beginning of each run, if the set limit is reached.

List item	Function
	<ul style="list-style-type: none"> • Only at Start up: at the beginning of each run, if the set limit is reached.
User Interface	Set the version of the Run menu: <ul style="list-style-type: none"> • Standard: version with all buttons • Quick Start: version with less buttons
Prep-Ahead Waiting Pos.	For special application only.
Log	Set if the autosampler must send the logs to HTA Monitor. <ul style="list-style-type: none"> • Enable: the autosampler sends logs to HTA Monitor, • Enable w/o warnings: the autosampler sends logs to HTA Monitor, but does not show warnings in case the PC is not connected to the autosampler. • Disable: the autosampler does not send logs to HTA Monitor.
Log Level	Set the type of logs the autosampler sends to the HTA Monitor. <ul style="list-style-type: none"> • Normal: the autosampler sends the standard logs, • Troubleshooting: the autosampler sends more detailed logs to HTA Monitor.
Protocol version	Set the communication protocol version used by the autosampler. Note: Do not change unless: <ul style="list-style-type: none"> ○ you operate by CDS through HTA software connectors ○ required by your third party software documentation ○ required by your Customer Representative
Exit	Exit the page and return to the Configuration menu.

9.7.2 Input Output

✱ Run menu > Settings > Setup > Configuration > Input Output



List item	Function
Ready On	Set the ON state of the Ready signal. This signal is used to start the autosampler and must match the signal given by the analyzer. <ul style="list-style-type: none"> • Low: the autosampler starts the injection when the analyzer provides a low logical level of voltage for the Ready signal. • High: the autosampler starts the injection when the analyzer provides a high logical level of voltage (+5 V) for the Ready signal.
Start On	Set the ON state of the Start signal from the autosampler to the analyzer during the injection. This parameter has to be set according to the analyzer. <ul style="list-style-type: none"> • Open: during the injection, the pins 3 (SAMPINS- NO) and 6 (SAMPINS COM) are


List item	Function
	<p>open for about one second. They are closed before and after the injection.</p> <ul style="list-style-type: none"> • Close: during the injection, the pins 3 (SAMPINS- NO) and 6 (SAMPINS COM) are closed for about one second (short- circuit). They are open before and after the injection.
Prep Mode On	<p>Set the ON state of the Prep Mode signal from the autosampler to the analyzer during the run. This parameter has to be set according to the analyzer.</p> <ul style="list-style-type: none"> • Open: before the processing of the sample vial, the pins 11 (ENDSAMP - NO) and 14 (ENDSAMP- COM) are open for about one second. If within one minute the Ready signal is not received, this signal is activated again. This operation is repeated until the Ready signal is received. • Close: before the processing of the sample vial, the pins 11 (ENDSAMP - NO) and 14 (ENDSAMP- COM) are closed for about one second (short- circuit). If within one minute the Ready signal is not received, this signal is activated again. This operation is repeated until the Ready signal is received.
Exit	Exit the page and return to the Configuration menu.

9.7.3 Maintenance

✱ *Run menu > Settings > Setup > Configuration > Maintenance*

Setup: Maintenance	
Plunger Strokes Counter: 36	↑
Plunger Strokes Limit: 25000	
Injection Counter: 19	OK
Injection Limit: Off	
System Integrity test: Disabled	↓
Penetration Counter: 14	


Setup: Maintenance	
Penetration Counter: 14	↑
Penetration Limit: 10000	
Syringe expire within [months]: ---	OK
Syringe Expiry [months]: 12	
Current Date: Start PC to use AI features.	↓
Injector Septa Testing: Enabled	



List item	Description
Plunger Strokes Counter	Number of plunger strokes since the last reset. It allows to monitor the syringe wear- and- tear. This counter is automatically reset after a new syringe installation.
Plunger Strokes Limit	Limit of plunger strokes. Once reached, the autosampler shows a maintenance warning message. Off: no limit
Injection Counter	Number of injections performed in the analyzer since the last reset. It allows to predict when the analyzer septum or injection liner should be replaced.
Injection Limit	Limit of injection number. Once reached, the autosampler shows a maintenance warning message. Off: no limit
System Integrity Test	Procedure to check syringe status before to execute a sequence. <ul style="list-style-type: none"> • Enable • Disable To be performed, it requires the 'Syringe integrity tool' option.
Penetration Counter	Total number of perforation of the sample vial septa, and of the injector. It allows to predict when the syringe should be replaced, due to needle wear.
Penetration Limit	Limit of perforation of the sample vial septa, and of the injector. Once reached, the autosampler shows a maintenance warning message. Off: no limit
Syringe expiry within [months]	Only if connected to the HTA Monitor. Remaining months to the syringe expiration. It allows to predict when the syringe should be replaced: it is particularly relevant when the autosampler performs a few injection per day.
Syringe Expiry [months]	Only if connected to the HTA Monitor. Month in which the syringe expires. Once reached, the autosampler shows a maintenance warning message. Off: no limit
Current Date	Only if connected to the HTA Monitor. Current Date.
Injector Septa Testing	Procedure to avoid damaging the syringe needle during the injection. <ul style="list-style-type: none"> • Enable • Disable
Turret Testing	Initializing procedure performed by the autosampler every times it is switched on. It checks the correct movement of the turret to verify the good functioning of the autosampler. <ul style="list-style-type: none"> • Enable • Disable

List item	Description
Execute Turret Baseline	Perform the baseline necessary to perform the Turret Testing. To be performed only on Customer Representative request.
Execute Turret Testing	Perform a Turret Testing.
Exit	Exit the page and return to the Configuration menu.

9.7.4 BCR

 *Run menu > Settings > Setup > Configuration > BCR*

List item	Description
Enable	Read the sample vial barcode and save it in the autosampler memory so that the information contained can be read by the available software: <ul style="list-style-type: none"> • Yes • No
On Reading Error	Behaviour if a reading error occurs: <ul style="list-style-type: none"> • Skip: skip the run for the vial. • Process: perform the run anyway.
Test	Test the alignment of the turret in front of the barcode reader for correct reading.
Exit	Exit the page and return to the Configuration menu.

9.8 The Run page

9.8.1 Injection Sychro

✱ Run menu > Settings > Setup > Configuration > Run > Injection Sychro

Note: Prep Mode is available only for some analyzer. Check in your autosampler mounting kit instruction manual if this feature is supported on your analyzer.




Note: If using **Prep Mode**, it is mandatory to start sample list on analyzer software before to start the run on the autosampler, otherwise the analysis could be not executed and/or not saved properly.

Injection Sychro	Sample preparation start	Sample injection	Start signal to the analyzer
Normal	At the reception of the Ready signal from the analyzer.	After ending sample preparation the autosampler checks again the presence of the Ready signal from the analyzer and then starts the injection.	At the beginning of the syringe plunger movement.
Trigger			
A-Start			When the autosampler touches the injector.
D-Start			At the end of the syringe plunger movement.
EA	At the end of the analysis time set in the autosampler method.	After the ending sample preparation.	At the beginning of the syringe plunger movement.
Normal w/o Rdy			
DIN	At the reception of the Ready signal from the analyzer.	No injection is carried out.	At the beginning of the syringe plunger movement, when needle has reached the proper depth inside the vial.
Prep Mode	At the end of the analysis time set in the autosampler method, the autosampler sends a PREP signal to the analyzer and wait for the Ready signal from the analyzer.	After the autosampler has checked that the Ready signal is still present.	When the syringe plunger starts to move with the syringe in the injector.

9.9 The Method page

9.9.1 Method

✱ *Run menu > Settings > Methods > (method) OK*



List item	Function
General	Access the page for setting the general parameters of the method. For details, see 'General' on page 269 .
Advanced	Access the page for setting the advanced parameters of the method. For details, see 'Advanced' on page 269 .
Preparation	Access the page for setting the preparation parameters of the method. For details, see 'Preparation' on page 270 .
Sample	Access the page for setting the sample parameters of the method. For details, see 'Sample' on page 270 .
Injection	Access the page for setting the injection parameters of the method. For details, see 'Injection' on page 271 .

List item	Function
Exit	Exit the menu and return to the Methods page.

9.10 The Method menu

9.10.1 General

✱ Run menu > Settings > Methods > (method) OK > General

Method 0 General	
Syringe Volume [ml]: 2.5	↑
Analysis Time [min]: 1	
Method Type: Constant Mode	OK
Preparation Time: 00:00:43	
Exit	↓

List item	Description
Syringe Volume [ml]	Volume of the syringe. The syringe volume affects the limits of other parameters such as speed and volume.
Analysis Time [min]	Time of the sample analysis.
Method Type	Parameter not editable.
Preparation Time	Display how long the autosampler needs to perform a cycle. It is calculated from the instrument and is not editable (Time needed to transfer the vial from the tray to the incubator, filling time, injection time, etc.).
Exit	Exit the page and return to the selected method page.

9.10.2 Advanced

✱ Run menu > Settings > Methods > (method) OK > Advanced

Method 0 Advanced	
Enrichment Cycle(s): 1	↑
Dwell Between Injections [min]: ---	
Pre-Flush: Off	OK
Exit	
	↓

List item	Function
Enrichment Cycle(s)	The number of enrichments in the injection cycle. The enrichment allows the user to inject a higher quantity of sample. Note: It is only possible to enrich the sample by multiple injections if the gas chromatograph can support this feature by a proper trapping system.
Dwell Between Injections [min]	Delay time between any enrichment, to allow the headspace vial to equilibrate.
Pre-Flush	Let the flush gas to pass thorough the syringe immediately before to insert the needle in the vial; it is used to eject the room air from the syringe before to take the

List item	Function
	sample.
Exit	Exit the page and return to the selected method page.

9.10.3 Preparation

✱ Run menu > Settings > Methods > (method) OK > Preparation

Method 0 Preparation	
Syringe Temperature [C]: Off	↑
Oven Temperature [C]: Off	
Incubation [min]: 0	OK
Shaker On [min]: 0.1	
Shaker Off [min]: 0.1	↓
Exit	

List item	Function
Syringe Temperature [°C]	Temperature of syringe conditioning. This value can be between 40°C and 150°C, or off (inactive conditioning).
Oven Temperature [°C]	Temperature of oven conditioning. This value can be between 40°C and 170°C, or Off (inactive conditioning).
Incubation Time [min]	Sample incubation time inside the oven.
Shaker On [min]	Interval time during incubation time in which the sample is shaken.
Shaker Off [min]	Interval time during incubation time in which the sample is not shaken.
Exit	Exit the page and return to the selected method page.

9.10.4 Sample

✱ Run menu > Settings > Methods > (method) OK > Sample

Method 0 Sample	
Sample Volume [ml]: 2.25	↑
Fill Volume [ml]: 1.00	
Pull Up Strokes: 0	OK
Equilibration Delay [s]: 0	
Sample Speed [ml/min]: 34.0	↓
Syringe Prefill: No	


Method 0 Sample	
Fill Volume [ml]: 1.00	↑
Pull Up Strokes: 0	
Equilibration Delay [s]: 0	OK
Sample Speed [ml/min]: 34.0	
Syringe Prefill: No	↓
Exit	

List item	Function
Sample Volume [ml]	Sample volume to be injected in the analyzer.
Fill Volume [ml]	Sample volume to be aspirated for the syringe prefill and the pull up strokes.
Pull Up Strokes	Up and down movements of the plunger to homogenize the head space sample in the vial.

List item	Function
Equilibration Delay [s]	Delay time between sample pick up and the syringe moving from the sample.
Sample Speed [ml/min]	Speed of the plunger during the syringe filling operation and pull up operation.
Syringe Prefill	<p>Syringe-filling mode.</p> <ul style="list-style-type: none"> • Yes: the syringe plunger has to be pulled up... • No: the syringe plunger has to be pulled down... <p>...when the needle comes in the vial, before the septum has been pierced.</p>
Exit	Exit the page and return to the selected method page.

9.10.5 Injection

✱ *Run menu > Settings > Methods > (method) OK > Injection*



List item	Function
Injection speed [ml/min]	Syringe dispensation speed while injecting the sample in the analyzer.
Pre Injection Dwell [s]	Waiting time of the needle inside the injector before injecting the sample in the analyzer. This allows the syringe temperature to level with the analyzer injector temperature.
Post Injection Dwell [s]	Waiting time of the needle inside the injector after having injected the sample in the analyzer. This allows the complete evaporation of the sample.
Flush Time [min]	Time during which the syringe is flushed with gas after injection.
Exit	Exit the page and return to the selected method page.

10. Screen description in Liquid mode

10.1 Screen elements

Note: The interface may be updated in future. In this case, please follow the screen instructions.




10.1.1 Interacting with the touch screen

The touch screen can be tapped with a finger or a capacitive touch screen stylus.

10.1.2 Types of screens


The screens of the touch screen are mainly of two types: menus, with icons, or pages, with a list.

10.1.3 Common buttons

Part	Function
	Move up and down in a list and select an item.
	<ul style="list-style-type: none"> • Move the needle/plunger up and down. • Open/Close the tray. • Increase and decrease the visualized parameter.
	Move the turret left and right.
OK	<ul style="list-style-type: none"> • Confirm the selection. • Confirm a parameter value. • Commute a signal status.
Save	Save the modifications.
Next	Go to the next page.
Exit	Exit the menu.
Cancel	Exit the page without saving the modifications.
Continue	Confirm and go to the next page.
Abort	Stop the autosampler immediately.
Stop	Pause the running operation.
Menu	Access additional options during the automatic run.
Resume	Continue a paused sequence.
Terminate current	Stop the sequence. The sequence stops only after having completed the injection in progress in the analyzer.


10.2 Main menus

10.2.1 Screen flowchart



10.2.2 Run menu

Quick Start version




Standard version



List item	Function
Start	Only in the quick start version. Process all the sample vials listed in all the sequence steps stored. After two consecutive empty positions, the autosampler stops.
Single Injection	Only in the standard version. Process a single sample vial. For details, see 'Single Injection' on page 277 .
Automatic Run	Only in the standard version. Perform a sequence of injections from different sample vials. For details, see 'Automatic Run' on page 278 .
Syringe Wash	Allow an extra wash of the syringe with the solvent/s. For details, see 'Syringe Wash' on page 278 .
Load Sample/Solvent	Open the tray to load/unload the solvent vials, sample vials and racks.
Settings	Access the Settings menu. For details, see 'Settings' on page 279 .
Injector Septa Test	Check if the syringe can pierce the injector septum easily and correctly.
Generate QR code	Open the Generate QR code wizard on the PC in the HTA Monitor.
Change Mode	Change the operation mode (Headspace or SPME).

10.2.3 Automatic Run menu

 (during the automatic run) Menu






List item	Function
Methods	Show the methods menu.
Sequences	Show the sequences menu.
Priority	Set a priority injection during the automatic run in progress.
Load Sample Tray	Open the tray to add sample vials.
Exit	Exit the menu and return to the summary page.

10.3 The Run menu (main menu)


10.3.1 Single Injection

✱ *Run menu > Single Injection*





List item	Description
Start Injection	Start the single injection set.
Method	From 0 to 9.
Sample	Position of the sample vial to be injected in the analyzer.
Internal Standard	Position of the internal standard vial.
Pre Wash Solvent	Position of the solvent vial/s to be used for the pre-washing operations.
Post Wash Solvent	Position of the solvent vial/s to be used for the post-washing operations.
Exit	Exit the page and return to the Run menu.

10.3.2 Automatic Run


 Run menu > Automatic Run

Run: Automatic Run



Start Autorun	
First Step: 1	
Last Step: 1	OK
Run-All-Samples: Yes	
Exit	

List item	Function
Start Autorun	Start the set sequence.
First Step	Set the first step of the sequence.
Last Step	Set the last step of the sequence.
Run-All-Samples	Set the behaviour of the autosampler when it detects two consecutive empty positions: <ul style="list-style-type: none">• Yes: the autosampler stops.• No: the autosampler searches for other sample vials.
Exit	Exit the page and return to the Run menu.

10.3.3 Syringe Wash


 Run menu > Syringe Wash


Run: Syringe Wash

Start Wash	
Solvent Vial: A	
Wash Volume [µl]: 10.0	OK
Cycle(s): 1	
Exit	

List item	Description
Start Wash	Run the syringe wash.
Solvent Vial	Position of the solvent vial
Wash Volume [µl]	Amount of solvent to be used for each wash cycle
Cycle(s)	Number of wash cycles
Exit	Exit the page and return to the Run menu.

10.3.4 Settings

 Run menu > Settings





List item	Function
Methods	Set the method parameters. For details, see 'Methods' on page 280 .
Sequences	Set the sequence parameters. For details, see 'Sequences' on page 281 .
Tray Type	Set the rack type to be used. For details, see 'Choose the type of rack to be used' on page 172 .
View Parameters	Set the autosampler parameters. For details, see 'View parameters' on page 281 .
Setup	Access the menu for autosampler setup. For details, see 'Setup' on page 282 .
Exit	Exit the menu and return to the Run menu.

10.4 The Settings menu

10.4.1 Methods


✳️ *Run menu > Settings > Methods*



List item	Function
Method	Select the method to be modified. For details, see 'The method menu' on page 300 .
Tools	Manage the method. For details, see 'Tools' on page 280 .
Exit	Exit the page and return to the Settings menu.

10.4.2 Tools

✳️ *Run menu > Settings > Methods > Tools*




List item	Function
Restore Default	Load the default settings for a method.

List item	Function
Copy Methods	Copy the parameters from one method to another.
Exit	Exit the page and return to the Methods menu.

10.4.3 Sequences



 *Run menu > Settings > Sequences*



List item	Function
Step	Select the desired step. For details, see 'Step' on page 288 .
Tools	Manage the sequence steps. For details, see 'Tools' on page 289 .
Exit	Exit the page and return to the Settings menu.

10.4.4 View Parameters


 *Run menu > Settings > View Parameters*






List item	Description
Syringe Volume [µl]	Volume of the installed syringe.
Tray Type	Tray type installed, with number of vial positions, dimension of sample vials and needle draw depth inside the vial.
Analyzer Brand	Brand of the analyzer to which the autosampler is connected.
Analyzer Model	Model of the analyzer to which the autosampler is connected.
Front Injector Depth [mm]	Needle insertion depth into the front injector.
Front Injector Insertion Speed	Needle insertion speed into the front injector.
Rear Injector Depth [mm]	Only with rear injector. Needle insertion depth into the rear injector.
Rear Injector Insertion Speed	Only with rear injector. Needle insertion speed into the rear injector.
IP Address	LAN connection data.
Subnet Mask	
Gateway	
Exit	Exit the page and return to the Settings menu.

10.4.5 Setup

 Run menu > Settings > Setup




List item	Function
Alignment	Access the autosampler alignment procedures to adjust the working positions. For details, see 'Alignment' on page 284 .

List item	Function
Syringe Installation	Install the syringe.
Manual Operation	<p>Move the plunger, needle, turret and tray for diagnostic or emergency needs. Status:</p> <ul style="list-style-type: none"> • RUN: moving • IDLE: not moving • BUSY: an error impedes the movement • NEED ENC ALIGN: the motor has not been aligned yet (only for service purpose) <p>For details, see 'Manual operation' on page 284.</p>
Configuration	<p>Edit the autosampler general settings.</p> <p>For details, see 'Configuration' on page 286.</p>
Service	Password protected. Allow service operations on the autosampler.
Exit	Exit the menu and return to the Run menu.

10.5 The Setup menu

10.5.1 Alignment

✱ *Run menu > Settings > Setup > Alignment*







List item	Function
Analyzer Default	Select the analyzer on which the autosampler is mounted.
Front Injector	Enabled only if the analyzer is selected. Align the vial locator to the front injector.
Rear Injector	Enabled only if the analyzer is selected. Align the vial locator to the rear injector.
Touch Plunger Zero	Perform the autosampler alignment procedure.
Tray Install	List the racks needed for the sequence with the related draw depth (mm) and manage the list.
Exit	Exit the menu and return to the Setup menu.

10.5.2 Manual Operation

✱ *Run menu > Settings > Setup > Manual Operation*

This page is useful for troubleshooting.






List item	Description
Plunger	<ul style="list-style-type: none"> Show the status sent by the plunger motor. If selected, it allows to move the plunger motor.
Needle	<ul style="list-style-type: none"> Show the status sent by the needle motor. If selected, it allows to move the needle motor.
Turret	<ul style="list-style-type: none"> Show the status sent by the turret motor. If selected, it allows to move the turret motor.
Tray	<ul style="list-style-type: none"> Show the status sent by the tray motor. If selected, it allows to move the tray motor.
Touch sensor	Status of the touch sensor: <ul style="list-style-type: none"> On: the vial locator has detected the object. Off: the vial locator has not detected the object.
Cover	<ul style="list-style-type: none"> Show the status sent by cover motor. If selected, it allows to move the cover motor.
Shaker	Allow to select a defined speed and, by pressing OK button, to activate it at the defined speed.


List item	Description
	It is used for service purpose.
Shaker sensor	Status of the shaker sensor: <ul style="list-style-type: none"> • On: the shaker sensor is on. • Off: the shaker sensor is off. It is used for service purpose.
Ready	Signal sent by the analyzer to communicate that it has finished the analysis. <ul style="list-style-type: none"> • On • Off
Start	Signal sent by the autosampler to communicate that it has started the run. <ul style="list-style-type: none"> • On • Off Select OK to switch the signal.
Sync-out	Status of the sync- out signal. Select OK to switch the signal.
Sync-in	Status of the sync- in signal. <ul style="list-style-type: none"> • On • Off
Prep	Signal sent by the autosampler to communicate to quit the PREP status and to go in READY status. <ul style="list-style-type: none"> • On • Off Select OK to switch the signal.
Exit	Exit the page and return to the Setup menu.

10.5.3 Configuration

✱ *Run menu > Settings > Setup > Configuration*

This menu contains general settings of the autosampler.





List item	Function
Sound	Change the sound effects of the autosampler: <ul style="list-style-type: none"> • Level 2: emit a sound when the screen is selected • Level 1: emit a sound in the following cases: <ul style="list-style-type: none"> ○ The single injection is finished. ○ The sample has been injected in the analyzer. ○ Save or Store have been selected. ○ The autosampler memory has been erased. • Level 0: no sound
Light	Change the syringe illumination settings: <ul style="list-style-type: none"> • On: syringe LED on • Off: syringe LED off • Injection: syringe LED on during injection only
Run	Change the settings for the injection sequence. For details, see 'Run' on page 290 .
Input Output	Change the input/output parameters to/from the analyzer. For details, see 'Input Output' on page 292 .
Maintenance	Show information and allow to change the settings of the preventive maintenance counters. For details, see 'Maintenance' on page 293 .
Advanced	Password protected. Edit the advanced parameters of the autosampler.
Activation	Insert the password for the optional features.
BCR	Show information about the sample vial barcode and change the barcode reader settings. For details, see 'BCR' on page 294 .
SyringeID	Show information about the syringe ID and change the SyringeID system settings. For details, see 'SyringeID' on page 295 .
Exit	Exit the menu and return to the Setup menu.

10.6 The Sequences page

10.6.1 Step

✱ *Run menu > Settings > Sequences > Step*

Run: Sequences

Step: 1	↑
Tools	
Exit	OK
	↓

Step: 1

Tray Type: 121 - 2ml	↑
Method: 0	
First Sample: A.1	OK
Last Sample: A.1	
Cycle(s): 1	↓
Internal Standard: A.1	

Step: 1


Internal Standard: A.1	↑
Pre Wash Solvent: C	
Post Wash Solvent: C	OK
Injection Mode: Default	
Name:	↓
Exit	

List item	Description
Tray Type	Type of rack to be used.
Method	Method to be used.
First Sample	Position of the first sample vial to be processed with the selected method.
Last Sample	Position of the last sample vial to be processed with the selected method. All the sample vials from the position set in First Sample to this one will be processed.
Cycle(s)	Number of injections of the same sample in the analyzer.
Internal Standard	Sample vial or solvent vial to be used as the internal standard.
Pre Wash Solvent	Solvent vial/s to be used for the pre- washing operations.
Post Wash Solvent	Solvent vial/s to be used for the post- washing operations.

List item	Description
Injection Mode	<ul style="list-style-type: none"> • Default: the autosampler injects in the injector defined in 'Run' on page 290. • Front: the autosampler injects only in the front injector. • Rear: the autosampler injects only in the rear injector. • Confirmation: the autosampler injects the same sample in both injectors. • High throughput: the autosampler injects different samples alternatively in the front and rear injector.
Name	Name of the sequence step. Editable from HTA Autosampler Manager.
Exit	Exit the page and return to the Sequences page.

10.6.2 Tools


✱ *Run menu > Settings > Sequences > Tools*






List item	Function
Add New	Add a new sequence step.
Copy as New	Create a new sequence step by copying the parameters of an existing ones.
Reorder	Change the order of the sequence steps.
Remove	Remove a sequence step.
Exit	Exit the page and return to the Sequences page.

10.7 The Configuration menu

10.7.1 Run

 *Run menu > Settings > Setup > Configuration > Run*


This page allows to change the autosampler behaviour during the run.




List item	Function
Tray Behaviour	Set if the tray must stay open or closed during the run.
Power-on Restart	Set the autosampler restart after a power failure during the run: <ul style="list-style-type: none">• On: the run restarts.• Off: the run does not restart and the Run menu appears.
Injection Mode	Only if two injectors are present. Set the default injector in which to inject the sample. <ul style="list-style-type: none">• Front: the autosampler injects only in the front injector.• Rear: the autosampler injects only in the rear injector.• Confirmation: the autosampler injects the same sample in both injectors.• High throughput: the autosampler injects different samples alternatively in the

List item	Function
	front and rear injector.
Injection Synchro	Set the autosampler behaviour in relation to the analyzer. For details, see ' Injection Synchro ' on page 296.
Empty Position Synchro	Set the autosampler behaviour when a vial is missing: <ul style="list-style-type: none"> • Yes: the autosampler sends the Start signal to the analyzer anyway. • No: the autosampler does not send the Start signal to the analyzer and proceeds with the next vial.
Maintenance Warning	Set when the maintenance warning messages can appear: <ul style="list-style-type: none"> • Normal: at Start up and at the beginning of each run, if the set limit is reached. • Only at Start up: at the beginning of each run, if the set limit is reached.
User Interface	Set the version of the Run menu: <ul style="list-style-type: none"> • Standard: version with all buttons • Quick Start: version with less buttons
Prep-Ahead Waiting Pos.	Only if Injection Synchro is set to Trigger . Set the position of the autosampler while it is waiting the signal for injecting the sample. <ul style="list-style-type: none"> • In vial: the autosampler waits with the needle in the vial. • On Injector: the autosampler waits with the turret over the injector.
Vial locator type	Set the kind of vial locator installed. Default value: A.
Log	Set if the autosampler must send the logs to HTA Monitor. <ul style="list-style-type: none"> • Enable: the autosampler sends logs to HTA Monitor, • Enable w/o warnings: the autosampler sends logs to HTA Monitor, but does not show warnings in case the PC is not connected to the autosampler. • Disable: the autosampler does not send logs to HTA Monitor.
Log Level	Set the type of logs the autosampler sends to the HTA Monitor. <ul style="list-style-type: none"> • Normal: the autosampler sends the standard logs, • Troubleshooting: the autosampler sends more detailed logs to HTA Monitor.
Protocol version	Set the communication protocol version used by the autosampler. Note: Do not change unless: <ul style="list-style-type: none"> ○ you operate by CDS through HTA software connectors ○ required by your third party software documentation ○ required by your Customer Representative
Exit	Exit the page and return to the Configuration menu.


10.7.2 Input Output

 Run menu > Settings > Setup > Configuration > Input Output





List item	Function
Ready On	<p>Set the ON state of the Ready signal. This signal is used to start the autosampler and must match the signal given by the analyzer.</p> <ul style="list-style-type: none">• Low: the autosampler starts the injection when the analyzer provides a low logical level of voltage for the Ready signal.• High: the autosampler starts the injection when the analyzer provides a high logical level of voltage (+5 V) for the Ready signal.
Start On	<p>Set the ON state of the Start signal from the autosampler to the analyzer during the injection. This parameter has to be set according to the analyzer.</p> <ul style="list-style-type: none">• Open: during the injection, the pins 3 (SAMPINS- NO) and 6 (SAMPINS COM) are open for about one second. They are closed before and after the injection.• Close: during the injection, the pins 3 (SAMPINS- NO) and 6 (SAMPINS COM) are closed for about one second (short- circuit). They are open before and after the injection.
Prep Mode On	<p>Set the ON state of the Prep Mode signal from the autosampler to the analyzer during the run. This parameter has to be set according to the analyzer.</p> <ul style="list-style-type: none">• Open: before the processing of the sample vial, the pins 11 (ENDSAMP - NO) and 14 (ENDSAMP- COM) are open for about one second. If within one minute the Ready signal is not received, this signal is activated again. This operation is repeated until the Ready signal is received.• Close: before the processing of the sample vial, the pins 11 (ENDSAMP - NO) and 14 (ENDSAMP- COM) are closed for about one second (short- circuit). If within one minute the Ready signal is not received, this signal is activated again. This operation is repeated until the Ready signal is received.
Exit	<p>Exit the page and return to the Configuration menu.</p>



10.7.3 Maintenance

 *Run menu > Settings > Setup > Configuration > Maintenance*



Setup: Maintenance

Plunger Strokes Counter: 9	
Plunger Strokes Limit: 25000	
Injection Counter: 13	OK
Injection Limit: Off	
Penetration Counter: 3	
Penetration Limit: 10000	

Setup: Maintenance

Penetration Limit: 10000	
Syringe expire within [months]: ---	
Syringe Expiry [months]: 3	OK
Syringe Testing: Enable	
Current Date: Start PC to use AI features.	
Execute Plunger Baseline	


Setup: Maintenance

Execute Plunger Baseline	
Injector Septa Testing: Enabled	
Turret Testing: Enabled	OK
Execute Turret Baseline	
Execute Turret Testing	
Exit	

List item	Description
Plunger Strokes Counter	Number of plunger strokes since the last reset. It allows to monitor the syringe wear- and- tear. This counter is automatically reset after a new syringe installation.
Plunger Strokes Limit	Limit of plunger strokes. Once reached, the autosampler shows a maintenance warning message. Off: no limit
Injection Counter	Number of injections performed in the analyzer since the last reset. It allows to predict when the analyzer septum or injection liner should be replaced.
Injection Limit	Limit of injection number. Once reached, the autosampler shows a maintenance warning message. Off: no limit
Penetration Counter	Total number of perforation of the solvent and sample vial septa, and of the injector. It allows to predict when the syringe should be replaced, due to

List item	Description
	needle wear.
Penetration Limit	Limit of perforation of the solvent and sample vial septa, and of the injector. Once reached, the autosampler shows a maintenance warning message. Off: no limit
Syringe expiry within [months]	Only if connected to the HTA Monitor. Remaining months to the syringe expiration. It allows to predict when the syringe should be replaced: it is particularly relevant when the autosampler performs a few injection per day.
Syringe Expiry [months]	Only if connected to the HTA Monitor. Month in which the syringe expires. Once reached, the autosampler shows a maintenance warning message. Off: no limit
Syringe Testing	Initializing procedure performed by the autosampler every times it is switched on. It checks the correct movement and the syringe plunger to verify the good functioning of the syringe. Available only for syringes PN 1.23.103, 1.23.104, 1.23.106. <ul style="list-style-type: none"> • Enable • Disable
Current Date	Only if connected to the HTA Monitor. Current Date.
Execute Plunger Baseline	Perform the baseline necessary to perform the Syringe Testing. To be performed only on Customer Representative request.
Injector Septa Testing	Procedure for the syringe PN 1.23.103 to avoid damaging the syringe needle during the injection. <ul style="list-style-type: none"> • Enable • Disable
Turret Testing	Initializing procedure performed by the autosampler every times it is switched on. It checks the correct movement of the turret to verify the good functioning of the autosampler. <ul style="list-style-type: none"> • Enable • Disable
Execute Turret Baseline	Perform the baseline necessary to perform the Turret Testing. To be performed only on Customer Representative request.
Execute Turret Testing	Perform a Turret Testing.
Exit	Exit the page and return to the Configuration menu.

10.7.4 BCR

 *Run menu > Settings > Setup > Configuration > BCR*

List item	Description
Enable	Read the sample vial barcode and save it in the autosampler memory so that the information contained can be read by the available software: <ul style="list-style-type: none"> • Yes • No
On Reading Error	Behaviour if a reading error occurs: <ul style="list-style-type: none"> • Skip: skip the run for the vial.

List item	Description
	<ul style="list-style-type: none"> • Process: perform the run anyway.
Test	Test the alignment of the turret in front of the barcode reader for correct reading.
Exit	Exit the page and return to the Configuration menu.

10.7.5 SyringeID

 *Run menu > Settings > Setup > Configuration > SyringeID*

List item	Description
Enable	Read the syringe ID and save it in the autosampler memory so that the information contained can be read by the available software: <ul style="list-style-type: none"> • Yes • No
Reading	Read the information contained in the RFID: <ul style="list-style-type: none"> • syringe serial number • syringe part number • syringe volume • number of plunger strokes performed • number of penetration performed
Exit	Exit the page and return to the Configuration menu.

10.8 The Run page

10.8.1 Injection Synchro

✱ Run menu > Settings > Setup > Configuration > Run > Injection Synchro

Note: **Prep Mode** is available only for some analyzer. Check in your autosampler mounting kit instruction manual if this feature is supported on your analyzer.

Note: If using **Prep Mode**, it is mandatory to start sample list on analyzer software before to start the run on the autosampler, otherwise the analysis could be not executed and/or not saved properly.



Injection Synchro	Sample preparation start	Sample injection	Start signal to the analyzer
Normal	When the autosampler receives the Ready signal from the analyzer after the end of the analysis time set in the autosampler method.	After the autosampler has checked that the Ready signal is still present.	When the syringe plunger starts to move with the syringe in the injector.
Trigger	When the autosampler receives the Ready signal from the analyzer after the end of the analysis time set in the autosampler method. After ending sample preparation it activates the Sync-Out signal.	when the autosampler receives a Sync-In signal, it starts the injection and deactivates the Sync-Out signal when the Start signal is activated.	
A-Start	When the autosampler receives the Ready signal from the analyzer after the end of the analysis time set in the autosampler method.	After the autosampler has checked that the Ready signal is still present.	When the autosampler touches the injector.
D-Start			When the syringe plunger finishes to move with the syringe in the injector.
EA	When the autosampler receives the Ready signal from the analyzer and the Sync- in signal after the end of the analysis time set in the autosampler method.		When the syringe plunger starts to move with the syringe in the injector.
Normal w/o Rdy	At the end of the analysis time set in the autosampler method	At the end of the sample preparation.	<ul style="list-style-type: none"> • The first Start signal after having injected the flush volume and waited for the viscosity delay. • The second Start signal after having emptied the syringe and waited for the viscosity delay.
D-Start w/Flush	When the autosampler receives the Ready signal from the analyzer after the end of the analysis time set in the autosampler method.	After the autosampler has checked that the Ready signal is still present. The autosampler injects the sample in two phases: <ul style="list-style-type: none"> • first phase: it injects the flush volume • second phase: when Flush Injection dwell is elapsed, it empties the syringe. 	

Injection Synchro	Sample preparation start	Sample injection	Start signal to the analyzer
Prep Mode	At the end of the analysis time set in the autosampler method, the autosampler sends a PREP signal to the analyzer and wait for the Ready signal from the analyzer.	After the autosampler has checked that the Ready signal is still present.	When the syringe plunger starts to move with the syringe in the injector.

10.9 The Method page

10.9.1 Method

✱ Run menu > Settings > Methods > (method) OK



List item	Function
General	Access the page for setting the general parameters of the method. For details, see 'General' on page 300 .
Solvent Wash	Access the page for setting the parameters of the washing solvents of the method. For details, see 'Solvent Wash' on page 300 .
Internal Standard	Access the page for setting the internal standard parameters of the method. For details, see 'Internal Standard' on page 301 .
Sample	Access the page for setting the sample parameters of the method. For details, see 'Sample' on page 301 .
Injection	Access the page for setting the injection parameters of the method. For details, see 'Injection' on page 302 .

List item	Function
Exit	Exit the menu and return to the Methods page.

10.10 The Method menu

10.10.1 General

✱ *Run menu > Settings > Methods > (method) OK > General*

Method 0 General	
Syringe Volume [μl]: 10	↑
Analysis Time [min]: 0	
Fill Speed [μl/s]: 1.0	OK
Viscosity Delay [s]: 5	
Exit	↓

List item	Description
Syringe Volume [μl]	Volume of the syringe.
Analysis Time [min]	Time between the injection and the moment in which the analyzer is ready again.
Fill Speed [μl/s]	Syringe aspiration speed.
Viscosity Delay [s]	Time the needle remains in the sample vial after the plunger has been raised. It allows viscous samples time to fill the syringe.
Exit	Exit the page and return to the selected method page.

10.10.2 Solvent Wash

✱ *Run menu > Settings > Methods > (method) OK > Solvent Wash*

Method 0 Solvent Wash	
Pre Wash Cycle(s): 3	↑
Pre Wash Volume [μl]: 10.0	
Pre Wash Mode: inject	OK
Post Wash Cycle(s): 3	
Post Wash Volume [μl]: 10.0	↓
Post Wash Mode: inject	

Method 0 Solvent Wash	
Pre Wash Volume [μl]: 10.0	↑
Pre Wash Mode: inject	
Post Wash Cycle(s): 3	OK
Post Wash Volume [μl]: 10.0	
Post Wash Mode: inject	↓
Exit	

List item	Description
Pre Wash Cycle(s)	Number of syringe washes before injecting the sample in the analyzer
Pre Wash Volume [μl]	Quantity of solvent to be used for each pre- washing operation
Pre Wash Mode	Frequency of the pre- washing operations: <ul style="list-style-type: none"> • Inject: wash at every injection • Sample: wash for each new vial • Step: wash for each new step of the sequence
Post Wash Cycle(s)	Number of syringe washes after having injected the sample in the analyzer
Post Wash Volume [μl]	Quantity of solvent for each post- washing operation

List item	Description
Post Wash Mode	Frequency of the post- washing operation: <ul style="list-style-type: none"> • Inject: wash at every injection • Sample: wash for each new vial • Step: wash for each new step of the sequence
Exit	Exit the page and return to the selected method page.

10.10.3 Internal Standard

✱ *Run menu > Settings > Methods > (method) OK > Internal Standard*

Method 0 Internal Standard	
Enable: No	↑
Internal Standard Volume [μl]: ---	
Air Gap mode: Post	OK
Air Gap Volume [μl]: ---	
Exit	↓

List item	Function
Enable	Enable or disable the internal standard.
Internal Standard Volume [μl]	Volume of the Internal standard
Air Gap Mode	Manage the air gap in the syringe. <ul style="list-style-type: none"> • Post: one air gap between the internal standard and the sample • Double: two air gaps, one air gap between the internal standard and the sample, the second one before the internal standard
Air Gap Volume [μl]	Volume of the air gap
Exit	Exit the page and return to the selected method page.

Note: See ['Use of an internal standard' on page 35](#)

10.10.4 Sample


✱ *Run menu > Settings > Methods > (method) OK > Sample*

Method 0 Sample	
Sample Volume [μl]: 1.0	↑
Fill Volume [μl]: 5.0	
Air Volume [μl]: 1.0	OK
Sample Washes: 1	
Pull Up Strokes: 1	↓
Exit	

List item	Function
Sample Volume [μl]	Sample volume to be injected in the analyzer.
Fill Volume [μl]	Sample volume to be aspirated for the sample wash and the pull up strokes.
Air Volume [μl]	Quantity of air to be aspirated after the sample has been drawn. This prevents volatile samples evaporating from the needle.
Sample Washes	Number of sample washes from the sample vial to the waste vial.
Pull Up Strokes	Rapid up and down movements of the plunger to eliminate air bubbles from the syringe.
Exit	Exit the page and return to the selected method page.

10.10.5 Injection

✱ *Run menu > Settings > Methods > (method) OK > Injection*



List item	Function
Injection speed [μl/s]	Syringe dispensation speed while injecting the sample in the analyzer.
Pre Injection Dwell [s]	Waiting time of the needle inside the injector before injecting the sample in the analyzer. This allows the syringe temperature to level with the analyzer injector temperature.
Post Injection Dwell [s]	Waiting time of the needle inside the injector after having injected the sample in the analyzer. This allows the complete evaporation of the sample.
Exit	Exit the page and return to the selected method page.

11. Screen description in SPME mode

11.1 Screen elements

Note: The interface may be updated in future. In this case, please follow the screen instructions.




11.1.1 Interacting with the touch screen

The touch screen can be tapped with a finger or a capacitive touch screen stylus.

11.1.2 Types of screens


The screens of the touch screen are mainly of two types: menus, with icons, or pages, with a list.

11.1.3 Common buttons


Part	Function
	Move up and down in a list and select an item.
	<ul style="list-style-type: none"> • Move the needle/plunger up and down. • Open/Close the tray. • Increase and decrease the visualized parameter.
	Move the turret left and right.
OK	<ul style="list-style-type: none"> • Confirm the selection. • Confirm a parameter value. • Commute a signal status.
Save	Save the modifications.
Next	Go to the next page.
Exit	Exit the menu.
Cancel	Exit the page without saving the modifications.
Continue	Confirm and go to the next page.
Abort	Stop the autosampler immediately.
Stop	Pause the running operation.
Menu	Access additional options during the automatic run.
Resume	Continue a paused sequence.
Terminate current	Stop the sequence. The sequence stops only after having completed the injection in progress in the analyzer.

11.2 Main menus


11.2.1 Screen flowchart



11.2.2 Run menu
Quick Start version





Standard version



List item	Function
Start	Only in the quick start version. Process all the sample vials listed in all the sequence steps stored. After two consecutive empty positions, the autosampler stops.
Single Injection	Only in the standard version. Process a single sample vial. For details, see ' Single Injection ' on page 308.
Automatic Run	Only in the standard version. Perform a sequence of injections from different sample vials. For details, see ' Automatic Run ' on page 308.
Fiber Cleaning Conditioning	Allow an extra fiber cleaning or conditioning. For details, see ' Fiber Cleaning Conditioning ' on page 309.
Load Sample Tray	Open the tray to load/unload the sample vials and racks.
Settings	Access the Settings menu. For details, see ' Settings ' on page 309.
Generate QR code	Open the Generate QR code wizard on the PC in the HTA Monitor .
Change Mode	Change the operation mode (Headspace or Liquid).

11.2.3 Automatic Run menu

 (during the automatic run) Menu





List item	Function
Methods	Show the methods menu.
Sequences	Show the sequences menu.
Load Sample Tray	Open the tray to add sample vials.
Immediate Injection	Skip the remaining conditioning time and inject immediately.
Skip Next Vial	Skip the remaining conditioning time and unload the vial, without to perform the injection.
Sample Info	View the remaining conditioning time of the samples in the oven.
Run Info	View important details of the run.
Exit	Exit the menu and return to the summary page.

11.3 The Run menu (main menu)

11.3.1 Single Injection


☀ *Run menu > Single Injection*



List item	Description
Start Injection	Start the single injection set.
Method	From 0 to 9.
Sample	Position of the sample vial to be injected in the analyzer.
Derivatizer	Position of the vial used as derivatizer (if enabled in the method).
Exit	Exit the page and return to the Run menu.

11.3.2 Automatic Run


☀ *Run menu > Automatic Run*



List item	Function
Start Autorun	Start the set sequence.
First Step	Set the first step of the sequence.
Last Step	Set the last step of the sequence.
Run-All-Samples	Set the behaviour of the autosampler when it detects two consecutive empty positions: <ul style="list-style-type: none"> • Yes: the autosampler stops. • No: the autosampler searches for other sample vials.
Exit	Exit the page and return to the Run menu.

11.3.3 Fiber Cleaning Conditioning


✱ *Run menu > Fiber Cleaning Conditioning*



List item	Description
Start	Start the cleaning/conditioning process.
Position	Allow to choose where to expose the fiber.
Remaining time [min]	Allow to choose for how many minutes to expose the fiber.
Exit	Exit the page and return to the selected method page.

11.3.4 Settings

✱ *Run menu > Settings*





List item	Function
Methods	Set the method parameters. For details, see 'Methods' on page 311 .

List item	Function
Sequences	Set the sequence parameters. For details, see 'Sequences' on page 312 .
Tray Type	Set the rack type to be used. For details, see 'Choose the type of rack to be used' on page 209 .
Eco Savings	Set the oven stand by temperature. St-by Oven: allow to set the oven temperature in stand-by conditions.
View Parameters	Set the autosampler parameters. For details, see 'View parameters' on page 312 .
Setup	Access the menu for autosampler setup. For details, see 'Setup' on page 314 .
Exit	Exit the menu and return to the Run menu.

11.4 The Settings menu

11.4.1 Methods


✱ *Run menu > Settings > Methods*



List item	Function
Method	Select the method to be modified. For details, see 'The method menu' on page 330 .
Tools	Manage the method. For details, see 'Tools' on page 311 .
Exit	Exit the page and return to the Settings menu.

11.4.2 Tools

✱ *Run menu > Settings > Methods > Tools*




List item	Function
Restore Default	Load the default settings for a method.

List item	Function
Copy Methods	Copy the parameters from one method to another.
Exit	Exit the page and return to the Methods menu.

11.4.3 Sequences


 *Run menu > Settings > Sequences*



List item	Function
Step	Select the desired step. For details, see 'Step' on page 320 .
Tools	Manage the sequence steps. For details, see 'Tools' on page 321 .
Exit	Exit the page and return to the Settings menu.

11.4.4 View Parameters

 *Run menu > Settings > View Parameters*





Run: View Parameters 2/3	
Front Injector Insertion Speed:	Very low
Rear Injector Depth [mm]:	0
Rear Injector Insertion Speed:	Low
Cleaning Device:	No
Cleaning Device Temperature [C]:	Off
Exit	Next

Run: View Parameters 3/3	
IP Address:	192.168.0.235
Subnet Mask:	255.255.0.0
Gateway:	192.168.0.1
Exit	Next

List item	Description
Fiber Length	Length and type of the installed fiber.
Tray Type	Tray type installed, with number of vial positions, dimension of sample vials and needle draw depth inside the vial.
Analyzer Brand	Brand of the analyzer to which the autosampler is connected.
Analyzer Model	Model of the analyzer to which the autosampler is connected.
Front Injector Depth [mm]	Needle insertion depth into the front injector.
Front Injector Insertion Speed	Needle insertion speed into the front injector.
Rear Injector Depth [mm]	Only with rear injector. Needle insertion depth into the rear injector.
Rear Injector Insertion Speed	Only with rear injector. Needle insertion speed into the rear injector.
Cleaning Device	Cleaning device enabled or not.
Cleaning Device Temperature [C]	Temperature at which the fiber will be exposed in the cleaning device.
IP Address	LAN connection data.
Subnet Mask	
Gateway	
Exit	Exit the page and return to the Settings menu.

11.4.5 Setup

 *Run menu > Settings > Setup*




List item	Function
Alignment	Access the autosampler alignment procedures to adjust the working positions. For details, see 'Alignment' on page 315 .
Syringe Installation	Install the syringe.
Manual Operation	Move the plunger, needle, turret and tray for diagnostic or emergency needs. Status: <ul style="list-style-type: none">• RUN: moving• IDLE: not moving• BUSY: an error impedes the movement• NEED ENC ALIGN: the motor has not been aligned yet (only for service purpose) For details, see 'Manual operation' on page 315 .
Configuration	Edit the autosampler general settings. For details, see 'Configuration' on page 318 .
Service	Password protected. Allow service operations on the autosampler.
Exit	Exit the menu and return to the Run menu.

11.5 The Setup menu

11.5.1 Alignment

✱ *Run menu > Settings > Setup > Alignment*







List item	Function
Analyzer Default	Select the analyzer on which the autosampler is mounted.
Front Injector	Enabled only if the analyzer is selected. Align the vial locator to the front injector.
Rear Injector	Enabled only if the analyzer is selected. Align the vial locator to the rear injector.
Touch Plunger Zero	Perform the autosampler alignment procedure.
Tray Install	List the racks needed for the sequence with the related draw depth (mm) and manage the list.
Exit	Exit the menu and return to the Setup menu.

11.5.2 Manual Operation

✱ *Run menu > Settings > Setup > Manual Operation*

This page is useful for troubleshooting.







List item	Description
Plunger	<ul style="list-style-type: none"> Show the status sent by the plunger motor. If selected, it allows to move the plunger motor.
Needle	<ul style="list-style-type: none"> Show the status sent by the needle motor. If selected, it allows to move the needle motor.
Turret	<ul style="list-style-type: none"> Show the status sent by the turret motor.

List item	Description
	<ul style="list-style-type: none"> If selected, it allows to move the turret motor.
Tray	<ul style="list-style-type: none"> Show the status sent by the tray motor. If selected, it allows to move the tray motor.
Touch sensor	Status of the touch sensor: <ul style="list-style-type: none"> On: the vial locator has detected the object. Off: the vial locator has not detected the object.
Cover	<ul style="list-style-type: none"> Show the status sent by cover motor. If selected, it allows to move the cover motor.
Shaker	Allow to select a defined speed and, by pressing OK button, to activate it at the defined speed. It is used for service purpose.
Shaker sensor	Status of the shaker sensor: <ul style="list-style-type: none"> On: the shaker sensor is on. Off: the shaker sensor is off. It is used for service purpose.
Fan	Allow to select a defined speed and, by pressing OK button, to activate it at the defined speed. It is used for service purpose.
Ready	Signal sent by the analyzer to communicate that it has finished the analysis. <ul style="list-style-type: none"> On Off
Start	Signal sent by the autosampler to communicate that it has started the run. <ul style="list-style-type: none"> On Off Select OK to switch the signal.
Flush	Status of the Flush valve. Select OK to switch the Flush valve status.
Sync-out	Status of the sync- out signal. Select OK to switch the signal.
Sync-in	Status of the sync- in signal. <ul style="list-style-type: none"> On Off
Prep	Signal sent by the autosampler to communicate to quit the PREP status and to go in READY status. <ul style="list-style-type: none"> On Off Select OK to switch the signal.
Exit	Exit the page and return to the Setup menu.

11.5.3 Configuration

 *Run menu > Settings > Setup > Configuration*

This menu contains general settings of the autosampler.



List item	Function
Sound	Change the sound effects of the autosampler: <ul style="list-style-type: none">• Level 2: emit a sound when the screen is selected• Level 1: emit a sound in the following cases:<ul style="list-style-type: none">○ The single injection is finished.○ The sample has been injected in the analyzer.○ Save or Store have been selected.○ The autosampler memory has been erased.• Level 0: no sound
Light	Change the syringe illumination settings: <ul style="list-style-type: none">• On: syringe LED on• Off: syringe LED off• Injection: syringe LED on during injection only
Run	Change the settings for the injection sequence. For details, see 'Run' on page 322 .
Input Output	Change the input/output parameters to/from the analyzer. For details, see 'Input Output' on page 324 .
Maintenance	Show information and allow to change the settings of the preventive maintenance counters. For details, see 'Maintenance' on page 325 .
Advanced	Password protected. Edit the advanced parameters of the autosampler.

List item	Function
Activation	Insert the password for the optional features.
BCR	Show information about the sample vial barcode and change the barcode reader settings. For details, see 'BCR' on page 326 .
Cleaning Device	NA.
Exit	Exit the menu and return to the Setup menu.

11.6 The Sequences page



11.6.1 Step

 *Run menu > Settings > Sequences > Step*

Step: 1


Tray Type: 42 positions	
Method: 0	
First Sample: A.1	OK
Last Sample: A.1	
Derivatizer: A.1	
Desorption Mode: Default	


Step: 1

First Sample: A.1	
Last Sample: A.1	
Derivatizer: A.1	OK
Desorption Mode: Default	
Name:	
Exit	

List item	Description
Tray Type	Type of rack to be used.
Method	Method to be used.
First Sample	Position of the first sample vial to be processed with the selected method.
Last Sample	Position of the last sample vial to be processed with the selected method. All the sample vials from the position set in First Sample to this one will be processed.
Derivatizer	Position of the vial used as derivatizer (if enabled in the method).
Desorption Mode	<ul style="list-style-type: none">• Default: the autosampler injects in the injector defined in 'Run' on page 322.• Front: the autosampler injects only in the front injector.• Rear: the autosampler injects only in the rear injector.
Name	Name of the sequence step. Editable from HTA Autosampler Manager.
Exit	Exit the page and return to the Sequences page.

11.6.2 Tools


 *Run menu > Settings > Sequences > Tools*






List item	Function
Add New	Add a new sequence step.
Copy as New	Create a new sequence step by copying the parameters of an existing ones.
Reorder	Change the order of the sequence steps.
Remove	Remove a sequence step.
Exit	Exit the page and return to the Sequences page.

11.7 The Configuration menu

11.7.1 Run

 *Run menu > Settings > Setup > Configuration > Run*


This page allows to change the autosampler behaviour during the run.




List item	Function
Tray Behaviour	Set if the tray must stay open or closed during the run.
Power-on Restart	Set the autosampler restart after a power failure during the run: <ul style="list-style-type: none">• On: the run restarts.• Off: the run does not restart and the Run menu appears.
Desorption Mode	Only if two injectors are present. Set the default injector in which to inject the sample. <ul style="list-style-type: none">• Front: the autosampler injects only in the front injector.• Rear: the autosampler injects only in the rear injector.
Desorption Synchro	Set the autosampler behaviour in relation to the analyzer. For details, see ' Desorption Synchro ' on page ³²⁷ .

List item	Function
Clean Mode	Only if two injectors are present. Set the default position in which to clean the fiber. <ul style="list-style-type: none"> • Front: the autosampler executes the post-cleaning procedure in the front injector. • Rear: the autosampler executes the post-cleaning procedure in the rear injector.
Shaker Speed	Define the shaker speed used during the run. <ul style="list-style-type: none"> • Very Low • Low • Normal • High • Very High
Maintenance Warning	Set when the maintenance warning messages can appear: <ul style="list-style-type: none"> • Normal: at Start up and at the beginning of each run, if the set limit is reached. • Only at Start up: at the beginning of each run, if the set limit is reached.
User Interface	Set the version of the Run menu: <ul style="list-style-type: none"> • Standard: version with all buttons • Quick Start: version with less buttons
Fiber Protection	Allow to enable or not the shaker during the extraction procedure. <ul style="list-style-type: none"> • No: the shaker is enabled. • Yes: the shaker is disabled.
Log	Set if the autosampler must send the logs to HTA Monitor. <ul style="list-style-type: none"> • Enable: the autosampler sends logs to HTA Monitor, • Enable w/o warnings: the autosampler sends logs to HTA Monitor, but does not show warnings in case the PC is not connected to the autosampler. • Disable: the autosampler does not send logs to HTA Monitor.
Log Level	Set the type of logs the autosampler sends to the HTA Monitor. <ul style="list-style-type: none"> • Normal: the autosampler sends the standard logs, • Troubleshooting: the autosampler sends more detailed logs to HTA Monitor.
Protocol version	Set the communication protocol version used by the autosampler. Note: Do not change unless: <ul style="list-style-type: none"> ○ you operate by CDS through HTA software connectors ○ required by your third party software documentation ○ required by your Customer Representative
Exit	Exit the page and return to the Configuration menu.


11.7.2 Input Output

 *Run menu > Settings > Setup > Configuration > Input Output*



List item	Function
Ready On	<p>Set the ON state of the Ready signal. This signal is used to start the autosampler and must match the signal given by the analyzer.</p> <ul style="list-style-type: none">• Low: the autosampler starts the injection when the analyzer provides a low logical level of voltage for the Ready signal.• High: the autosampler starts the injection when the analyzer provides a high logical level of voltage (+5 V) for the Ready signal.
Start On	<p>Set the ON state of the Start signal from the autosampler to the analyzer during the injection. This parameter has to be set according to the analyzer.</p> <ul style="list-style-type: none">• Open: during the injection, the pins 3 (SAMPINS- NO) and 6 (SAMPINS COM) are open for about one second. They are closed before and after the injection.• Close: during the injection, the pins 3 (SAMPINS- NO) and 6 (SAMPINS COM) are closed for about one second (short- circuit). They are open before and after the injection.
Prep Mode On	<p>Set the ON state of the Prep Mode signal from the autosampler to the analyzer during the run. This parameter has to be set according to the analyzer.</p> <ul style="list-style-type: none">• Open: before the processing of the sample vial, the pins 11 (ENDSAMP - NO) and 14 (ENDSAMP- COM) are open for about one second. If within one minute the Ready signal is not received, this signal is activated again. This operation is repeated until the Ready signal is received.• Close: before the processing of the sample vial, the pins 11 (ENDSAMP - NO) and 14 (ENDSAMP- COM) are closed for about one second (short- circuit). If within one minute the Ready signal is not received, this signal is activated again. This operation is repeated until the Ready signal is received.
Exit	<p>Exit the page and return to the Configuration menu.</p>

11.7.3 Maintenance


 *Run menu > Settings > Setup > Configuration > Maintenance*

Setup: Maintenance	
Injection Counter: 16	↑
Injection Limit: Off	
Penetration Counter: 66	OK
Penetration Limit: 100000	
Current Date: Start PC to use AI features.	↓
Turret Testing: Disabled	

Setup: Maintenance	
Penetration Limit: 100000	↑
Current Date: Start PC to use AI features.	
Turret Testing: Disabled	OK
Execute Turret Baseline	
Execute Turret Testing	↓
Exit	

List item	Description
Injection Counter	Number of injections performed in the analyzer since the last reset. It allows to predict when the analyzer septum or injection liner should be replaced.
Injection Limit	Limit of injection number. Once reached, the autosampler shows a maintenance warning message. Off : no limit
Penetration Counter	Total number of perforation of the sample vial septa, and of the injector. It allows to predict when the syringe should be replaced, due to needle wear.
Penetration Limit	Limit of perforation of the sample vial septa, and of the injector. Once reached, the autosampler shows a maintenance warning message. Off : no limit
Current Date	Only if connected to the HTA Monitor. Current Date.
Turret Testing	Initializing procedure performed by the autosampler every times it is switched on. It checks the correct movement of the turret to verify the good functioning of the autosampler. <ul style="list-style-type: none"> • Enable • Disable
Execute Turret Baseline	Perform the baseline necessary to perform the Turret Testing. To be performed only on Customer Representative request.
Execute Turret Testing	Perform a Turret Testing.
Exit	Exit the page and return to the Configuration menu.

11.7.4 BCR

 *Run menu > Settings > Setup > Configuration > BCR*

List item	Description
Enable	Read the sample vial barcode and save it in the autosampler memory so that the information contained can be read by the available software: <ul style="list-style-type: none"> • Yes • No
On Reading Error	Behaviour if a reading error occurs: <ul style="list-style-type: none"> • Skip: skip the run for the vial. • Process: perform the run anyway.
Test	Test the alignment of the turret in front of the barcode reader for correct reading.
Exit	Exit the page and return to the Configuration menu.

11.8 The Run page

11.8.1 Desorption Synchro

✱ Run menu > Settings > Setup > Configuration > Run > Desorption Synchro

Note: Prep Mode is available only for some analyzer. Check in your autosampler mounting kit instruction manual if this feature is supported on your analyzer.




Note: If using **Prep Mode**, it is mandatory to start sample list on analyzer software before to start the run on the autosampler, otherwise the analysis could be not executed and/or not saved properly.

Injection Synchro	Sample preparation start	Sample injection	Start signal to the analyzer
Normal	At the reception of the Ready signal from the analyzer.	After ending sample preparation the autosampler checks again the presence of the Ready signal from the analyzer and then starts the injection.	At the beginning of the plunger movement.
A-Start			When the autosampler touches the injector.
D-Start			When needle has been completely removed from the injector, after the desorption.
Normal w/o Rdy	At the end of the analysis time set in the autosampler method.	After the ending sample preparation.	At the beginning of the plunger movement.
Prep Mode	At the end of the analysis time set in the autosampler method, the autosampler sends a PREP signal to the analyzer and wait for the Ready signal from the analyzer.	After the autosampler has checked that the Ready signal is still present.	At the beginning of the plunger movement.

11.9 The Method page

11.9.1 Method

✱ Run menu > Settings > Methods > (method) OK



List item	Function
General	Access the page for setting the general parameters of the method. For details, see 'General' on page 330 .
Conditioning	Access the page for setting the advanced parameters of the method. For details, see 'Conditioning' on page 330 .
Wash&Clean	Access the page for setting the preparation parameters of the method. For details, see 'Wash&Clean' on page 331 .
Extraction	Access the page for setting the sample parameters of the method. For details, see 'Extraction' on page 331 .
Desorption	Access the page for setting the injection parameters of the method. For details, see 'Desorption' on page 332 .

List item	Function
Exit	Exit the menu and return to the Methods page.

11.10 The Method menu

11.10.1 General

✱ *Run menu > Settings > Methods > (method) OK > General*

Method 0 General	
Vial Type: 20ml	↑
Analysis Time [min]: 1	
Preparation Time: 00:00:49	OK
Exit	
	↓

List item	Description
Vial Type	Type of the vial to be used.
Analysis Time [min]	Time between the injection and the moment in which the analyzer is ready again.
Preparation Time	Display how long the autosampler needs to perform a cycle. It is calculated from the instrument and is not editable (Time needed to transfer the vial from the tray to the incubator, filling time, injection time, etc.).
Exit	Exit the page and return to the selected method page.

11.10.2 Conditioning

✱ *Run menu > Settings > Methods > (method) OK > Conditioning*

Method 0 Conditioning	
Oven Temperature [C]: Off	↑
Incubation Time [min]: 0	
Shaker On [min]: 0.1	OK
Shaker Off [min]: 0.1	
Exit	↓

List item	Function
Oven Temperature [°C]	Temperature of oven conditioning. This value can be between 40°C and 170°C, or Off (inactive conditioning).
Incubation Time [min]	Sample incubation time inside the oven.
Shaker On [min]	Interval time during incubation time in which the sample is shaken.
Shaker Off [min]	Interval time during incubation time in which the sample is not shaken.
Exit	Exit the page and return to the selected method page.

11.10.3 Wash&Clean

✱ *Run menu > Settings > Methods > (method) OK > Wash&Clean*

Method 0 Wash & Clean	
Pre-Washing: No	↑
Pre-Washing Time [min]: 0.2	
Post-Cleaning: Yes	OK
Post-Cleaning Time [min]: 0.2	
Exit	↓

List item	Function
Pre-Washing	Enable/disable the pre-washing operation. The fiber pre-washing is performed after the extraction (before the desorption) by immersion inside a proper solution in the wash vial.
Pre-Washing Time [min]	Set the time of the pre-washing operation.
Post-Cleaning	Enable/disable the post-cleaning operation. The fiber post-cleaning is the fiber cleaning after the desorption operation. The cleaning is performed in the location specified as ' Clean mode ' in the Setup Configuration > Run .
Post-Cleaning Time [min]	Set the time of the post-cleaning operation.
Exit	Exit the page and return to the selected method page.

11.10.4 Extraction

✱ *Run menu > Settings > Methods > (method) OK > Extraction*


Method 0 Extraction	
Extraction Time [min]: 0.2	↑
Extraction Depth [mm]: 15	
Derivatization: None	OK
Derivatization Time [min]: 0.2	
Derivatization Depth [mm]: 5	↓
Exit	

List item	Function
Extraction Time [min]	Time the fiber is exposed inside the vial to absorb the analytes.
Extraction Depth [mm]	The depth of the fiber inside the vial during the extraction phase.
Derivatization	It can have the values: <ul style="list-style-type: none"> • Pre: derivatization before the extraction. • Post: derivatization after the extraction. • None: no derivatization.
Derivatization Time	Time the fiber is exposed to the derivatization agent.

List item	Function
[min]	
Derivatization Depth [mm]	Depth of the fiber inside the vial during the derivatization.
Exit	Exit the page and return to the selected method page.

11.10.5 Desorption

✱ Run menu > Settings > Methods > (method) OK > Desorption



List item	Function
Desorption [min]	Time the fiber remains in the injection port to desorpt the analytes.
Exit	Exit the page and return to the selected method page.

12. Maintenance

12.1 Preventive maintenance

12.1.1 Basic concepts

Maintenance requires the substitution of consumables that are subject to wear or stress (syringe, solvent vial septa, analyzer injector septa, injection liner). For this reason, the autosampler has three types of preventive maintenance features:

- **Maintenance counters:** allow to set the maximum number of times in which each consumable can be used and provides counters for each of them. Once the set limit has been reached, a warning message appears on the screen to remind that the consumable needs to be substituted. Counters are automatically reset when the warning message is cleared.
- **Self-diagnosis tests:** allow to verify the status of components subject to wear. Two self-diagnosis tests are available:

Self-diagnosis tests	Function	When is performed if enabled
Syringe testing	Test the syringe and in particular the syringe plunger. Available only in Liquid mode and only for syringes PN 1.23.103, 1.23.104.	<ul style="list-style-type: none"> • At each autosampler switch on • At each new syringe installation
Turret testing	Test the turret.	<ul style="list-style-type: none"> • At each autosampler switch on

- **Additional tests:**

Additional tests	Function	When is performed if enabled
Injector septa test	Available only in Headspace and Liquid mode and only for syringes PN 1.23.103, 1.23.104. Allow to prevent damaging of the syringe needle during the injection.	When user executes the test: it is suggested to execute it after installing a new injector septum.

12.2 Maintain the autosampler

12.2.1 Test the injector septum perforation

This procedure needs to be performed after installing a new injector septum to avoid damaging the syringe needle during the injection. It is highly important to execute this test with the injector at room temperature and without pressure.

1. In the **Run** menu, select **Injector Septa Test**: the autosampler verifies the capacity of the syringe of perforating the injector septum. If **Injector Septa Test** is not available, select **Settings > Setup > Configuration > Maintenance** and enable it.
2. Select the injector(s) where to perform the test.
3. Select **Continue**.
4. Do what follows:

If...	then...
the test passes	select Continue .
the test fails	<ol style="list-style-type: none"> 1. Slightly unscrew the injector septum/septa. 2. Select Retry to repeat the Injector septa test. 3. If the problem persist, change the injector septa/septum and retry the Injector septa test. 4. If the problem persist, contact your Customer Representative.

12.2.2 Set the preventive maintenance counters

✱ *Run menu > Settings > Setup > Configuration > Maintenance*

Select the desired counters and set their limit values.

12.2.3 Reset a counter value manually

✱ *Run menu > Settings > Setup > Configuration > Maintenance*

Select the desired counter and select **OK**: a confirmation message appears. Select **Erase > Continue** to confirm.


12.2.4 Replace septa and caps of solvent and waste vials

Change septa and caps on periodical basis. The septa should be replaced more frequently than the caps.

12.3 Clean the autosampler

12.3.1 Clean the external surfaces of the autosampler

1. Switch the autosampler OFF.
2. Disconnect the power cord from the rear.
3. Dampen a soft, lint- free or microfiber cloth with distilled water.
NOTICE: *The cloth must not be wet. Do not use solvents nor detergents.*
4. Wipe the touch screen and all fragile parts of the autosampler gently.
5. Wipe the antifriction sticker to remove any residual.



6. Dampen another soft, lint-free or microfiber cloth with a non- abrasive detergent solution and clean the other external parts of the autosampler.
NOTICE: *Do not use water, solvents, detergents or spray on electrical parts.*
Note: *Alternatively, it is also possible to clean the touch screen with a screen cleaner kit that includes antistatic wipes.*
7. Dry the autosampler with a dry lint-free cloth.
8. Switch the autosampler ON.

12.3.2 Clean the rack and the tray

1. Open the tray and remove the rack.
2. Clean the rack and the tray with a non- abrasive detergent solution.
3. Dry the rack and the tray with a dry lint- free cloth.
4. Reposition the rack inside the tray and close the tray.

12.4 Autosampler packing


12.4.1 Uninstall the autosampler

1. Remove the syringe and the solvent vials, the syringe integrity tool, the external pressure regulator and the swagelok adapter from the autosampler.
2. Open the tray and remove the sample standard rack.
3. Close the tray.
4. Clean the autosampler and its components and refer to the applicable safety regulations to decontaminate them from any residual hazardous substances.
5. Loosen the safety lock.
6. Raise the vial locator completely.
7. Tighten the safety lock.
8. Lower the sliding lid.
9. Raise the syringe holder to the top using the needle motor.
10. Switch the autosampler off.
11. Disconnect the power cord from the wall socket and the power supply from the power supply socket of the autosampler.
12. Disconnect the interface cable from the autosampler port and from the analyzer.
13. Disconnect the Ethernet cable from the autosampler and from the PC or of wall jack.
14. If present, disconnect the purge gas from the connector from the autosampler rear panel.
15. Open the locking hook of the mounting plate.
16. Lift the autosampler and place it on a flat surface in a vertical position. See '[Lift the autosampler](#)' on [page 24](#).
17. To proceed, see '[Pack the autosampler for shipping](#)' on [page 337](#).

12.4.2 Pack the autosampler for shipping


Note: The autosampler can only be shipped in its original packaging.

1. Repack the accessories in the accessory boxes. Seal the accessory boxes with a suitable adhesive tape.
2. Reposition the turret protection.





3. Reposition the polythene bag and the protection belt.


NOTICE: Lift up the sliding lid before to apply the protection belt. The protection belt must not be applied on the sliding lid.




4. Position the autosampler on the cardboard, then reposition the lower package cushioning.




5. Reposition the sealed accessory boxes to the turret sides.




6. Reposition the upper package cushioning.





7. Close the packaging with a suitable adhesive tape.



8. Reposition the autosampler in its original box.



9. Reposition the free package cushioning.



10. Close and seal the box with the suitable adhesive tape.

13. Troubleshooting and errors in Headspace mode

13.1 Troubleshooting

13.1.1 Analyzer reproducibility issue

Note: the indications given below are for gas chromatography applications and may not be suitable for other analysis techniques.

Problem	Cause	What to do
Peak area or peak retention times are not reproducible.	Injection depth	Ensure the needle is set to the correct depth for the analyzer. This is particularly important if the unit has just been installed or moved.
	Leak in the inlet septum	<ol style="list-style-type: none"> 1. If there is a leak in the inlet septum, change it. 2. If the inlet septum has experienced less than 100 injections, ensure what follows to avoid premature septum wear: <ul style="list-style-type: none"> ○ the syringe needle is straight ○ the septum retainer nut is not too tight ○ the syringe has been installed correctly
	Syringe dirty or worn	<ol style="list-style-type: none"> 1. Check the pressure at the external pressure regulator; 2. Check gas flow at needle tip during the flush phase. 3. Increase the syringe flush time in the method. 4. Set the syringe temperature 10°C above the oven temperature. 5. Replace syringe glass and syringe plunger.
	Vacuum created in sample vial	Reduce sample volume.
	Method parameters are not correct	Check each parameter of the method to be sure they fits your application.
Contamination or ghost peak issue	Purge gas is dirty	Check purge gas for impurities.
Peak size issue	Injection depth	Ensure the needle is set to the correct depth for the analyzer. This is particularly important if the unit has just been installed or moved.
	Method parameters are not correct	Check each parameter of the method to be sure they fits your application.
	Leak in the analyzer system	<ol style="list-style-type: none"> 1. Replace the injector septum and check the fitting to avoid leaks. 2. If the septum has experienced less than 100 injections, ensure what follows to avoid premature septum problems: <ul style="list-style-type: none"> ○ the syringe needle is straight ○ the septum retainer nut is not too tight ○ the syringe has been installed correctly
	GC injection liner not suitable	Check the liner geometry according to the recommendations of the GC manufacturer.
	GC related parameters	Any parameter like injector, oven or detector temperature can contribute to bad peak shape. Check the column type and connections (square cut connections). Check the troubleshooting guide of the GC manufacturer.

Problem	Cause	What to do
Carry over issue	Syringe dirty or worn	<ol style="list-style-type: none"> 1. Check the pressure at the external pressure regulator; 2. Check gas flow at needle tip during the flush phase. 3. Increase the syringe flush time in the method. 4. Set the syringe temperature 10°C above the oven temperature. 5. Replace syringe glass and syringe plunger.
No signal/no peak	Syringe is not working properly	Replace syringe glass and syringe plunger.
	Injection depth too high	Ensure the needle goes deep enough into the analyzer injector.
	Injection in the wrong injector	Select the correct injection mode in the Setup menu and Sequences page.

13.1.2 Communication with analyzer issues

Problem	Cause	What to do
The autosampler starts also if the analyzer is "Not Ready".	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The autosampler doesn't start.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The analyzer doesn't start the analysis when the autosampler injects.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The analyzer starts the analysis before that autosampler injects.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.

13.1.3 Usage

Below is a list of possible problems happening and their potential solutions. Some of these are hypothetical as these problems have not occurred.

Problem	Cause	What to do
During the instrument research procedure, no instrument appears on the list	Firewall is blocking HTA Monitor	<p>Be sure that the firewall is allowing HTA Monitor to work. In details, be sure that HTA Monitor could communicate on both Public and Private net.</p> <p>For example, with Windows 10:</p> <ol style="list-style-type: none"> 1. Reach Windows Defender Firewall (Control Panel) 2. Select Allow an application... 3. Select both Public and Private for HTA Monitor app

Problem	Cause	What to do
		4. Retry

13.1.4 Configuration issues

Problem	Cause	What to do
Bad alignment on the injector position(s)	Injector adapter is missing.	Install the proper injector adapter.
	Injector alignment is needed.	Execute Injector alignment . See 'Realign injectors' on page 145 .
The autosampler injects in the wrong injector	Injection Mode in Step is wrong.	In the Run menu, select Settings > Sequences > Step to set the correct injector. See 'Step' on page 259 .
	Injection Mode in Run is wrong.	In the Run menu, select Settings > Setup > Configuration > Run to set the correct injector. See 'Run' on page 261 .
Missing: front injector alignment	The front injector alignment has not been performed after having changed the analyzer brand or model.	In the Run menu, select Settings > Setup > Alignment > Front Injector to align the autosampler to the front injector. See 'Realign injectors' on page 145 .
Missing: rear injector alignment	The rear injector alignment has not been performed after having changed the analyzer brand or model.	In the Run menu, select Settings > Setup > Alignment > Rear Injector to align the autosampler to the rear injector. See 'Realign injectors' on page 145 .
Missing: touch & plunger zero	The Touch & Plunger Zero procedure has not been performed.	In the Run menu, select Settings > Setup > Alignment > Touch & Plunger Zero to perform the alignment procedure. See 'Reset the position of the vial' on page 147 .
WARNING Check Touch Sensor then execute Touch & Plunger Zero.	A syringe with a longer needle is installed and touch and plunger zero procedure is not performed.	1. Select Continue . 2. Select Touch Plunger Zero . See 'Reset the position of the vial' on page 147 .

13.1.5 Consumables/syringe issues

Element	Problem	Cause	What to do
Syringe	Syringe needle damaged.	Use of incorrect vials or injection septum or incorrect syringe installation/configuration.	<ol style="list-style-type: none"> 1. Identify the cause of the damage and the position in which it occurs. 2. Ensure the alignment on the oven cover and on the injector/s is correct. 3. Check that the septa of sample vials and injectors are suitable and correctly pierced. 4. Check that there are no obstacles in the vials, such as an insert.

Element	Problem	Cause	What to do
			5. Replace the syringe glass. See 'Syringe management' on page 131 .
	Needle breaking.	Needle not set up correctly.	Align vial locator to needle tip (see 'Align the vial locator to the needle tip' on page 139) and execute Touch & Plunger Zero (see 'Reset the position of the vial' on page 147).
	Needle obstructed.	Pieces of septum inside needle.	Check quality of septum of sample vial and injector.
		Pieces of septum inside needle.	Check that sample vials are not overcrimped: this may result in septum damaging.
		Sample full of dirt.	Check sample quality or modify tray depth to draw at higher level.
Sample vials	Unit picks up and drops vials.	Safety lock and/or needle height regulator aren't installed properly.	Align the vial locator to the needle tip and execute Touch & Plunger Zero . See 'Reset the position of the vial' on page 147 .
		Vial locator and its two steel bars can't move along the syringe needle.	Check if the safety lock and the needle height regulator are fixed correctly; they should not be fixed in a non-parallel position.
		Septum is too tight.	Check if crimping or screwing is performed correctly.

13.2 Start up errors

13.2.1 Generic start-up errors

Error message	Cause	Solution
<ul style="list-style-type: none"> Plunger error: Code 1 1 Plunger error: Code 1 2 Plunger error: Code 1 4 Plunger error: Code 2 1 Plunger error: Code 2 2 Plunger error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>Syringe is damaged or not installed properly.</p>	<p>Remove the obstacle and press Retry.</p> <p>Remove the syringe and press Retry; then, execute the syringe installation procedure properly. See 'Syringe management' on page 131.</p>
<ul style="list-style-type: none"> Needle error: Code 1 1 Needle error: Code 1 2 Needle error: Code 1 4 Needle error: Code 2 1 Needle error: Code 2 2 Needle error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>Instrument is in shipping position.</p>	<p>Remove the obstacle and press Retry.</p> <p>Press Disengage, set the autosampler in operating position and press Exit.</p>
<ul style="list-style-type: none"> Turret error: Code 1 1 Turret error: Code 1 2 Turret error: Code 1 4 Turret error: Code 2 1 Turret error: Code 2 2 Turret error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p>	<p>Remove the obstacle and press Retry.</p>
<ul style="list-style-type: none"> Tray error: Code 1 1 Tray error: Code 1 2 Tray error: Code 1 4 Tray error: Code 2 1 Tray error: Code 2 2 Tray error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>A sample vial fell between frontal panel and rack support.</p> <p>The rack is not positioned properly.</p>	<p>Remove the obstacle and press Retry.</p> <p>Press Disengage to open the tray and remove the sample vial; press Exit to retry the start-up.</p> <p>Press Disengage to open the tray and reposition the rack correctly; press Exit to retry the start-up.</p>
Unit power up in Service mode: a PIN code is required.	Problem during instrument start-up.	Switch off and on the autosampler to retry the start-up.
Error: "Timeout communication from peripheral: plunger / needle / tray / turret".	Problem during instrument start-up.	Switch off and on the autosampler to retry the start-up.

13.2.2 Safety lock errors

Error message	Cause	Solution
Is safety lock in operating position? If NO, switch off and check user manual.	The syringe location has not been prepared.	<ol style="list-style-type: none"> Switch the autosampler off. See 'Prepare the syringe location' on page 62.

Error message	Cause	Solution
Is safety lock in operating position? If NO, switch off and check user manual.	The touch sensor is not working correctly.	<ol style="list-style-type: none"> 1. Ensure to have prepared the syringe location properly. See 'Prepare the syringe location' on page 62. 2. If the problem persists, select Yes and then Run menu > Settings > Setup > Manual Operation. 3. Select the needle and raise it. 4. Ensure that the vial locator can move up and down freely and that the safety lock is higher than the needle height regulator. 5. Select Exit.

13.2.3 Vial locator errors

Error message	Cause	Solution
Attention! Vial locator in critical position! Press Disengage and move the vial locator to a safe position, then press Exit.	Unexpected switching off of the autosampler while the vial locator was picking up/replacing a vial.	<ol style="list-style-type: none"> 1. Select Disengage: the Manual Operation page appears. 2. Select the needle and raise it. 3. If present, remove the vial from the vial locator and place it back to its original position. 4. Check that all the vials are correctly positioned in the rack. 5. Select Exit.

13.2.4 Memory errors

Error message	Cause	Solution
Error Reading Border mark-out!	The autosampler cannot read the Border Mark Out parameters.	<ol style="list-style-type: none"> 1. Select Continue: the Service page appears. 2. Contact your Customer Representative to solve the problem.
Error Reading Methods!	The autosampler cannot read the parameters of the saved methods.	<ol style="list-style-type: none"> 1. Select Continue. 2. Switch the autosampler off and on. 3. If the error message appears again, contact your Customer Representative to solve the problem.
Error Reading Sequence!	The autosampler cannot read the parameters of the saved sequences.	
Error Reading Features!	The autosampler cannot read the features.	<ol style="list-style-type: none"> 1. Select Continue. 2. Switch the autosampler off and on. 3. If the error message appears again, contact your Customer Representative to solve the problem.

13.3 Run issues / errors

13.3.1 Generic run errors

Error message	Cause	Solution
<ul style="list-style-type: none"> Plunger error: Code 4 2 Plunger error: Code 4 4 Plunger error: Code 4 8 Plunger error: Code 10 7 	Plunger locker is in incorrect position.	Install correctly the plunger locker and execute the Syringe Installation procedure. See 'Syringe management' on page 131 .
	The syringe was installed without to perform completely the Syringe Installation procedure.	Execute completely the Syringe Installation procedure. See 'Syringe management' on page 131 .
	Syringe is not installed correctly.	Install correctly the syringe and perform Syringe Installation procedure. See 'Syringe management' on page 131 .
	Damaged syringe or worn out syringe.	Change the syringe and perform Syringe Installation procedure. See 'Syringe management' on page 131 .
<ul style="list-style-type: none"> Needle error: Code 4 2 Needle error: Code 4 4 Needle error: Code 4 8 Needle error: Code 10 7 	Safety lock or needle height regulator loosen/not fasten properly.	Fasten completely the loosen part and press Retry ; to align the vial locator at the needle tip may be required. See 'Align the vial locator to the needle tip' on page 139 .
	Injector depth too high, not reachable with the syringe needle in use.	Change injector depth according to the syringe needle in use or change the syringe with a syringe with longer needle.
<ul style="list-style-type: none"> Turret error: Code 4 2 Turret error: Code 4 4 Turret error: Code 4 8 Turret error: Code 10 7 	An obstacle does not allow the motor to run for all its possible movement.	Remove the obstacle and press Retry .
	The anti-friction sticker needs to be cleaned.	Wipe the anti-friction sticker with a soft, lint-free or microfiber cloth dampen with distilled water.
<ul style="list-style-type: none"> Tray error: Code 4 2 Tray error: Code 4 4 Tray error: Code 4 8 Tray error: Code 10 7 	A sample vial fell between frontal panel and rack support.	Press Disengage to open the tray and remove the sample vial.
	The rack is not positioned properly.	Press Disengage to open the tray and correct rack position.
<ul style="list-style-type: none"> Cover error: Code 7 1 	A vial has fallen during vial transport and blocks the cover to move.	Remove the vial, and check 'Vial transport problem' on page 349 .
<ul style="list-style-type: none"> Shaker error: Code 7 1 	Something is blocking the shaker to move.	Remove the obstacle and press Retry .
<ul style="list-style-type: none"> Oven error: Code 8 1 Oven error: Code 8 2 Oven error: Code 8 3 Oven error: Code 8 4 	The oven cannot reach or keep the desired temperature.	Contact your customer representative.
<ul style="list-style-type: none"> Syringe error: Code 8 1 Syringe error: Code 8 2 Syringe error: Code 8 3 Syringe error: Code 8 4 	The syringe warmer is not installed or not installed properly.	Install the syringe warmer, if already installed try to remove it and install it again.
	The syringe cannot reach or keep the desired temperature.	Contact your customer representative.

13.3.2 Touch sensor error

Error message	Cause	Solution
WARNING Check Touch Sensor!	The touch sensor is off and/or it is not working correctly.	<ol style="list-style-type: none"> 1. Select Disengage: the Manual Operation page appears. 2. Select the needle and raise it. 3. Ensure that the vial locator can move up and down freely and that the safety lock is higher than the needle height regulator. 4. Select Exit.

13.3.3 Vial transport problem

Error	Cause	Solution
Falling vials	Sample vials not prepared properly.	<ol style="list-style-type: none"> 1. Check the external status of vial, cap and septa. 2. Be sure they are not overcrimped.
	Vial locator damaged.	Contact customer representative.

13.3.4 Vial detection problem

Error message	Cause	Solution						
Warning: Obstacle found over [x]! Try to remove the obstacle	The autosampler found an obstacle over the indicated position.	<div>1. Remove the obstacle.</div> <div>2. To open the Manual Operation page, select Disengage.</div> <div>3. Select the desired item and move its motor and unlock the system.</div> <div>4. To restart the automatic run from the point of interruption select Retry. Otherwise, select Exit.</div> <div>5. If error persists, execute Touch & Plunger Zero. See 'Reset the position of the vial' on page 147.</div>						
Warning: [x] missing! Insert the vial, then press Retry	Missing vial in the indicated position.	<div>Do what follows:</div> <table><tr><th>If you...</th><th>Then...</th></tr><tr><td>can access the desired position</td><td><div>1. Insert the missing vial.</div><div>2. Select Retry.</div></td></tr><tr><td>cannot access the desired position</td><td><div>1. To open the Manual Operation page, select Disengage.</div><div>2. Select the desired item, move its motor and unlock the system.</div><div>3. Insert the missing vial.</div><div>4. Select Exit: automatic run restarts from the point of interruption.</div></td></tr></table> <div>If error persists, execute Touch & Plunger Zero. See 'Reset the position of the vial' on page 147.</div>	If you...	Then...	can access the desired position	<div>1. Insert the missing vial.</div> <div>2. Select Retry.</div>	cannot access the desired position	<div>1. To open the Manual Operation page, select Disengage.</div> <div>2. Select the desired item, move its motor and unlock the system.</div> <div>3. Insert the missing vial.</div> <div>4. Select Exit: automatic run restarts from the point of interruption.</div>
If you...	Then...							
can access the desired position	<div>1. Insert the missing vial.</div> <div>2. Select Retry.</div>							
cannot access the desired position	<div>1. To open the Manual Operation page, select Disengage.</div> <div>2. Select the desired item, move its motor and unlock the system.</div> <div>3. Insert the missing vial.</div> <div>4. Select Exit: automatic run restarts from the point of interruption.</div>							

Error message	Cause	Solution
Bad alignment on the sample vials OR oven cover OR oven locations: <ul style="list-style-type: none"> • autosampler touches more times the vial before recognizing it • Obstacle Found error appears 	Touch & Plunger Zero is needed.	Execute Touch & Plunger Zero . See ' Reset the position of the vial' on page 147 .

13.4 System Integrity test errors

Error	Cause	Solution
Integrity test has failed. Check for clogged needle.	Needle is clogged	Replace the syringe glass.
Integrity test has failed. Check for syringe or System Integrity Tool.	System Integrity Tool has not been placed in its site	Place the System Integrity Tool.
	Syringe damaged	Change the syringe plunger (more probable) or syringe glass.
	System Integrity Tool with damaged septum	Replace system integrity septum.

13.5 Mismatch errors

Error message	Cause	Solution								
Change [x] parameter. Unknown value.	Inconsistency between the parameter and the Protocol version set.	Contact your Customer Representative.								
Warning: Run Error parameters Step [i]: Set Meth [x] Anatetime>0 or change Injection Synchro	Inconsistency between the Analysis Time and the Injection Synchro set.	1. Select Abort . 2. In the Run menu select Settings > Methods and modify the methods setting Analysis Time to a value higher than 0.								
Warning: Run Error parameters Step [x]: Injection Mode mismatch	Inconsistency between the injection mode set in Settings > Sequences > Step and the injector type/number (visible in Settings > View Parameters).	1. Select Exit . 2. In the Run menu, select Settings > Sequences > Step and modify the injection mode. Otherwise, select Settings > Setup > Alignment and modify the injector type/number. See 'Realign injectors' on page 145 .								
Warning: Run Error Parameters Method [x] Syringe Volume mismatch in Method [x]	Inconsistency between the syringe volume set during installation and the syringe volume set in the method.	<table><tr><td colspan="2">1. Select Exit.2. Do what follows:</td></tr><tr><th>If the installed syringe...</th><th>Then...</th></tr><tr><td>is correct</td><td>1. in the Run menu select Settings > Methods > (method) OK > General and check the value of Syringe Volume. If not correct, modify it and save the method. 2. in the Run menu select Settings > View Parameter and check the value of Syringe Volume. If not correct, execute Syringe Installation procedure taking care to select the correct syringe volume. See 'Syringe management' on page 131.</td></tr><tr><td>is not correct</td><td>Replace the syringe. See 'Syringe management' on page 131.</td></tr></table>	1. Select Exit .2. Do what follows:		If the installed syringe...	Then...	is correct	1. in the Run menu select Settings > Methods > (method) OK > General and check the value of Syringe Volume. If not correct, modify it and save the method. 2. in the Run menu select Settings > View Parameter and check the value of Syringe Volume. If not correct, execute Syringe Installation procedure taking care to select the correct syringe volume. See 'Syringe management' on page 131 .	is not correct	Replace the syringe. See 'Syringe management' on page 131 .
1. Select Exit .2. Do what follows:										
If the installed syringe...	Then...									
is correct	1. in the Run menu select Settings > Methods > (method) OK > General and check the value of Syringe Volume. If not correct, modify it and save the method. 2. in the Run menu select Settings > View Parameter and check the value of Syringe Volume. If not correct, execute Syringe Installation procedure taking care to select the correct syringe volume. See 'Syringe management' on page 131 .									
is not correct	Replace the syringe. See 'Syringe management' on page 131 .									

13.6 Preventive maintenance warnings

Error message	Cause	Solution
FAIL. [Injector position] injector septum/septa is/are not suitable for injection. Check septum/septa and retry.	Injector septa test failed. Septa too tightened or too hard for the injection.	<ol style="list-style-type: none"> 1. Slightly unscrew the injector septum/septa. 2. Select Retry to repeat the Injector septa test. 3. If the problem persist, change the injector septa/septum and retry the Injector septa test. 4. If the problem persist, contact your Customer Representative.
	Injector septa test failed: the test was execute with injector not at room temperature and/or with pressure inside the injector.	<ol style="list-style-type: none"> 1. Cool down the injector to room temperature. 2. Set injector pressure to OFF/0 bar. 3. Select Retry to repeat the Injector septa test. 4. If the problem persists, change the injector septa/septum and retry the Injector septa test. 5. If the problem persists, contact your Customer Representative.
Injection limit reached.	Injection counter reached injection limit set.	Select Ok and replace the injector septum.
Penetrations limit reached. Replace your syringe.	Penetration counter reached penetration limit set.	Select Replace and replace the syringe.
Plunger strokes limit reached. Replace your syringe.	Plunger strokes counter reached plunger strokes limit set.	Select Replace and replace the syringe.
Syringe is expired. Replace your syringe.	Syringe reached syringe expire date set.	Select Replace and replace the syringe.
Turret assembly requires maintenance. Contact your Service Rep. This is a predictive maintenance warning. You can continue to operate the device safely in the meantime. ([Code])	Turret testing failed. Turret worn.	Contact your Customer Representative indicating the error message and the error code.
Missing Turret baseline.	Turret baseline is missing (required to perform turret testing).	<ol style="list-style-type: none"> 1. Select Execute. 2. Select Continue.
Switch off and on the autosampler to perform the preventive	The autosampler has not been turned off and on for more than 24 hours.	Switch the autosampler off and then switch it on, or select Cancel to close the message.

Error message	Cause	Solution
maintenance task.		

14. Troubleshooting and errors in Liquid mode

14.1 Troubleshooting

14.1.1 Analyzer reproducibility issue

Note: the indications given below are for gas chromatography applications and may not be suitable for other analysis techniques.

Problem	Cause	What to do
Peak area or peak retention times are not reproducible	Injection depth	Ensure the needle is set to the correct depth for the analyzer. This is particularly important if the unit has just been installed or moved.
	Leak in the inlet septum	<ol style="list-style-type: none"> 1. If there is a leak in the inlet septum, change it. 2. If the inlet septum has experienced less than 100 injections, ensure what follows to avoid premature septum wear: <ul style="list-style-type: none"> ○ the syringe needle is straight ○ the septum retainer nut is not too tight ○ the syringe has been installed correctly
	Syringe dirty or worn	Clean the syringe according to the manufacturer's cleaning instructions.
	Sample volume in the vial is not correct.	Adjust the sample volume so that it corresponds to half of the vial nominal volume.
	Sample vials are not closed properly.	Ensure the vial caps are not loosen.
	Sample instability	Check the stability of the samples, and ensure they are not sensitive to light or heat.
	Sample size variation and syringe probably imprecise or worn	Change the syringe. See 'Syringe management' on page 173 .
	Air bubbles inside the syringe	<ol style="list-style-type: none"> 1. If the samples is viscous, dilute it in a proper low-viscosity solvent. 2. In the Run menu, select Settings > Methods > (method) OK > General and modify the fill speed and the viscosity delay. 3. Select Exit: the selected method page appears. 4. Select Injection and modify the injection speed. 5. Increase the pull-ups.
Contamination or ghost peak issue	Drop on waste vial cap or on the vial locator	Change the waste vial septum with a pre- slit septum or with a low maintenance septum.
	Small fragments of the solvent cap septum dissolved in the sample.	<ol style="list-style-type: none"> 1. Inject several blanks to check for ghost peaks persistence or elimination. 2. Check what follows: <ul style="list-style-type: none"> ○ the syringe needle has not burrs that can cause septum cutting ○ the septum resistance to the solvent in use ○ the septum is flat and not damaged
	Contamination in sample vials	Use new or clean vials and store new vials in a contaminant- free place.

Problem	Cause	What to do
	Sample instability	Check the stability of the samples, and ensure they are not sensitive to light or heat. Ensure they are stored correctly.
	Drop on waste vial cap or on the vial locator	Change waste vial septum with a pre- slit septum or with a low maintenance septum.
Peak size issue	Leak in the analyzer system	<ol style="list-style-type: none"> 1. Replace the injector septum and check the fitting to avoid leaks. 2. If the septum has experienced less than 100 injections, ensure what follows to avoid premature septum problems: <ul style="list-style-type: none"> ○ the syringe needle is straight ○ the septum retainer nut is not too tight ○ the syringe has been installed correctly
	Sample instability	Check the stability of the samples, and ensure they are not sensitive to light or heat. Ensure they are stored correctly.
	The vials are not closed properly	Ensure the vial caps are not loosed.
	Drop on waste vial cap or on the vial locator	Change waste vial septum with a pre- slit septum or with a low maintenance septum.
Carry over issue	Type or number of washes are inappropriate.	Check if the washing modalities are suitable for your specific application.
	Solvent almost finished	<ul style="list-style-type: none"> • Top up the solvent vial. • If the waste vial is almost full, replace it or empty it.
	Syringe dirty or worn	Clean the syringe according to the manufacturer's cleaning instructions.
	Samples are immiscible.	In this case the washes should not rinse properly the syringe. Increase the wash number or use solvents that can rinse different types of samples.
	Drop on waste vial cap or on the vial locator	Change waste vial septum with a pre- slit septum or with a low maintenance septum.
	Damaged solvent septa: not compatible with solvent (e.g. acetone).	Use low maintenance septa.
No signal/no peak	Syringe plunger is not working properly.	<ol style="list-style-type: none"> 1. Ensure the plunger is secured by the plunger locker. 2. Ensure the syringe needle is not plugged. If needed, replace or clean the syringe.
	Sample volume is insufficient	Top up the sample volume in the vial or change the needle draw depth inside the sample vial. See 'Set the characteristics of the rack to be used' on page 172 .
	Sample is viscous	In the Run menu, select Settings > Methods > (method) OK > General and modify the fill speed and the viscosity delay or dilute the sample in a proper low viscosity solvent.
	Injection depth too high	Ensure the needle goes deep enough into the analyzer injector.
	Injection in the wrong injector	Select the correct injection mode in the Setup menu and Sequences page.

Problem	Cause	What to do
	Drop on waste vial cap or on the vial locator	Change waste vial septum with a pre- slit septum or with a low maintenance septum.

14.1.2 Communication with analyzer issues

Problem	Cause	What to do
The autosampler starts also if the analyzer is "Not Ready".	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The autosampler doesn't start.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The analyzer doesn't start the analysis when the autosampler injects.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The analyzer starts the analysis before that autosampler injects.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.

14.1.3 Usage

Below is a list of possible problems happening and their potential solutions. Some of these are hypothetical as these problems have not occurred.

Problem	Cause	What to do
During the instrument research procedure, no instrument appears on the list	Firewall is blocking HTA Monitor	Be sure that the firewall is allowing HTA Monitor to work. In details, be sure that HTA Monitor could communicate on both Public and Private net. For example, with Windows 10: 1. Reach Windows Defender Firewall (Control Panel) 2. Select Allow an application... 3. Select both Public and Private for HTA Monitor app 4. Retry

14.1.4 Configuration issues

Problem	Cause	What to do
Bad alignment on the injector position(s)	Injector adapter is missing.	Install the proper injector adapter.
	Injector alignment is needed.	Execute Injector alignment . See 'Realign injectors' on page 181 .
The autosampler injects in the wrong injector	Injection Mode in Step is wrong.	In the Run menu, select Settings > Sequences > Step to set the correct injector. See 'Step' on

Problem	Cause	What to do
		page ²⁸⁸ .
	Injection Mode in Run is wrong.	In the Run menu, select Settings > Setup > Configuration > Run to set the correct injector. See 'Run' on page ²⁹⁰ .
Missing: front injector alignment	The front injector alignment has not been performed after having changed the analyzer brand or model.	In the Run menu, select Settings > Setup > Alignment > Front Injector to align the autosampler to the front injector. See 'Realign injectors' on page ¹⁸¹ .
Missing: rear injector alignment	The rear injector alignment has not been performed after having changed the analyzer brand or model.	In the Run menu, select Settings > Setup > Alignment > Rear Injector to align the autosampler to the rear injector. See 'Realign injectors' on page ¹⁸¹ .
Missing: touch & plunger zero	The Touch & Plunger Zero procedure has not been performed.	In the Run menu, select Settings > Setup > Alignment > Touch & Plunger Zero to perform the alignment procedure. See 'Reset the position of the vial' on page ¹⁸³ .
WARNING Check Touch Sensor then execute Touch & Plunger Zero.	A syringe with a longer needle is installed and touch and plunger zero procedure is not performed.	<ol style="list-style-type: none"> 1. Select Continue. 2. Select Touch Plunger Zero. 'Reset the position of the vial' on page ¹⁸³ .

14.1.5 Consumables/syringe issues

Element	Problem	Cause	What to do
Syringe	Syringe needle damaged.	Use of incorrect vials, syringes, needles or injection septum or incorrect syringe installation/configuration.	<ol style="list-style-type: none"> 1. Identify the cause of the damage and the position in which it occurs. 2. Ensure the alignment on the sample vials and on the injector/s is correct. 3. Check that the septa of sample vials and injectors are suitable and correctly pierced. 4. Check that there are no obstacles in the vials, such as an insert. 5. Replace the syringe. See 'Syringe management' on page ¹⁷³.
	Syringe is sucking up air.	Too much sample is being drawn at one time and is causing cavitation.	Use pre-slit septa. Use multiple injections if possible.
		Fill speed too high.	In the Run menu, select Settings > Methods > (method) OK > General and modify the fill speed.
		Tray depth not proper (needle tip is not below liquid level).	Correct tray depth. See 'Set the characteristics of the rack to be used' on page ¹⁷² .
	Needle breaking.	Needle not set up	Align vial locator to needle tip (see 'Align

Element	Problem	Cause	What to do
		correctly.	the vial locator to the needle tip' on page ¹⁷⁶ and execute Touch & Plunger Zero (see 'Reset the position of the vial' on page ¹⁸³).
	Needle obstructed.	Pieces of septum inside needle.	Check quality of septum of sample vial and injector.
		Pieces of septum inside needle.	Check that sample vials are not overcrimped/overtighten: this may result in septum damaging.
		Sample full of dirt.	Check sample quality or modify tray depth to draw at higher level.
	Plunger does not touch the bottom of the syringe barrel. If the plunger does not rest at the bottom of the syringe barrel the unit will not operate correctly and will inject air into the system (the volume of sample may also be inaccurate).	Plunger locker inserted wrong way.	Insert plunger locker in the correct way and execute Touch & Plunger Zero . See 'Reset the position of the vial' on page ¹⁸³ .
		Incorrect Touch & Plunger Zero / Syringe Installation .	Repeat Touch & Plunger Zero or Syringe Installation .
		Incompatible syringe.	Change syringe.
Sample vials	Unit picks up and drops vials.	Safety lock and/or needle height regulator aren't installed properly.	Align the vial locator to the needle tip and execute Touch & Plunger Zero . See 'Reset the position of the vial' on page ¹⁸³ .
		Vial locator and its two steel bars can't move along the syringe needle.	Check if the safety lock and the needle height regulator are fixed correctly; they should not be fixed in a non-parallel position.
		Septum is too tight.	Check if crimping or screwing is performed correctly.
	Vial locator wet and/or liquid on solvent caps.	Waste is without septum or with not pre-slit septum.	Use pre-slit septum.
Solvent vials	Solvent vials caps do not stay in the correct position.	Solvent vial septa have been put upside down.	Remount septa. Red layer should face down.
	Plastic parts inside solvent vials.	Solvent vial septa have been used for too many times.	Replace solvent vial septa.
	Needle touches vial bottom during wash.	Missing septa.	Remount septa.
		Touch & Plunger Zero missing.	Execute Touch & Plunger Zero . See 'Reset the position of the vial' on page ¹⁸³ .
		Waste vial in different configuration (septum/cap) compared to solvent vial	Install the same type of septum / cap on Solvent and waste vial and execute Touch & Plunger Zero . See 'Reset the position of the vial' on page ¹⁸³ .

Element	Problem	Cause	What to do
		configuration.	
		Brand/type of vial has changed since last execution of touch plunger zero.	Execute Touch & Plunger Zero . See ' Reset the position of the vial ' on page 183.
	Vial locator wet and/or liquid on solvent caps.	Waste is without septum or with not pre-slit septum.	Use pre-slit septum.

14.2 Start up errors

14.2.1 Generic start-up errors

Error message	Cause	Solution
<ul style="list-style-type: none"> Plunger error: Code 1 1 Plunger error: Code 1 2 Plunger error: Code 1 4 Plunger error: Code 2 1 Plunger error: Code 2 2 Plunger error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>Syringe is damaged or not installed properly.</p>	<p>Remove the obstacle and press Retry.</p> <p>Remove the syringe and press Retry; then, execute the syringe installation procedure properly. See 'Syringe management' on page 173.</p>
<ul style="list-style-type: none"> Needle error: Code 1 1 Needle error: Code 1 2 Needle error: Code 1 4 Needle error: Code 2 1 Needle error: Code 2 2 Needle error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>Instrument is in shipping position.</p>	<p>Remove the obstacle and press Retry.</p> <p>Press Disengage, set the autosampler in operating position and press Exit.</p>
<ul style="list-style-type: none"> Turret error: Code 1 1 Turret error: Code 1 2 Turret error: Code 1 4 Turret error: Code 2 1 Turret error: Code 2 2 Turret error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p>	<p>Remove the obstacle and press Retry.</p>
<ul style="list-style-type: none"> Tray error: Code 1 1 Tray error: Code 1 2 Tray error: Code 1 4 Tray error: Code 2 1 Tray error: Code 2 2 Tray error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>A sample vial fell between frontal panel and rack support.</p> <p>The rack is not positioned properly.</p>	<p>Remove the obstacle and press Retry.</p> <p>Press Disengage to open the tray and remove the sample vial; press Exit to retry the start-up.</p> <p>Press Disengage to open the tray and reposition the rack correctly; press Exit to retry the start-up.</p>
Unit power up in Service mode: a PIN code is required.	Problem during instrument start-up.	Switch off and on the autosampler to retry the start-up.
Error: "Timeout communication from peripheral: plunger / needle / tray / turret".	Problem during instrument start-up.	Switch off and on the autosampler to retry the start-up.

14.2.2 Safety lock errors

Error message	Cause	Solution
Is safety lock in operating position? If NO, switch off and check user manual.	The syringe location has not been prepared.	<ol style="list-style-type: none"> Switch the autosampler off. See 'Prepare the syringe location' on page 62.

Error message	Cause	Solution
Is safety lock in operating position? If NO, switch off and check user manual.	The touch sensor is not working correctly.	<ol style="list-style-type: none"> 1. Ensure to have prepared the syringe location properly. See 'Prepare the syringe location' on page 62. 2. If the problem persists, select Yes and then Run menu > Settings > Setup > Manual Operation. 3. Select the needle and raise it. 4. Ensure that the vial locator can move up and down freely and that the safety lock is higher than the needle height regulator. 5. Select Exit.

14.2.3 Vial locator errors

Error message	Cause	Solution
Attention! Vial locator in critical position! Press Disengage and move the vial locator to a safe position, then press Exit.	Unexpected switching off of the autosampler while the vial locator was picking up/replacing a vial.	<ol style="list-style-type: none"> 1. Select Disengage: the Manual Operation page appears. 2. Select the needle and raise it. 3. If present, remove the vial from the vial locator and place it back to its original position. 4. Check that all the vials are correctly positioned in the rack. 5. Select Exit.

14.2.4 Memory errors

Error message	Cause	Solution
Error Reading Border mark-out!	The autosampler cannot read the Border Mark Out parameters.	<ol style="list-style-type: none"> 1. Select Continue: the Service page appears. 2. Contact your Customer Representative to solve the problem.
Error Reading Methods!	The autosampler cannot read the parameters of the saved methods.	<ol style="list-style-type: none"> 1. Select Continue. 2. Switch the autosampler off and on. 3. If the error message appears again, contact your Customer Representative to solve the problem.
Error Reading Sequence!	The autosampler cannot read the parameters of the saved sequences.	
Error Reading Features!	The autosampler cannot read the features.	<ol style="list-style-type: none"> 1. Select Continue. 2. Switch the autosampler off and on. 3. If the error message appears again, contact your Customer Representative to solve the problem.

14.3 Run issues / errors

14.3.1 Generic run errors

Error message	Cause	Solution
<ul style="list-style-type: none"> Plunger error: Code 4 2 Plunger error: Code 4 4 Plunger error: Code 4 8 Plunger error: Code 10 7 	Plunger locker is in incorrect position.	Install correctly the plunger locker and execute the Syringe Installation procedure. See ' Syringe management ' on page 173.
	The syringe was installed without to perform completely the Syringe Installation procedure.	Execute completely the Syringe Installation procedure. See ' Syringe management ' on page 173.
	Syringe is not installed correctly.	Install correctly the syringe and perform Syringe Installation procedure. See ' Syringe management ' on page 173.
	Damaged syringe or worn out syringe.	Change the syringe and perform Syringe Installation procedure. See ' Syringe management ' on page 173.
<ul style="list-style-type: none"> Needle error: Code 4 2 Needle error: Code 4 4 Needle error: Code 4 8 Needle error: Code 10 7 	Safety lock or needle height regulator loosen/not fasten properly.	Fasten completely the loosen part and press Retry ; to align the vial locator at the needle tip may be required. See ' Align the vial locator to the needle tip ' on page 176.
	Microinserts are present inside the sample vials and vial depth is too high.	Install a new tray with the proper vial depth or change the vial depth in the tray in use. See ' Racks management ' on page 172.
	The vial locator is not aligned with needle tip, therefore the syringe hits the vial locator when inserting the needle in solvent vial or in the injector.	Align the vial locator at the needle tip. See ' Align the vial locator to the needle tip ' on page 176.
	Injector depth too high, not reachable with the syringe needle in use.	Change injector depth according to the syringe needle in use or change the syringe with a syringe with longer needle.
<ul style="list-style-type: none"> Turret error: Code 4 2 Turret error: Code 4 4 Turret error: Code 4 8 Turret error: Code 10 7 	An obstacle does not allow the motor to run for all its possible movement.	Remove the obstacle and press Retry .
	The anti-friction sticker needs to be cleaned.	Wipe the anti-friction sticker with a soft, lint-free or microfiber cloth dampen with distilled water.
<ul style="list-style-type: none"> Tray error: Code 4 2 Tray error: Code 4 4 Tray error: Code 4 8 Tray error: Code 10 7 	A sample vial fell between frontal panel and rack support.	Press Disengage to open the tray and remove the sample vial.
	The rack is not positioned properly.	Press Disengage to open the tray and correct rack position.

14.3.2 Touch sensor error

Error message	Cause	Solution
WARNING Check Touch Sensor!	The touch sensor is off and/or it is not working correctly.	<ol style="list-style-type: none"> 1. Select Disengage: the Manual Operation page appears. 2. Select the needle and raise it. 3. Ensure that the vial locator can move up and down freely and that the safety lock is higher than the needle height regulator. 4. Select Exit.

14.3.3 Vial detection problem

Error message	Cause	Solution						
Warning: Obstacle found over [x]! Try to remove the obstacle	The autosampler found an obstacle over the indicated position.	<div>1. Remove the obstacle.</div> <div>2. To open the Manual Operation page, select Disengage.</div> <div>3. Select the desired item and move its motor and unlock the system.</div> <div>4. To restart the automatic run from the point of interruption select Retry. Otherwise, select Exit.</div> <div>5. If error persists, execute Touch & Plunger Zero. See 'Reset the position of the vial' on page 183.</div>						
Warning: [x] missing! Insert the vial, then press Retry	Missing vial in the indicated position.	<div>Do what follows:</div> <table><tr><th>If you...</th><th>Then...</th></tr><tr><td>can access the desired position</td><td><div>1. Insert the missing vial.</div><div>2. Select Retry.</div></td></tr><tr><td>cannot access the desired position</td><td><div>1. To open the Manual Operation page, select Disengage.</div><div>2. Select the desired item, move its motor and unlock the system.</div><div>3. Insert the missing vial.</div><div>4. Select Exit: automatic run restarts from the point of interruption.</div></td></tr></table> <div>If error persists, execute Touch & Plunger Zero. See 'Reset the position of the vial' on page 183.</div>	If you...	Then...	can access the desired position	<div>1. Insert the missing vial.</div> <div>2. Select Retry.</div>	cannot access the desired position	<div>1. To open the Manual Operation page, select Disengage.</div> <div>2. Select the desired item, move its motor and unlock the system.</div> <div>3. Insert the missing vial.</div> <div>4. Select Exit: automatic run restarts from the point of interruption.</div>
If you...	Then...							
can access the desired position	<div>1. Insert the missing vial.</div> <div>2. Select Retry.</div>							
cannot access the desired position	<div>1. To open the Manual Operation page, select Disengage.</div> <div>2. Select the desired item, move its motor and unlock the system.</div> <div>3. Insert the missing vial.</div> <div>4. Select Exit: automatic run restarts from the point of interruption.</div>							
Bad alignment on the sample vials OR solvent vials OR waste vialOR oven cover OR oven locations: <div><div>• autosampler touches more times the vial before recognizing it</div><div>• Obstacle Found error appears</div></div>	Touch & Plunger Zero is needed.	Execute Touch & Plunger Zero . See 'Reset the position of the vial' on page 183 .						

Error message	Cause	Solution
Vial not found. Insert proper vial in position C, then continue	The vial in the oven is missing or not inserted in position C.	Insert the vial in position C.
	The vial is a 10 ml or 6 ml, but the spacer used is the one for liquid mode (the black one).	Get the proper spacers kit to work with 10 ml vials or 6 ml vials.

14.4 Mismatch errors

Error message	Cause	Solution						
Change [x] parameter. Unknown value.	Inconsistency between the parameter and the Protocol version set.	Contact your Customer Representative.						
Please add a volume for Fill Volume or set to 0 Pull Up Strokes	Inconsistency between the number of pull up strokes and the fill volume set.	<div><div>1. Select Continue: the Sample page appears.</div><div>2. Do what follows:<table><tr><th>If you want to...</th><th>Then...</th></tr><tr><td>avoid the pull up strokes</td><td>set Pull Up Strokes to 0.</td></tr><tr><td>perform the pull up strokes</td><td>set Fill Volume to a value different from 0.</td></tr></table></div><div>3. Select Exit > Exit and Save if necessary.</div></div>	If you want to...	Then...	avoid the pull up strokes	set Pull Up Strokes to 0 .	perform the pull up strokes	set Fill Volume to a value different from 0 .
If you want to...	Then...							
avoid the pull up strokes	set Pull Up Strokes to 0 .							
perform the pull up strokes	set Fill Volume to a value different from 0 .							
Please add a volume for Fill Volume or set to 0 Sample Washes	Inconsistency between the number of sample washes and the fill volume set.	<div><div>1. Select Continue: the Sample page appears.</div><div>2. Do what follows:<table><tr><th>If you want to...</th><th>Then...</th></tr><tr><td>avoid the sample washes</td><td>set Sample Washes to 0.</td></tr><tr><td>perform the sample washes</td><td>set Fill Volume to a value different from 0.</td></tr></table></div><div>3. Select Exit > Exit and Save if necessary.</div></div>	If you want to...	Then...	avoid the sample washes	set Sample Washes to 0 .	perform the sample washes	set Fill Volume to a value different from 0 .
If you want to...	Then...							
avoid the sample washes	set Sample Washes to 0 .							
perform the sample washes	set Fill Volume to a value different from 0 .							
Please add a volume for Internal Standard or disable Internal Standard	The internal standard is enabled, but its volume is not set.	<div><div>1. Select Continue: the Internal Standard page appears.</div><div>2. Do what follows:<table><tr><th>If you want to...</th><th>Then...</th></tr><tr><td>set the volume</td><td>select Internal Standard Volume and set it to a value different from 0.</td></tr><tr><td>disable the internal standard</td><td>select Enable and set it to NO.</td></tr></table></div><div>3. Select Exit > Exit and Save if necessary.</div></div>	If you want to...	Then...	set the volume	select Internal Standard Volume and set it to a value different from 0 .	disable the internal standard	select Enable and set it to NO .
If you want to...	Then...							
set the volume	select Internal Standard Volume and set it to a value different from 0 .							
disable the internal standard	select Enable and set it to NO .							
Warning: Run Error parameters Step [i]: Set Meth [x] Anatetime>0 or change Injection Synchro	Inconsistency between the Analysis Time and the Injection Synchro set.	<div><div>1. Select Abort.</div><div>2. In the Run menu select Settings > Methods and modify the methods setting Analysis Time to a value higher than 0.</div></div>						
Warning: Run Error parameters Step [x]: Injection Mode mismatch	Inconsistency between the injection mode set in Settings > Sequences > Step and the injector type/number (visible in Settings > View Parameters).	<div><div>1. Select Exit.</div><div>2. In the Run menu, select Settings > Sequences > Step and modify the injection mode. Otherwise, select Settings > Setup > Alignment and modify the injector type/number. See 'Realign injectors' on page 181.</div></div>						
Warning: Run Error Parameters	Inconsistency between the syringe volume set during installation and the syringe volume set in the method.	<div><div>1. Select Exit.</div><div>2. Do what follows:</div></div>						

Error message	Cause	Solution	
Method [x] Syringe Volume mismatch in Method [x]		If the installed syringe...	Then...
		is correct	<ol style="list-style-type: none"> 1. in the Run menu select Settings > Methods > (method) OK > General and check the value of Syringe Volume. If not correct, modify it and save the method. 2. in the Run menu select Settings > View Parameter and check the value of Syringe Volume. If not correct, execute Syringe Installation procedure taking care to select the correct syringe volume. See 'Syringe management' on page 173.
		is not correct	Replace the syringe. See ' Syringe management ' on page 173 .
Step [i]: Change [x] parameter. Unknown value	The parameter is not compatible with the 'protocol' parameter.	Contact your Customer Representative to solve the problem.	
Step [i]: Tray in sequence is not installed	The tray chosen in the sequence menu is not the installed one.	Change the tray in Tray Type menu or change the tray in Sequence menu.	
Step [i]: Syringe Volume mismatch in Method [Y]	The method is created for a syringe different than the installed one.	Install the proper syringe or change the syringe volume in the method.	
Rear injector not installed	The autosampler is set to work just with one (front) injector and it can't inject in the second (rear) injector.	Change the injector in the sequence or add and align the rear injector.	

14.5 Preventive maintenance warnings

Error message	Cause	Solution
FAIL. [Injector position] injector septum/septa is/are not suitable for injection. Check septum/septa and retry.	Injector septa test failed. Septa too tightened or too hard for the injection.	<ol style="list-style-type: none"> 1. Slightly unscrew the injector septum/septa. 2. Select Retry to repeat the Injector septa test. 3. If the problem persist, change the injector septa/septum and retry the Injector septa test. 4. If the problem persist, contact your Customer Representative.
	Injector septa test failed: the test was execute with injector not at room temperature and/or with pressure inside the injector.	<ol style="list-style-type: none"> 1. Cool down the injector to room temperature. 2. Set injector pressure to OFF/0 bar. 3. Select Retry to repeat the Injector septa test. 4. If the problem persists, change the injector septa/septum and retry the Injector septa test. 5. If the problem persists, contact your Customer Representative.
Injection limit reached.	Injection counter reached injection limit set.	Select Ok and replace the injector septum.
Penetrations limit reached. Replace your syringe.	Penetration counter reached penetration limit set.	Select Replace and replace the syringe.
Plunger strokes limit reached. Replace your syringe.	Plunger strokes counter reached plunger strokes limit set.	Select Replace and replace the syringe.
Syringe is expired. Replace your syringe.	Syringe reached syringe expire date set.	Select Replace and replace the syringe.
Syringe requires attention. Syringe plunger hard to move. Clean or replace syringe.	Syringe testing failed. Syringe plunger hard to move.	<ol style="list-style-type: none"> 1. Remove the syringe and check it for plunger stickiness or binding. Slide the syringe plunger up and down a few times. It should move smoothly without sticking and binding. 2. If the syringe plunger is sticky, remove the plunger and clean it with solvent. 3. Select Replace and replace the syringe if necessary. 4. If the problem persists, contact your Customer Representative.
Turret assembly requires maintenance. Contact your Service Rep. This is a predictive maintenance warning. You can	Turret testing failed. Turret worn.	Contact your Customer Representative indicating the error message and the error code.

Error message	Cause	Solution
continue to operate the device safely in the meantime. ([Code])		
Missing Plunger baseline.	Plunger baseline is missing (required to perform syringe testing).	<ol style="list-style-type: none"> 1. Select Execute. 2. Remove the syringe if necessary. 3. Select Continue. 4. Select Continue.
Missing Turret baseline.	Turret baseline is missing (required to perform turret testing).	<ol style="list-style-type: none"> 1. Select Execute. 2. Select Continue.
Switch off and on the autosampler to perform the preventive maintenance task.	The autosampler has not been turned off and on for more than 24 hours.	Switch the autosampler off and then switch it on, or select Cancel to close the message.

15. Troubleshooting and errors in SPME mode

15.1 Troubleshooting

15.1.1 Analyzer reproducibility issue

Note: the indications given below are for gas chromatography applications and may not be suitable for other analysis techniques.

Problem	Cause	What to do
Peak area or peak retention times are not reproducible.	Injection depth	Ensure the fiber is set to the correct depth for the analyzer. This is particularly important if the unit has just been installed or moved.
	Leak in the inlet septum	<ol style="list-style-type: none"> 1. If there is a leak in the inlet septum, change it. 2. If the inlet septum has experienced less than 100 injections, ensure what follows to avoid premature septum wear: <ul style="list-style-type: none"> ○ the fiber needle is straight ○ the septum retainer nut is not too tight ○ the fiber and fiber holder have been installed correctly
	Fiber dirty or worn	<ol style="list-style-type: none"> 1. Replace the fiber. 2. Increase the desorption time (autosampler method). 3. Increase the post cleaning time (autosampler method). 4. Increase the injector temperature (GC setup).
	Vacuum created in sample vial	Reduce sample volume.
	Method parameters are not correct	Check each parameter of the method to be sure they fits your application.
Peak size issue	Injection depth	Ensure the fiber is set to the correct depth for the analyzer. This is particularly important if the unit has just been installed or moved.
	Method parameters are not correct	Check each parameter of the method to be sure they fits your application.
	Leak in the analyzer system	<ol style="list-style-type: none"> 1. Replace the injector septum and check the fitting to avoid leaks. 2. If the septum has experienced less than 100 injections, ensure what follows to avoid premature septum problems: <ul style="list-style-type: none"> ○ the fiber needle is straight ○ the septum retainer nut is not too tight ○ the fiber and fiber holder have been installed correctly
	GC injection liner not suitable	Check the liner geometry according to the recommendations of the GC manufacturer.
	GC related parameters	Any parameter like injector, oven or detector temperature can contribute to bad peak shape. Check the column type and connections (square cut connections). Check the troubleshooting guide of the GC manufacturer.
Carry over issue	Fiber dirty or worn	<ol style="list-style-type: none"> 1. Replace syringe glass and syringe plunger.

Problem	Cause	What to do
		2. Increase the desorption time (autosampler method). 3. Increase the post cleaning time (autosampler method). 4. Increase the injector temperature (GC setup).
No signal/no peak	Syringe is not working properly	Replace syringe glass and syringe plunger.
	Injection depth too high	Ensure the fiber goes deep enough into the analyzer injector.
	Injection in the wrong injector	Select the correct injection mode in the Setup menu and Sequences page.

15.1.2 Communication with analyzer issues

Problem	Cause	What to do
The autosampler starts also if the analyzer is "Not Ready".	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The autosampler doesn't start.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The analyzer doesn't start the analysis when the autosampler injects.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.
The analyzer starts the analysis before that autosampler injects.	Interface cable is not properly connected.	Check if cable is installed properly.
	Wrong analyzer or autosampler settings.	Check installation manual for the specific analyzer.

15.1.3 Usage

Below is a list of possible problems happening and their potential solutions. Some of these are hypothetical as these problems have not occurred.

Problem	Cause	What to do
During the instrument research procedure, no instrument appears on the list	Firewall is blocking HTA Monitor	Be sure that the firewall is allowing HTA Monitor to work. In details, be sure that HTA Monitor could communicate on both Public and Private net. For example, with Windows 10: 1. Reach Windows Defender Firewall (Control Panel) 2. Select Allow an application... 3. Select both Public and Private for HTA Monitor app 4. Retry

15.1.4 Configuration issues

Problem	Cause	What to do
Bad alignment on the injector position(s)	Injector adapter is missing.	Install the proper injector adapter.
	Injector alignment is needed.	Execute Injector alignment . See ' Realign injectors' on page 234 '.
The autosampler injects in the wrong injector	Injection Mode in Step is wrong.	In the Run menu, select Settings > Sequences > Step to set the correct injector. See ' Step' on page 320 '.
	Injection Mode in Run is wrong.	In the Run menu, select Settings > Setup > Configuration > Run to set the correct injector. See ' Run' on page 322 '.
Missing: front injector alignment	The front injector alignment has not been performed after having changed the analyzer brand or model.	In the Run menu, select Settings > Setup > Alignment > Front Injector to align the autosampler to the front injector. See ' Realign injectors' on page 234 '.
Missing: rear injector alignment	The rear injector alignment has not been performed after having changed the analyzer brand or model.	In the Run menu, select Settings > Setup > Alignment > Rear Injector to align the autosampler to the rear injector. See ' Realign injectors' on page 234 '.
Missing: touch & plunger zero	The Touch & Plunger Zero procedure has not been performed.	In the Run menu, select Settings > Setup > Alignment > Touch & Plunger Zero to perform the alignment procedure. See ' Reset the position of the vial' on page 236 '.
WARNING Check Touch Sensor then execute Touch & Plunger Zero.	A syringe with a longer needle is installed and touch and plunger zero procedure is not performed.	1. Select Continue . 2. Select Touch Plunger Zero . See ' Reset the position of the vial' on page 236 '.

15.1.5 Consumables/fiber issues

Element	Problem	Cause	What to do
Fiber	Fiber damaged.	Use of incorrect vials or injection septum or incorrect syringe installation/configuration.	1. Identify the cause of the damage and the position in which it occurs. 2. Ensure the alignment on the oven cover and on the injector/s is correct. 3. Check that the septa of sample vials and injectors are suitable and correctly pierced. 4. Check that there are no obstacles in the vials, such as an insert. 5. Replace the fiber. See ' Fiber management' on page 210 '.
	Fiber breaking.	Fiber not set up correctly.	Align vial locator to fiber tip (see ' Align the vial locator to the needle tip' on page 217 ') and execute Touch & Fiber Zero (see ' Reset the position of the vial' on page 236 ').

Element	Problem	Cause	What to do
		Pieces of septum inside needle.	Check quality of septum of sample vial and injector.
		Pieces of septum inside needle.	Check that sample vials are not overcrimped: this may result in septum damaging.
		Sample full of dirt.	Check sample quality or modify extraction depth to sample at higher level.
Sample vials	Unit picks up and drops vials.	Safety lock and/or needle height regulator aren't installed properly.	Align the vial locator to the fiber tip and execute Touch & Fiber Zero . See ' Reset the position of the vial' on page 236 '.
		Vial locator and its two steel bars can't move along the fiber needle.	Check if the safety lock and the needle height regulator are fixed correctly; they should not be fixed in a non-parallel position.
		Septum is too tight.	Check if crimping or screwing is performed correctly.

15.2 Start up errors

15.2.1 Generic start-up errors

Error message	Cause	Solution
<ul style="list-style-type: none"> Plunger error: Code 1 1 Plunger error: Code 1 2 Plunger error: Code 1 4 Plunger error: Code 2 1 Plunger error: Code 2 2 Plunger error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>Fiber is damaged or not installed properly.</p>	<p>Remove the obstacle and press Retry.</p> <p>Remove the fiber and fiber holder and press Retry; then, execute the syringe installation procedure properly. See 'Fiber management' on page 210.</p>
<ul style="list-style-type: none"> Needle error: Code 1 1 Needle error: Code 1 2 Needle error: Code 1 4 Needle error: Code 2 1 Needle error: Code 2 2 Needle error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>Instrument is in shipping position.</p>	<p>Remove the obstacle and press Retry.</p> <p>Press Disengage, set the autosampler in operating position and press Exit.</p>
<ul style="list-style-type: none"> Turret error: Code 1 1 Turret error: Code 1 2 Turret error: Code 1 4 Turret error: Code 2 1 Turret error: Code 2 2 Turret error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p>	<p>Remove the obstacle and press Retry.</p>
<ul style="list-style-type: none"> Tray error: Code 1 1 Tray error: Code 1 2 Tray error: Code 1 4 Tray error: Code 2 1 Tray error: Code 2 2 Tray error: Code 2 4 	<p>An obstacle does not allow the motor to run for all its possible movement.</p> <p>A sample vial fell between frontal panel and rack support.</p> <p>The rack is not positioned properly.</p>	<p>Remove the obstacle and press Retry.</p> <p>Press Disengage to open the tray and remove the sample vial; press Exit to retry the start-up.</p> <p>Press Disengage to open the tray and reposition the rack correctly; press Exit to retry the start-up.</p>
Unit power up in Service mode: a PIN code is required.	Problem during instrument start-up.	Switch off and on the autosampler to retry the start-up.
Error: "Timeout communication from peripheral: plunger / needle / tray / turret".	Problem during instrument start-up.	Switch off and on the autosampler to retry the start-up.

15.2.2 Safety lock errors

Error message	Cause	Solution
Is safety lock in operating position? If NO, switch off and check user manual.	The syringe location has not been prepared.	<ol style="list-style-type: none"> Switch the autosampler off. See 'Prepare the syringe location' on page 62.

Error message	Cause	Solution
Is safety lock in operating position? If NO, switch off and check user manual.	The touch sensor is not working correctly.	<ol style="list-style-type: none"> 1. Ensure to have prepared the syringe location properly. See 'Prepare the syringe location' on page 62. 2. If the problem persists, select Yes and then Run menu > Settings > Setup > Manual Operation. 3. Select the needle and raise it. 4. Ensure that the vial locator can move up and down freely and that the safety lock is higher than the needle height regulator. 5. Select Exit.

15.2.3 Vial locator errors

Error message	Cause	Solution
Attention! Vial locator in critical position! Press Disengage and move the vial locator to a safe position, then press Exit.	Unexpected switching off of the autosampler while the vial locator was picking up/replacing a vial.	<ol style="list-style-type: none"> 1. Select Disengage: the Manual Operation page appears. 2. Select the needle and raise it. 3. If present, remove the vial from the vial locator and place it back to its original position. 4. Check that all the vials are correctly positioned in the rack. 5. Select Exit.

15.2.4 Memory errors

Error message	Cause	Solution
Error Reading Border mark-out!	The autosampler cannot read the Border Mark Out parameters.	<ol style="list-style-type: none"> 1. Select Continue: the Service page appears. 2. Contact your Customer Representative to solve the problem.
Error Reading Methods!	The autosampler cannot read the parameters of the saved methods.	<ol style="list-style-type: none"> 1. Select Continue. 2. Switch the autosampler off and on. 3. If the error message appears again, contact your Customer Representative to solve the problem.
Error Reading Sequence!	The autosampler cannot read the parameters of the saved sequences.	
Error Reading Features!	The autosampler cannot read the features.	<ol style="list-style-type: none"> 1. Select Continue. 2. Switch the autosampler off and on. 3. If the error message appears again, contact your Customer Representative to solve the problem.

15.3 Run issues / errors

15.3.1 Generic run errors

Error message	Cause	Solution
<ul style="list-style-type: none"> Plunger error: Code 4 2 Plunger error: Code 4 4 Plunger error: Code 4 8 Plunger error: Code 10 7 	Plunger locker is not present.	Install correctly the plunger locker and execute the Fiber/Seal Tip Installation procedure. See 'Fiber management' on page 210 .
	The fiber and fiber holder were installed without to perform completely the Syringe Installation procedure.	Execute completely the Syringe Installation procedure. See 'Fiber management' on page 210 .
	Syringe is not installed correctly.	Install correctly the syringe and perform Syringe Installation procedure. See 'Fiber management' on page 210 .
	Damaged syringe or worn out syringe.	Change the syringe and perform Fiber/Seal Tip Installation procedure. See 'Fiber management' on page 210 .
<ul style="list-style-type: none"> Needle error: Code 4 2 Needle error: Code 4 4 Needle error: Code 4 8 Needle error: Code 10 7 	Safety lock or needle height regulator loosen/not fasten properly.	Fasten completely the loosen part and press Retry ; to align the vial locator at the needle tip may be required. See 'Align the vial locator to the needle tip' on page 217 .
<ul style="list-style-type: none"> Turret error: Code 4 2 Turret error: Code 4 4 Turret error: Code 4 8 Turret error: Code 10 7 	An obstacle does not allow the motor to run for all its possible movement.	Remove the obstacle and press Retry .
	The anti-friction sticker needs to be cleaned.	Wipe the anti-friction sticker with a soft, lint-free or microfiber cloth dampen with distilled water.
<ul style="list-style-type: none"> Tray error: Code 4 2 Tray error: Code 4 4 	A sample vial fell between frontal panel and rack support.	Press Disengage to open the tray and remove the sample vial.
<ul style="list-style-type: none"> Tray error: Code 4 8 Tray error: Code 10 7 	The rack is not positioned properly.	Press Disengage to open the tray and correct rack position.
<ul style="list-style-type: none"> Cover error: Code 7 1 	A vial has fallen during vial transport and blocks the cover to move.	Remove the vial, and check 'Vial transport problem' on page 349 .
<ul style="list-style-type: none"> Shaker error: Code 7 1 	Something is blocking the shaker to move.	Remove the obstacle and press Retry .
<ul style="list-style-type: none"> Oven error: Code 8 1 Oven error: Code 8 2 Oven error: Code 8 3 Oven error: Code 8 4 	The oven cannot reach or keep the desired temperature.	Contact your customer representative.

15.3.2 Touch sensor error

Error message	Cause	Solution
WARNING Check Touch Sensor!	The touch sensor is off and/or it is not working correctly.	<ol style="list-style-type: none"> 1. Select Disengage: the Manual Operation page appears. 2. Select the needle and raise it. 3. Ensure that the vial locator can move up and down freely and that the safety lock is higher than the needle height regulator. 4. Select Exit.

15.3.3 Vial transport problem

Error	Cause	Solution
Falling vials	Sample vials not prepared properly.	<ol style="list-style-type: none"> 1. Check the external status of vial, cap and septa. 2. Be sure they are not overcrimped.
	Vial locator damaged.	Contact customer representative.

15.3.4 Vial detection problem

Error message	Cause	Solution						
Warning: Obstacle found over [x]! Try to remove the obstacle	The autosampler found an obstacle over the indicated position.	<div>1. Remove the obstacle.</div> <div>2. To open the Manual Operation page, select Disengage.</div> <div>3. Select the desired item and move its motor and unlock the system.</div> <div>4. To restart the automatic run from the point of interruption select Retry. Otherwise, select Exit.</div> <div>5. If error persists, execute Touch & Fiber Zero. See 'Reset the position of the vial' on page 236.</div>						
Warning: [x] missing! Insert the vial, then press Retry	Missing vial in the indicated position.	<div>Do what follows:</div> <table><tr><th>If you...</th><th>Then...</th></tr><tr><td>can access the desired position</td><td><div>1. Insert the missing vial.</div><div>2. Select Retry.</div></td></tr><tr><td>cannot access the desired position</td><td><div>1. To open the Manual Operation page, select Disengage.</div><div>2. Select the desired item, move its motor and unlock the system.</div><div>3. Insert the missing vial.</div><div>4. Select Exit: automatic run restarts from the point of interruption.</div></td></tr></table> <div>If error persists, execute Touch & Fiber Zero. See 'Reset the position of the vial' on page 236.</div>	If you...	Then...	can access the desired position	<div>1. Insert the missing vial.</div> <div>2. Select Retry.</div>	cannot access the desired position	<div>1. To open the Manual Operation page, select Disengage.</div> <div>2. Select the desired item, move its motor and unlock the system.</div> <div>3. Insert the missing vial.</div> <div>4. Select Exit: automatic run restarts from the point of interruption.</div>
If you...	Then...							
can access the desired position	<div>1. Insert the missing vial.</div> <div>2. Select Retry.</div>							
cannot access the desired position	<div>1. To open the Manual Operation page, select Disengage.</div> <div>2. Select the desired item, move its motor and unlock the system.</div> <div>3. Insert the missing vial.</div> <div>4. Select Exit: automatic run restarts from the point of interruption.</div>							

Error message	Cause	Solution
Bad alignment on the sample vials OR wash vials OR oven cover OR oven locations: <ul style="list-style-type: none"> • autosampler touches more times the vial before recognizing it • Obstacle Found error appears 	Touch & Plunger Zero is needed.	Execute Touch & Fiber Zero . See 'Reset the position of the vial' on page 236 .
Vial not found. Insert proper vial in position C, then continue	The vial in the oven is missing or not inserted in position C.	Insert the vial in position C.
	The vial is a 10 ml or 6 ml, but the spacer used is the one for liquid mode (the black one).	Get the proper spacers kit to work with 10 ml vials or 6 ml vials.

15.4 Mismatch errors

Error message	Cause	Solution
Change [x] parameter. Unknown value.	Inconsistency between the parameter and the Protocol version set.	Contact your Customer Representative.
Warning: Run Error parameters Step [i]: Set Meth [x] Anatime>0 or change Injection Synchro	Inconsistency between the Analysis Time and the Injection Synchro set.	1. Select Abort . 2. In the Run menu select Settings > Methods and modify the methods setting Analysis Time to a value higher than 0.
Warning: Run Error parameters Step [x]: Injection Mode mismatch	Inconsistency between the injection mode set in Settings > Sequences > Step and the injector type/number (visible in Settings > View Parameters).	1. Select Exit . 2. In the Run menu, select Settings > Sequences > Step and modify the injection mode. Otherwise, select Settings > Setup > Alignment and modify the injector type/number. See 'Realign injectors' on page 234 .

15.5 Preventive maintenance warnings


Error message	Cause	Solution
Injection limit reached.	Injection counter reached injection limit set.	Select Ok and replace the injector septum.
Penetrations limit reached. Replace your syringe.	Penetration counter reached penetration limit set.	Select Replace and replace the syringe.
Turret assembly requires maintenance. Contact your Service Rep. This is a predictive maintenance warning. You can continue to operate the device safely in the meantime. ([Code])	Turret testing failed. Turret worn.	Contact your Customer Representative indicating the error message and the error code.
Missing Turret baseline.	Turret baseline is missing (required to perform turret testing).	1. Select Execute . 2. Select Continue .
Switch off and on the autosampler to perform the preventive maintenance task.	The autosampler has not been turned off and on for more than 24 hours.	Switch the autosampler off and then switch it on, or select Cancel to close the message.

16. Manage unsolved errors

16.1 Deal with unsolved errors

For errors and problems that cannot be solved, contact your Customer Representative.

16.2 Send log files to your Customer Representative

1. In the taskbar, right click  and click **Generate QR Code**: the **QR Code Wizard** window opens.
2. Follow the wizard instructions: the QR Code appears.
3. Scan the QR code using a dedicated QR reader app on your smartphone: a pre-filled email appears.
4. Send the email.

16.3 Tray emergency release

If tray doesn't open automatically (and samples are inside), please proceed as follows:

1. Switch the autosampler OFF.
2. Wait 30 sec.
3. Switch it ON again.
4. Try to open the tray by selecting the tray icon in the home screen. If its not possible to access the home screen, try to access the manual operation screen.
5. If this fails, switch the autosampler OFF.
6. Disconnect it from main power.
7. Pull the tray out gently using both hands. The tray opening must be performed by hand only if it is strictly necessary.

Note: *Do not pull the tray with only one hand.*

Do not pull the tray with the power switched on.

17. Software related errors

17.1 HTA Monitor related errors

Error message	Cause	Solution
HTA Monitor is not active. Start PC and press Continue to use HTA Monitor.	PC is switched off	<ol style="list-style-type: none"> 1. Switch on the PC. Wait for the PC startup. 2. Select Continue.
	HTA Monitor is not running.	<ol style="list-style-type: none"> 1. Start HTA Monitor. 2. Select Continue.
	Communication problem	<ol style="list-style-type: none"> 1. Check that all the autosampler and PC cables are correctly connected and that the autosampler is ON. 2. Open HTA Monitor, select Install and check that instrument appears in the list: if not, press Add new and follow the wizard. 3. Contact your Customer Representative.

17.2 HTA Autosampler Manager related errors

Problem	Cause	What to do
HTA Autosampler Manager cannot connect to the autosampler.	Communication problem	<ol style="list-style-type: none"> 1. Check that all the autosampler and PC cables are correctly connected and that the autosampler is ON. 2. Check that PC is working with a Static (fixed) IP. Contact your network administrator if necessary. 3. Check that no firewall is blocking the HTA Autosampler Manager software nor the ports 20101 (TCP), 20102 (TCP) and 20201 (UDP). Contact your network administrator if necessary. 4. Click Modify network configuration of autosampler and Proceed: the new suggested network parameters for the autosampler appear. 5. Contact your network administrator to check the compatibility of the suggested network parameters with your network. 6. To start the automatic assignment of the suggested network parameters, click Next. 7. If the problem persists, contact your Customer representative.

18. Technical features

18.1 Technical data

18.1.1 Features

Dimensions	<ul style="list-style-type: none"> • Width (closed oven cover): 33 cm (13") • Width (open oven cover): 35 cm (13.8") • Depth with closed tray: 32 cm (12.6") • Depth with open tray: 52 cm (20.5") • Height: 64 cm (25.2")
Weight	10 kg (22 lb)
Clearance from the instrument	10 cm (3.9")

18.1.2 Moving parts

Turret	Maximum rotation: 320°
Needle	Up to 26.7 cm (10.51") vertically
Plunger	Up to 6.54 cm (2.57") in the syringe
Tray	Up to 22 cm (8.66") outside the autosampler

18.1.3 Hardware

Touch screen	4.3" TFT Capacitive Full Viewing Angle						
PC requirements (for HTA Monitor)	<ul style="list-style-type: none"> • Operating system: Windows 7, Windows 8.1, Windows 10, Windows 11. PC Edition only (excluding mobile devices and appliances). PC are expected to run Windows OS with the latest update installed (unless differently specified). • Administrator account (just for SW installation: then, any type of account will be able to run the software). • Additionally required software: Microsoft .NET Framework 4.5.2 • 1 Ethernet connector 10/100 available on PC or on a switch • Minimum screen resolution: 1024x768 • Internet: Internet functionality requires Internet access (fees may apply). • The software install the SQL Express component. Considering this component the PC requirements are: <table border="1"> <thead> <tr> <th>If you are using...</th><th>Then...</th></tr> </thead> <tbody> <tr> <td>Windows 7 or Windows 8.1 O.S. *</td><td> <ul style="list-style-type: none"> • RAM: 2 GB minimum, 4 GB suggested • Disk space: 6 GB • Processor: 1.0 GHz or faster • Further requirements are available at the Microsoft SQL Server 2014 Express webpage </td></tr> <tr> <td>Windows 10 or Windows 11*</td><td> <ul style="list-style-type: none"> • RAM: 2 GB minimum, 4 GB suggested • Disk space: 8 GB • Processor: x64 2.0 GHz or faster • Further requirements are available at the Microsoft </td></tr> </tbody> </table> 	If you are using...	Then...	Windows 7 or Windows 8.1 O.S. *	<ul style="list-style-type: none"> • RAM: 2 GB minimum, 4 GB suggested • Disk space: 6 GB • Processor: 1.0 GHz or faster • Further requirements are available at the Microsoft SQL Server 2014 Express webpage 	Windows 10 or Windows 11*	<ul style="list-style-type: none"> • RAM: 2 GB minimum, 4 GB suggested • Disk space: 8 GB • Processor: x64 2.0 GHz or faster • Further requirements are available at the Microsoft
If you are using...	Then...						
Windows 7 or Windows 8.1 O.S. *	<ul style="list-style-type: none"> • RAM: 2 GB minimum, 4 GB suggested • Disk space: 6 GB • Processor: 1.0 GHz or faster • Further requirements are available at the Microsoft SQL Server 2014 Express webpage 						
Windows 10 or Windows 11*	<ul style="list-style-type: none"> • RAM: 2 GB minimum, 4 GB suggested • Disk space: 8 GB • Processor: x64 2.0 GHz or faster • Further requirements are available at the Microsoft 						

If you are using...	Then...
---------------------	---------

SQL Server 2019 webpage


* Software is tested under English Operative System (ISO/IEC 8859-Part 1 "Latin-1 Western European"), 32-bit (for Windows 7) and 64-bit (Windows 8.1 or higher) updated to the last service pack and/or versions.

• OS general settings:

If you are using...	Then...
---------------------	---------


Windows 7, 8 and 10

1. Reach the following menu: **Control Panel > System and Security > Administrative tool > Task scheduler.**
2. Be sure to have enabled all tasks history:



Windows 11

1. Reach the following menu: **Control Panel > System and Security > Windows tool > Task scheduler.**
2. Be sure to have enabled all tasks history:



• OS special settings with multi-user installation, if you want to use HTA Monitor from any user of the PC:

If you are using...	Then...
---------------------	---------

Windows 10

1. Reach the following menu: **Settings > Accounts > Sign-in options.**
2. Set "Use my sign-in info to autostart finish setting up my device after an update or restart" to **OFF**.
3. Repeat this operation for any user you want to

	If you are using...	Then...
		use.
	Windows 11	<ol style="list-style-type: none"> 1. Reach the following menu: Settings > Accounts > Sign-in options. 2. Set "Use my sign-in info to autostart finish setting up my device after an update" to OFF.
PC requirements (for HTA Autosampler Manager)	<ul style="list-style-type: none"> • Operating system: Windows Vista, Windows 7, Windows 8.1, Windows 10, Windows 11. PC Edition only (excluding mobile devices and appliances). PC are expected to run Windows OS with the latest update installed (unless differently specified). Software is tested under English Operative System (ISO/IEC 8859-Part 1 "Latin-1 Western European"), 32-bit (for Windows 7 or lower) and 64-bit (Windows 8.1 or higher) updated to the last service pack and/or versions. • Administrator account • RAM: 1 GB minimum • Disk space: 100 MB (for installation) • Additionally required software: Microsoft .NET Framework 4.5.2 • 1 Ethernet connector 10/100 available on PC or on a switch • Minimum screen resolution: 1152x864 	

18.1.4 Electrical data

Power supply voltage	100- 240 V AC. Connection cable included.
Power supply voltage fluctuation	Not to exceed $\pm 10\%$ of nominal voltage
Clearance from the instrument	10 cm (3.9")
Frequency	50- 60 Hz
Power connection	1.4 A
Autosampler voltage	24 V DC
Autosampler power consumption	120 W
Overvoltage category	II
Relay specification	<ul style="list-style-type: none"> • $V_{\max} = 24 \text{ V}$ • $I_{\max} = 1 \text{ A}$

18.1.5 Interface specifications

To the analyzer	TTL
To the PC	Ethernet 10/100
Interface cable category	5 or higher
Factory Ethernet configuration of the	<ul style="list-style-type: none"> • IP Address: 192.168.0.207

autosampler	<ul style="list-style-type: none"> • Subnet Mask: 255.255.0.0 • Gateway: 192.168.0.1 • Ports: 20101 (TCP), 20102 (TCP), 20201 (UDP)
--------------------	--

18.1.6 Environmental storage conditions

Note: The following conditions are valid only if the product is stored intact and in its original packaging.

Temperature	From - 20 °C to 60 °C (from - 4 °F to 140 °F)
Relative humidity	From 5% to 80%

18.1.7 Environmental conditions of use

Ambience	Only for use indoors. Pollution degree 2, laboratories and normal indoor environment (IEC CEI 60664- 1).
Room temperature limits	From 15°C to 40°C (from 59°F to 104°F)
Ambient humidity limits	From 20% to 80% (non- condensable)
Altitude	Up to 2 000 m (6 561.7 ft)
Atmosphere	No ATEX. Avoid source of fire close to the instrument.
Permitted cleaning products	<ul style="list-style-type: none"> • For the touch screen and all fragile parts: distilled water. • Other parts of the autosampler: non- abrasive detergent solution
Sound pressure level	Maximum measured level: 60 dBA (below the limits of 85 dBA defined by the regulations in force)

18.1.8 Analytical parameters for Headspace technique

Volume increment/decrement	According to the syringe type
Pull-up strokes	Up to 15 strokes
Sampling withdrawal	<ul style="list-style-type: none"> • Sample volume, air volume and aspiration speed: according to the syringe type • Viscosity time: from 0 to 15 s
Injection speed	According to the syringe type
Waiting time before and after the injection	From 0 to 99 s
Enrichment cycles	Up to 15
Dwell time between enrichment cycles	From 0 to 100min
Syringe temperature	Off; From 40°C to 150°C
Oven positions	6
Oven temperature	Off; From 40°C to 170°C
Shaking method	Orbital shaker
Shaking speed	From very low to very high
Shaking cycles	On/Off, from 0 to 9.9min

18.1.9 Analytical parameters for Liquid technique


Volume increment/decrement	According to the syringe type
Number of pre and post washes	Up to 15
Wash mode	Each injection, each sampler or step
Sample wash	Up to 15
Bubble elimination	Up to 15 strokes
Sampling withdrawal	<ul style="list-style-type: none"> Sample volume, air volume and aspiration speed: according to the syringe type Viscosity time: from 0 to 15 s
Injection speed	According to the syringe type
Waiting time before and after the injection	From 0 to 99 s
Internal Standard Technique mode	No, post or double

18.1.10 Analytical parameters for SPME technique

Extraction mode	Liquid phase / Headspace vapors
Fiber type	10mm, 20mm, Sealed tip
Oven positions	6
Oven temperature	Off; From 40°C to 170°C
Incubation time	From 0 to 999min step 1min
Shaking method	Orbital shaker
Shaking speed	From very low to very high
Shaking cycles	On/Off, from 0 to 9.9min

18.1.11 Sample vial specifications for Headspace and SPME technique

Note: Correct sample vial dimensions are essential for proper operation.




Part	Description	Value
A	Vial height	<ul style="list-style-type: none"> 20ml vial: from 77 mm to 79 mm (from 3.03" to 3.11")

Part	Description	Value
		<ul style="list-style-type: none">• 10ml vial: from 48 mm to 50 mm (from 1.89" to 1.97")• 6ml vial: from 40.2 mm to 42.2 mm (from 1.58" to 1.66")
B	Maximum vial body height	<ul style="list-style-type: none">• 20ml vial: 67 mm (2.64")• 10ml vial: 38 mm (1.50")• 6ml vial: 30.2 mm (1.19")
C	Vial diameter	<ul style="list-style-type: none">• 20ml vial: from 22.4 mm to 23.0 mm (from 0.88" to 0.91")• 10ml vial: from 22.4 mm to 23.0 mm (from 0.88" to 0.91")• 6ml vial: from 21.5 mm to 23.0 mm (from 0.85" to 0.91")
D	Vial top part	<ul style="list-style-type: none">• 20ml vial: 10 mm (0.39")• 10ml vial: 10 mm (0.39")• 6ml vial: 10 mm (0.39")
E	Vial rounded bottom	<ul style="list-style-type: none">• 20ml vial: 2 mm (0.08")• 10ml vial: 2 mm (0.08")• 6ml vial: 2 mm (0.08")

18.1.12 Sample vial specifications for Liquid technique

Note: Correct sample vial dimensions are essential for proper operation.



Part	Description	Value
A	Vial without cap height	From 31.5 mm to 32.5 mm (from 1.24" to 1.28")
B	Maximum vial body height	25 mm (0.98")
C	Vial diameter	From 11.5 mm to 11.8 mm (from 0.45" to 0.46")
D	Vial with cap height	From 33 mm to 35 mm (from 1.30" to 1.38")
E	Cap diameter	From 11.5 mm to 12.5 mm (from 0.45" to 0.49")
F	Cap height	From 5 mm to 6.5 mm (from 0.20" to 0.26")

18.1.13 Cap specifications for Headspace and SPME technique

The autosampler uses sample vials with crimp caps.

The Minimum diameter of the pierceable area is 9mm.

In general, do not use crimp caps more than once for Headspace and SPME analysis.

18.1.14 Septa specifications for Headspace technique

Septum material	Compatible with	Incompatible with	Thickness
PTFE/silicone rubber	PTFE resistance until punctured, then septa will have the compatibility of silicone (alcohol, acetone, ether, DMF, DMSO)	ACN, THF, benzene, chloroform, pyridine, toluene, hexane, heptane	2-3mm

For more detailed compatibility data refer to the manufacturer's recommendations. In general, do not use septa more than once for headspace analysis.

18.1.15 Septa specifications for SPME technique

Septum material	Compatible with	Incompatible with	Thickness
PTFE/silicone rubber	PTFE resistance until punctured, then septa will have the compatibility of silicone (alcohol, acetone, ether, DMF, DMSO)	ACN, THF, benzene, chloroform, pyridine, toluene, hexane, heptane	1.5mm ± 0.2mm

For more detailed compatibility data refer to the manufacturer's recommendations. In general, do not use septa more than once for SPME analysis.

18.1.16 Purge line specifications for Headspace technique

Gas connection	<ul style="list-style-type: none"> Use inlet 1/8" on the rear of the autosampler. The inlet gas has to be regulated by using a pressure regulator (only provided as an optional accessory). The gas should be at room temperature.
Maximum pressure	2 bar (29psi or 200kPa)
Acceptable purge gas	<ul style="list-style-type: none"> Helium (Hydrocarbon free) Nitrogen <p>Note: Purge gas must be 99.9995% pure. It is recommended that high quality traps are used in order to remove hydrocarbons, water and oxygen.</p>


External pressure regulator

This optional accessory is used to regulate the pressure of the gas used for the syringe flushing. For more details, see ['Accessories installation' on page 109](#).

Parameter	Value
Input pressure	Lower than 10 bar (145 psi or 1000 kPa)
Output pressure	2 bar maximum (29psi or 200kPa)

18.2 Electrical connections


18.2.1 Structure of the interface cable connector




Pin	Name	Note
1	+5 Volts	Max. 20 mA
2	FREE3- IN	Signal input (true low); Low- level = 0÷0.5 V Hi- level = 4÷12 V
3	START - NO	Signal out = Relay contact (open by default)
4	GCRDY	Signal input (true high); Low- Level = 0÷0.5 V Hi- level = 4÷12 V
5	GND GC	Ground
6	START - COM	Relay common contact
7	FREE1- NC	Signal out = Relay contact (closed by default)
8	ENDSAMP- NC	Signal out = Relay contact (closed by default)
9	READY	Signal input (true low); Low- Level = 0÷0.5 V Hi- level = 4÷12 V
10	FREE1- NO	Signal out = Relay contact (open by default)
11	ENDSAMPL- NO	Signal out = Relay contact (open by default)
12	FREE2	Signal input (true low); Low- Level = 0÷0.5 V Hi- level = 4÷12 V
13	FREE1- COM	Relay common contact
14	ENDSAMPL- COM	Relay common contact
15	START - NC	Signal out = Relay contact (closed by default)

18.2.2 Circuit diagrams

Interface cable connector:



Typical connection with the analyzer:



18.3 Ordering information

18.3.1 Consumables for Headspace technique

Original HTA consumables are the following:

Name	Code
Vial 6ml: clear, round bottom (100pcs/pk)	1.29.949
Vial 10ml: clear, round bottom (100pcs/pk)	1.29.941
Vial 20ml: clear, round bottom (100pcs/pk)	1.29.942
SIL/PTFE Septa & Crimp caps for headspace vial (100pcs/pk) - (-40°C - +200°C)	1.29.940

18.3.2 Syringes for Headspace technique

The following syringes can be provided by HTA.

Syringe volume	Description	Code	Syringe type and warmer assembly
1 ml	Glass Barrel	3.21.027	Standard version 1.93.841
1 ml	Plunger	3.21.025	
2,5 ml	Glass Barrel	3.21.400	
2,5 ml	Plunger	3.21.300	

18.3.3 Syringe kit for Headspace technique

Name	Characteristic	Code
Syringe kit 1ml for Headspace autosampler	<p>Syringe kit for specific syringe code. See 'Syringes for Headspace technique' on page 402.</p> <p>It contains:</p> <ul style="list-style-type: none"> - mechanical parts to use 1 ml syringe; - N°1 syringe glass; - N°1 syringe plunger. 	1.91.280

18.3.4 Consumables for Liquid technique

Original HTA consumables are the following:

Name	Code
Caps 11 mm (10 pieces)	1.29.945
Sample vials 1.5 ml with snap cap (12x32) (100 pieces)	1.29.961
Sample vials 1.5 ml clear with screw cap and septum (100 pieces)	1.28.064
Caps for solvent vials and waste vials (20 pieces, without septa)	1.29.965
Solvent & waste vial septa - low maintenance (20 pieces)	1.93.882
Solvent vials 10 ml (100 pieces, without caps and septa)	1.29.941
Vial assembly low maintenance (solvent): clear vial + septum + cap (6pc/pk)	1.29.967

18.3.5 Syringes for Liquid technique

The following syringes can be provided by HTA.

Syringe size	Needle kind	Syringe holder	Plunger with PTFE cover	Needle length (mm)	Gauge	Notes	Code
0.5 µl	Removable, conical tip	Standard	No	50	G23	-	1.23.109
5 µl	Fixed, conical tip	Standard	No	50	G26	-	1.23.102
	Fixed, conical tip	Long needles syringe holder	No	80	G26 (reinforced)	-	1.23.110
10 µl	Fixed, conical tip	Standard	No	50	G26	Standard syringe	1.23.103
	Fixed, conical tip	Long needles syringe holder	No	80	G26	-	1.23.104
	Fixed, conical tip	Long needles syringe holder	No	80	G22	-	1.23.105
	Fixed, conical tip	Long needles syringe holder	No	80	G23	-	1.23.106
	Removable, conical tip	Standard	No	51	G26	-	1.23.165
	Removable, conical tip	Standard	No	51	G26	Standard syringe, with Syringe ID	1.96.089
	Fixed, conical tip	Long needles syringe holder	No	80	G26 (reinforced)	-	3.10.013
25 µl	Fixed, conical tip	Standard	Yes	50	G23	-	1.23.161
100 µl	Fixed, conical tip	Long needles syringe holder	No	80	G23	-	1.23.120

18.3.6 Needles for Liquid technique

The following needles can be provided by HTA.

Needle kind	Needle length (mm)	Gauge	Notes	HTA Code
Conical tip	51	G26	-	1.23.170
Conical tip	80	G23	-	1.23.116

18.3.7 Consumables for SPME technique

Original HTA consumables are the following:

Name	Code
20ml SPME vial clear, crimp neck, round bottom (100pcs)	1.29.948
Septa & Crimp Caps for SPME vial (100pcs)	1.29.947

18.3.8 Fibers for SPME technique

Type of used fibers	Fibers/extraction devices
Standard SPME fibers	<p>The supplied fiber holder must be used with suitable fibers/extraction devices. It is required that these fibers/extraction devices are original (i.e. of authorized production) as they are subject to industrial property rights.</p> <p>Refer to the Sigma Aldrich/Supelco/Merck local office for the SPME Fibers. Always use Fibers for automatic use. 2800T can handle SPME Fibers of 1cm or 2cm long.</p> <p>The metal Fibers are preferred if they suit your specific applications.</p>
Other suitable fibers/extraction devices	<p>The supplied fiber holder must be used with suitable fibers/extraction devices. It is required that these fibers/extraction devices are original (i.e. of authorized production) as they are subject to industrial property rights.</p> <p>2800T can handle Sealed tip fibers of 1.1mm diameter and of 1.5mm diameter.</p>

18.3.9 Accessories

Original HTA accessories are the following:

Name	Code
SPME kit: standard fibers	1.91.716
SPME kit: standard fibers (no fibers included)	1.91.737
SPME kit: fibers S1.1/1.5mm (no fibers included)	1.91.780
External Pressure Regulator	3.23.000
Swagelok adapter	1.93.822
Spacers kit for 10ml vial (6pc/pk) for Headspace and SPME technique	1.93.820
Spacers kit for 6ml vial (6pc/pk) for Headspace and SPME technique	1.93.828
System integrity kit (Headspace): verify the syringe integrity at the beginning of each batch in a heuristic approach.	1.93.824
Manual Crimper for 20mm Caps	C00028
Manual De-Crimper for 20mm Caps	C00029

18.3.10 Spare parts

Only use original spare parts provided by HTA.

Name	Code
Ethernet cable	1.93.802
Power cord	1.28.013
Power supply 120 W	1.97.060

Name	Code
Sample standard rack for Headspace and SPME technique	1.91.050
Sample standard rack for Liquid technique	1.90.050
Syringe pointer	1.97.138
Vial locator for Liquid technique (green o-ring)	1.97.161
Vial locator for Headspace technique (red o-ring)	3.10.005
Vial locator for Headspace technique and SPME for standard fibers technique (blue o-ring)	1.91.767
Vial locator for SPME for 1.1 fibers / 1.5 fibers technique (yellow o-ring)	1.91.772

18.3.11 HTA Autosampler Manager activation keys

Name	Notes	Code
Activation key standard version (full license)	Standard license to control the autosampler.	1.90.499
Activation key CFR 21 Part 11 version	In addition to the full license, enable the HTA Autosampler Manager software to work according to CFR 21 part 11 compliance.	1.90.503

18.3.12 HTA software connectors

Name	Notes	Code
Driver for 28x0T autosampler: Agilent OpenLAB CDS ChemStation for GC rev. C	Compatible with: <ul style="list-style-type: none"> Openlab CDS Chemstation edition Openlab Chemstation edition 	1.90.465
Driver for 28x0T autosampler: Agilent OpenLab CDS 2	Compatible with: <ul style="list-style-type: none"> Openlab CDS 	1.90.456
Driver for 28x0T autosampler: ICF (for Scion Compass CDS release 4.1 or higher)	Compatible with: <ul style="list-style-type: none"> Scion Compass CDS 	1.90.455

18.4 Disposal

18.4.1 Disposal liability



As specified by the European Directive 2012/19/EU, never dispose of the product as standard urban waste.

The owner is responsible for disposing both these products and any other electronic or electrical equipment via the specific waste collection centers as specified by the government or local public authorities, or by returning the products when purchasing new replacements.

Correct disposal and recycling will contribute to the prevention of potentially harmful consequences for the environment and personal health. Incorrect disposal of the product will be subject to local administrative fines.

For more detailed information on disposal, contact your Customer Representative, visit <https://www.hta-it.com/support/weee.html> or send an email to weee@hta-it.com

19. Conformity

19.1 Declarations

19.1.1 Declaration of conformity

Following is the list of EU Directives and standards the product meets:

Conformity	<ul style="list-style-type: none"> • 2014/35/EU (Low- Voltage Directive) • 2014/30/EU (Electromagnetic Compatibility Directive) • 2006/42/EC (MD - Machinery) • 2011/65/EU (RoHS - Hazardous substances in electrical and electronic equipment) • 2015/863/EU (Delegated Directive amending Annex II to Directive 2011/65/EU regards the list of restricted substances (Restriction of DEHP, BBP, DBP and DIBP))
Reference standards	<ul style="list-style-type: none"> • EN 61010- 1: 2010/A1:2019 (Safety Requirements for electrical equipment for measurement, control and laboratory use Part 1: General requirements) • EN 61010-2-010:2020 (Safety requirements for electrical equipment for measurement, control and laboratory use Part 2-010: Particular requirements for laboratory equipment for the heating of materials) • EN 61010- 2- 081:2020 (Safety Requirements for electrical equipment for measurement, control and laboratory use Part 2- 081: Particular requirements for automatic and semiautomatic laboratory equipment for analysis and other purposes) • EN 61326- 1:2021 (Electrical equipment for measurement, control and laboratory use. EMC Requirements Part 1: General Requirements) • EN 63000:2018 (Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances) • The instrument is designed and manufactured under a quality system registered to ISO 9001

Following is the list of the applicable UK regulations:

Regulations	<ul style="list-style-type: none"> • Electrical Equipment (Safety) Regulations 2016 • Electromagnetic Compatibility Regulations 2016 • The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012
-------------	--

The autosampler has been tested to and complies with IEC CB Scheme, including USA and Canadian national deviations. In particular, the following standards are covered:

Standards	<ul style="list-style-type: none"> • UL 61010- 1 • CSA C22.2 No. 61010- 1 • EN 61010- 1
-----------	--

The autosampler has been tested to and complies with the following Electromagnetic Compatibility (EMC) standards:

Standards	<ul style="list-style-type: none"> • EN 61326- 1, EMC Requirements for Electrical Equipment for Measurement, Control and Laboratory Use – Part 1: General Requirements. • FCC CFR 47 Part 15, Code of Federal Regulations, Title 47 Part 15 Radio Frequency Devices Subpart B, Unintentional radiation. • Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This autosampler has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.
-----------	---

The autosampler is CAN ICES- 3 (A)/NMB- 3(A) following the ICES- 003 regulation.

The autosampler is also compliant with China Administrative Measure on the Control of Pollution Caused by Electronic Information Products regulation (China RoHS):

危险物质名称和浓度表 Table of hazardous substances' name and concentration						
组件名称 Component name	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr6+)	多溴化联苯 (PBB)	多溴联苯醚 (PBDE)
金属零件 Metal parts	0	0	0	0	0	0
塑胶零件 Plastic parts	0	0	0	0	0	0
电子零件 Electronic parts	0	0	0	0	0	0
马达 Motors	0	0	0	0	0	0
电源供应 Power supply	0	0	0	0	0	0
电缆线 Cables	0	0	0	0	0	0
其他部分 Other parts	0	0	0	0	0	0

0: 表明此有毒或有害物质包含在本部分的所有均质材料中，根据 EIP- A, EIP- B, EIP- C 低於極限要求 JS/T11363- 2014

0: Indicates that this toxic or hazardous substance contained in all the homogeneous materials for this part, according to EIP- A, EIP- B, EIP- C, is below the limit requirement in JS/T11363- 2014.

X: 表明此有毒或有害物质包含在本部分的所有均质材料中，根据 EIP- A, EIP- B, EIP- C 超过限額要求 JS/T11363- 2014.

X: Indicates that this toxic or hazardous substance contained in all the homogeneous materials for this part, according to EIP- A, EIP- B, EIP- C, is above the limit requirement in JS/T11363- 2014.

Note: The full EU declaration of conformity text is in the autosampler package.

20. Warranty and service

20.1 Warranty

20.1.1 Terms of warranty

The product is covered by warranty. The full warranty text is available in the product paper in the autosampler package.

20.1.2 Limited Warranty

Due to technology constraints, the touch screen display may have a minor defect. A minor defect refers to a defect which is not considered to be able to substantially degrade the product application or a defect that deviates from existing standards almost unrelated to the effective use of the product or its operation, such as black spot, white spot, bright spot, pinhole, black line, white line, contrast variation, glass defect, polarizer defect, etc. The minor defects do not give right to repair or part substitution under warranty, unless the number of minor defects substantially degrades the product application.

20.2 Service

20.2.1 Customer service

In the event of malfunction, fault or to obtain information and return the product, contact your Customer Representative.

21. Glossary of terms

A

air gap

Air cushion that can be inserted in the Internal Standard technique. The air gap can be between the internal standard and the sample (Post mode), or after both the internal standard and the sample (Double mode).

C

confirmation mode

Among the available injection modes, the autosampler injects the same sample into two injectors to confirm the analysis. The analyzer is started on the first injection.

F

front injector

Among the available analyzer injector, injector that is closest to the user.

R

rear injector

Analyzer injector that is furthest to the user.

S

sequence

A sequence is an ordered series of sample vials to prepare and inject, including the method needed to prepare each vial.

syringe draw depth

Millimeters under the septum of the sample vial cap where the syringe needle must be positioned to aspirate the sample.

system integrity test

It is used to verify the syringe integrity. If this test is enabled, at the beginning of each batch the system performs a check to verify the syringe integrity through a heuristic procedure. This test is carried out using the system integrity tool. The septum of this tool should be replaced periodically. This check can verify if it is better to replace the syringe (barrel and plunger). Please note that this check does not replace the necessary periodic validations you have to program for the instrument qualification. If the test fails, a warning message appears.

T

touch sensor

Sensor in the vial locator that detects that an object has been touched by the vial locator.

H

head space technique

It is an indirect method to inject components with low boiling points from liquid or solid samples, based on the analysis of vapours released from a sample that has been incubated at a constant temperature for a given period of time. -this technique has some advantages compared to other techniques: quick preparation time compared to analysis time, possibility to analyse samples that cannot be analysed in liquid or in solid state, less cross contamination problems and higher sensitivity, as this procedure can be considered as a pre-concentration step before injection.

high throughput

Among the available injection modes, the autosampler injects different samples into two separate injection port. The analyzer is started on the first injection.

V

vial leakage check

This function can only be enabled by the HTA Autosampler Manager (Setup-Options). If this option is enabled, the pressure inside vials is monitored by a heuristic procedure in order to check against anomalous values that are indicative of a vial leakage problem. The Vial leakage test can discriminate samples correctly sealed vs. bad crimping or missing septa. Each sample is marked with PASS or FAIL information in HTA Autosampler Manager Windows.



HTA s.r.l.
Via del Mella, 21
25131 Brescia (BS)
Italy

www.hta-it.com
enquiry@hta-it.com
+39 030 3582920

2800T
Instruction Manual
2800T
Version: h EN 02-Apr-2025
© 2025 HTA s.r.l.